


Medtronic

Onyx Frontier™ DES



Engineered to deliver



Engineered
to deliver



Different
by design




Optimised for
complex PCI



Engineered to deliver

Onyx Frontier DES introduces an enhanced delivery system[†] designed to take the acute performance of Resolute Onyx[™] DES even further.

Enhanced delivery system[†] features:




Dual-flex balloon provides increased flexibility and is comprised of a unique blend of two layers²:

- Inner layer enhances flexibility²
- Outer layer maintains strength²

This results in a **thinner balloon** with the same rated burst pressure (RBP) as Resolute Onyx DES.²



The dual-flex balloon enables a **7.5% lower crossing profile** than Resolute Onyx DES.²



An updated manufacturing process applied to the outer shaft results in a **29% more flexible catheter** (compared to Resolute Onyx DES) for improved deliverability.²

At least

24%


more deliverable than competitive DES^{†3}

[†]Stent delivery system updates were implemented on the 2.0-4.0 mm Onyx Frontier DES diameter.


Different by design




Onyx Frontier DES builds off the legacy of Resolute Onyx DES, featuring the same stent design differentiators that provide the conformability,⁴ visibility,⁵ fast healing,⁶ and size matrix you've come to rely on.



In comparison to laser-cut drug-eluting stents, only Medtronic DES are made from a single wire to enable a fluid range of motion and the conformability needed for **superior strut apposition.**⁴



The platinum-iridium core within Onyx Frontier DES is **more visible** than competitive DES, while enabling greater radial strength with thin struts.^{#7}



The zotarolimus drug inhibits neointimal growth,⁸ **promotes faster healing.**⁶

Only Medtronic offers DES in **2.0 mm to 5.0 mm** sizes to treat the broadest range of coronary vessel diameters.

| Platform | Diameter (mm) | Stent length (mm) | | | | | | | | | | MSID ⁸ (mm) |
|---------------------|---------------|-------------------|----|----|----|----|----|----|----|----|---|------------------------|
| Small vessels | 2.00 | 8 | 12 | 15 | 18 | 22 | 26 | 30 | — | — | — | 3.50 |
| | 2.25 | 8 | 12 | 15 | 18 | 22 | 26 | 30 | 34 | 38 | — | 3.50 |
| | 2.50 | 8 | 12 | 15 | 18 | 22 | 26 | 30 | 34 | 38 | — | 3.50 |
| Medium vessels | 2.75 | 8 | 12 | 15 | 18 | 22 | 26 | 30 | 34 | 38 | — | 4.00 |
| | 3.00 | 8 | 12 | 15 | 18 | 22 | 26 | 30 | 34 | 38 | — | 4.00 |
| Large vessels | 3.50 | 8 | 12 | 15 | 18 | 22 | 26 | 30 | 34 | 38 | — | 5.00 |
| | 4.00 | 8 | 12 | 15 | 18 | 22 | 26 | 30 | 34 | 38 | — | 5.00 |
| Extra-large vessels | 4.50 | — | 12 | 15 | 18 | 22 | 26 | 30 | — | — | — | 6.00 |
| | 5.00 | — | 12 | 15 | 18 | 22 | 26 | 30 | — | — | — | 6.00 |

Four platforms specifically designed to meet the unique needs of each vessel size.



Optimised for complex PCI

An exclusive set of design features and meaningful clinical data inherited from Resolute Onyx DES provide support for your most challenging cases.

Bifurcation PCI



- Other DES feature irregular cell shapes, which may obstruct wire or catheter advancement through the cell's opening
- Round struts create a smooth passage when accessing the side branch, while lowering the propensity to catch⁵

Rounded struts



Onyx Frontier DES

(Cross-section of actual stents)



4.50-5.00 mm expand up to 6.00 mm⁵ while maintaining structural integrity

Left main and other extra-large vessel PCI (4.50-5.00 mm)

- Specifically designed with additional crowns and thicker struts to provide the radial strength needed for extra-large vessels⁵
- Maintains a low profile,⁵ allowing for 5 F compatibility
- ROLEX Registry showed Resolute Onyx DES was safe and effective in left main PCI in a complex patient population⁹

59%

multivessel disease



53%

acute coronary syndrome

30%

diabetic

Low 5.1% TLF, 2.0% TLR, and 1.1% ST at one year



2.00-2.50 mm expand up to 3.50 mm⁵ with minimal foreshortening for tapered and extra-small vessels⁵

Extra-small vessels (2.00-2.50 mm)

- 2.0 mm offers the lowest crossing profile of any DES⁵
- Demonstrated 2% target lesion revascularisation and 0% stent thrombosis at one year in a complex, small-vessel population¹⁰



Onyx ONE Global Trial¹¹ 1,003 Resolute Onyx DES patients studied

COMPLEX PATIENTS

33% AF patients
53% ACS patients
39% diabetic patients

COMPLEX LESIONS

46% moderate to severe calcified lesions
80% B2/C lesions
38 mm average stented length

- Indication is based on the results from the Onyx ONE Global Trial, which evaluated real-world, complex, HBR patients on 1-month DAPT treated with a Resolute Onyx DES or a BioFreedom[™] DCS
- The data is intended to better inform short-DAPT decisions in these patients, including those at high risk of thrombotic events¹¹
- Results showed that Resolute Onyx DES was safe and effective¹¹

[™]Third-party brands are trademarks of their respective owners.

¹Stent delivery system updates were implemented on the 2.0-4.0 mm Onyx Frontier DES diameter.

²Onyx Frontier DES has the same platinum-iridium core as Resolute Onyx DES.

³Stents should not be expanded to a diameter beyond the maximum labeled diameter listed per the IFU. Post-dilation required for overexpansion.

¹Based on bench test data on file at Medtronic. [D00339634 - Test Report for DES Competitive Comparison with Frontier test methods, Rev C, 05-May-2022] May not be indicative of clinical performance. N = 5 DES of each tested: Onyx Frontier DES, Orsiro Mission DES, Resolute Onyx DES, XIENCE Skypoint DES, SYNERGY DES, Ultimaster Tansel DES.

²Based on bench test data on file at Medtronic. [44RD21031-040047 Onyx Frontier Vs Resolute Onyx Balloon Extrusion, Version 1.0, 17-Feb-2022] May not be indicative of clinical performance.

³Based on bench test data on file at Medtronic. [D00339634 - Test Report for DES Competitive Comparison with Frontier test methods, Rev C, 05-May-2022] May not be indicative of clinical performance. N = 7 of each DES tested.

⁴Third-party modeling and analysis. [Mortier MDT-ON14-report-curved-v10-20150220_Onyx_Synergy] Data may not be indicative of clinical performance. Evaluated the following stent platforms: Resolute Onyx DES, Multi-Link 8[™] BMS, SYNERGY[™] DES, XIENCE Alpine[™] DES, and Multi-Link 8 platform.

⁵Based on bench test data on file at Medtronic. [University of Budapest Visibility Testing, V0.1, 28-Sep-2021] May not be indicative of clinical performance.

⁶Roleder T, Kedhi E, Berta B, et al. Short-term stent coverage of second-generation zotarolimus-eluting durable polymer stents: Onyx one-month optical coherence tomography study. *Adv Interv Cardiol.* 2019;15(2):143-150.

⁷Based on bench test data on file at Medtronic. [University of Budapest Visibility Testing, V0.1, 28-Sep-2021; 10166182DOC Competitive Analysis Test Report, Rev AC, 08-Jun-2021] May not be indicative of clinical performance. Stents tested include Resolute Onyx DES, SYNERGY DES, XIENCE Sierra[™] DES, and Orsiro DES.

⁸Yeh RW, Silber S, Chen L, et al. 5-Year Safety and Efficacy of Resolute Zotarolimus-Eluting Stent: The RESOLUTE Global Clinical Trial Program. *JACC Cardiovasc Interv.* February 2017;10(3):247-254.

⁹Tarantini G, et al. The ROLEX Registry (Revascularization Of Left Main With Resolute onyx). Presented at PCR 2022. Investigator-initiated study funded by Medtronic.

¹⁰Cuellas C, et al. Use of a Zotarolimus-eluting stent for small vessel disease (DISCO 9 Study). Presented at PCR 2021.

¹¹Windecker S, Latib A, Kedhi E, et al. Polymer-based or Polymer-free Stents in Patients at High Bleeding Risk. *N Engl J Med.* March 26, 2020;382(13):1208-1218.

This material is for Healthcare Professionals in countries with applicable health authority product registrations.

Important: Always refer to the Instructions For Use (IFU) packaged with the product/e-IFU for complete instructions, indications, contraindications, warnings, and precautions.

UC202210416a EE-onyx-frontier-brochure-en-we-7450604 © 2023 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.

cema.medtronic.com

Medtronic

Medtronic CEMA
Dubai Knowledge Park
Dubai,
Dubai 500638
UAE