

Comparison of NIRS Technology Performance During Acute Hypoxia and Hemorrhagic Shock

Katherine E. Liu, PhD¹, Ulf Borg²

¹Clinical Research, Patient Monitoring, Medtronic; ²Medical Science, Patient Monitoring, Medtronic

Background

Near infrared spectroscopy (NIRS) is used to measure regional oxygen saturation, including both cerebral and somatic oximetry as a marker of perfusion.¹ Understanding cerebral and tissue oxygenation can assist clinicians in managing patients' hemodynamics. For example, cerebral desaturation events during surgery can lead to cognitive decline, major organ morbidity and mortality, renal failure, and stroke.²⁻⁶ In cardiopulmonary bypass surgery, cerebral desaturation events are common, occurring in up to 61% of patients.⁷ Previous Medtronic-independent, randomized prospective studies have shown that use of the INVOS™ monitoring system quickly identifies cerebral desaturation during surgery, allowing for interventions which can reduce patient ICU length of stay, decrease major organ morbidity, and reduce the odds of post operative cognitive decline.^{2,5}

Importantly, although NIRS devices are commercially available from several manufacturers, multiple studies have reported that the devices' clinical performance differs. For example, in a Medtronic-sponsored prospective healthy volunteer study, Tomlin et al found that during controlled desaturation, there were significant differences between Medtronic, Nonin, and Edwards monitors, with the INVOS™ regional oximeter meeting a 20% relative decrease threshold, defined as cerebral

desaturation, in all 10 healthy volunteers, compared to 1 or 2 healthy volunteers with EQUANOX and FORE-SIGHT, respectively.⁸ The INVOS™ regional oximeter also detected a 10% relative decrease in cerebral oxygenation an average of 28 seconds earlier than EQUANOX, and 43 seconds earlier than FORE-SIGHT.⁸ Other Medtronic-independent, prospective, healthy volunteer and patient studies have confirmed that the rates of desaturation and resaturation are significantly different between Medtronic, Nonin, and Masimo NIRS devices.⁹⁻¹¹ Due to differences in sensor designs, algorithms, and device performance, multiple groups have concluded that clinicians should be cautious in generalizing the findings of one NIRS device to all NIRS devices.^{8,11-14} This is supported by multiple prospective, Medtronic-independent studies,^{2,11,14} and by the American Society for Enhanced Recovery and Perioperative Quality Initiative Joint Consensus.¹⁵

The purpose of this study was to compare the performance of multiple NIRS devices in a pre-clinical setting designed to reflect relevant clinical scenarios, including acute hypoxia and hemorrhagic shock.

Methods

Test Devices

4 NIRS devices were tested, including the INVOS™ 7100 Regional Oximeter (Medtronic), EQUANOX 7600 (Nonin), FORE-SIGHT ELITE (Edwards), and Root O3 (Masimo) monitors.

Procedure

Healthy porcine weighing between 30-35lb were assessed. Each animal underwent anesthesia during the procedures. The study was performed at a lab accredited by the Association for Assessment and Accreditation of Laboratory Animal Care, and the protocol was approved by the site Animal Care and Use Committee.

Monitoring Procedures

NIRS sensors were placed between the ears, in front of the Nuchal crest, and over the peri-renal area on the left side. During each clinical scenario, one adult sensor from the Medtronic system and one adult sensor from either the Nonin, Edwards, or Masimo system were placed on the head simultaneously. The somatic sensors were placed over the peri-renal area with the emitters of each sensor facing in opposite directions, to avoid interference.

Acute Hypoxia

To induce acute hypoxia, the ventilator FiO₂ was rapidly decreased from 1.0 (baseline) to 0.08. When the animal reached SpO₂ of 40% ± 10%, this was held for 2 minutes, during which CrSO₂ and RrSO₂ were noted, and then FiO₂ was returned to baseline. This procedure was repeated once for each animal. After the repeat hypoxia procedure, the animal was allowed to recover to baseline.

Hemorrhagic Shock

Hypovolemia was induced by removing a total of 60% of the animal's calculated blood volume, in 3 steps of 20% blood volume reduction. At each step, physiological parameters were noted. To resuscitate the animal, the procedure was reversed, with blood volume restored in the same stepwise manner. At each step, physiological parameters were recorded.

Results


Across the procedures, the Nonin device exhibited signal interference, and was excluded from the analysis.

Acute Hypoxia


Compared to the Edwards device, the Medtronic device measured a large, rapid decrease in CrSO₂ within the first 2 minutes of hypoxia, dropping from 44% to 15% (66% decrease from baseline). In contrast, the Edwards device measured a drop from 52% to 39% (25% decrease from baseline). A similar pattern was observed in RrSO₂ measurements. While the Edwards device detected a decrease in both cerebral and peri-renal oxygen saturation, the magnitude of the change was consistently smaller (Figure 1A). The same trend occurred when comparing the responses of the Masimo and Medtronic devices during acute hypoxia. Relative to the baseline reading for each system, the Medtronic device measured CrSO₂ and RrSO₂ decreases of 70% and 42%, while the Masimo device measured decreases of 32% and 30%, respectively (Figure 1B).

Figure 1. Cerebral and peri-renal oxygen saturation during hypoxia, measured using A) Edwards and Medtronic devices, and B) Masimo and Medtronic devices.

A)



B)




Hemorrhagic Shock


During induction of hemorrhagic shock, the Medtronic device had a more dynamic change in both CrSO₂ and RrSO₂ (24% and 10% from baseline, respectively), compared to the Edwards device (16% and 5% from baseline, respectively) (Figure 2A). The percent decrease in CrSO₂ was also larger for the Medtronic device (48%), compared to the Masimo device (27%) (Figure 2A). The Medtronic device also had a larger magnitude response in CrSO₂ during blood re-infusion, compared to either the Edwards (79% vs 37% increase from baseline, respectively) or Masimo (133% and 55% increase from baseline, respectively) devices (Figure 2B).

Figure 2. Percent change in cerebral oxygen saturation during A) shock, and B) during blood re-infusion.

A)



B)



Discussion

During acute hypoxia and hemorrhagic shock, the Medtronic, Edwards, and Masimo devices all detected changes in cerebral and peri-renal oxygen saturation. Importantly, the Medtronic INVOS™ 7100 regional oximeter showed greater responsiveness in regional oxygen saturation, including both CrSO₂ and RrSO₂, compared to the Edwards and Masimo devices. This is consistent with a previous Medtronic-sponsored animal study of acute hypoxia comparing the INVOS™ 5100C regional oximeter to NIRS devices from other manufacturers, where the Medtronic device demonstrated a faster response time, a greater magnitude change in regional oxygen saturation, and strong correlation to SpO₂ during the transition from normal oxygenation to hypoxia.¹⁶ The results of this study are also consistent with prospective healthy volunteer studies, both MDT-sponsored and MDT-independent, which have reported that the INVOS™ regional oximetry system detects deoxygenation faster and has a more dynamic range compared to other NIRS devices, such as those from Nonin, Edwards, and Masimo.^{8,11}

In circumstances where a patient experiences cerebral or tissue desaturation, the regional oximetry system needs to respond quickly in order to alert the clinician and allow adequate time for intervention. Combined with existing literature, this new data supports the conclusion that other NIRS devices do not respond the same as the INVOS™ regional oximetry system during critical physiological challenges, such as hypoxia or shock.^{9,13,17} Care should be taken when selecting a regional oximeter, to ensure clinicians can achieve the best outcomes for their patients.

The INVOS™ monitoring system should not be used as the sole basis for diagnosis or therapy and is intended only as an adjunct in patient assessment. Reliance on the INVOS™ system alone for detecting cerebral desaturation events is not recommended.

References

1. Triffterer L, Marhofer P, Lechner G, et al. An observational study of the macro- and micro-haemodynamic implications of epidural anaesthesia in children. *Anaesthesia*. Apr 2017;72(4):488-495. doi:10.1111/anae.13746
2. Murkin JM, Adams SJ, Novick RJ, et al. Monitoring brain oxygen saturation during coronary bypass surgery: a randomized, prospective study. *Anesth Analg*. Jan 2007;104(1):51-8. doi:10.1213/01.ane.0000246814.29362.f4
3. Murkin JM, Adams SJ, Pardy E, Quantz M, McKenzie FN, Guo L. Monitoring brain oxygen saturation during coronary bypass surgery improves outcomes in diabetic patients: a post hoc analysis. *Heart Surg Forum*. Feb 2011;14(1):E1-6. doi:10.1532/hcf98.20101065
4. Goldman S, Sutter F, Ferdinand F, Trace C. Optimizing intraoperative cerebral oxygen delivery using noninvasive cerebral oximetry decreases the incidence of stroke for cardiac surgical patients. *The heart surgery forum*. 2004;7(5):E376-E381. doi:10.1532/HSF98.20041062
5. Colak Z, Borojevic M, Bogovic A, Ivancan V, Biocina B, Majeric-Kogler V. Influence of intraoperative cerebral oximetry monitoring on neurocognitive function after coronary artery bypass surgery: a randomized, prospective study. *European journal of cardio-thoracic surgery : official journal of the European Association for Cardio-thoracic Surgery*. 2015;47(3):447-454. doi:10.1093/ejcts/ezu193
6. Slater JP, Guarino T, Stack J, et al. Cerebral Oxygen Desaturation Predicts Cognitive Decline and Longer Hospital Stay After Cardiac Surgery. Article. *The Annals of Thoracic Surgery*. 01/01/January 2009 2009;87(1):36-45. doi:10.1016/j.athoracsur.2008.08.070
7. Subramanian B, Nyman C, Fritock M, et al. A Multicenter Pilot Study Assessing Regional Cerebral Oxygen Desaturation Frequency During Cardiopulmonary Bypass and Responsiveness to an Intervention Algorithm. *Anesth Analg*. 2016;122(6):1786-1793. doi:10.1213/ANE.0000000000001275
8. Tomlin KL, Neitenbach A-M, Borg U. Detection of critical cerebral desaturation thresholds by three regional oximeters during hypoxia: a pilot study in healthy volunteers. *BMC Anesthesiol*. 2017;17(1):6. doi:10.1186/s12871-016-0298-7
9. Fellahi J-L, Butin G, Fischer M-O, Zamparini G, Gérard J-L, Hanouz J-L. Dynamic evaluation of near-infrared peripheral oximetry in healthy volunteers: A comparison between INVOS and EQUANOX. Article. *J Crit Care*. 10/01/October 2013 2013;28(5):881-881. doi:10.1016/j.jcrc.2013.05.004
10. Pisano A, Galdieri N, Iovino TP, Angelone M, Corcione A. Direct comparison between cerebral oximetry by INVOS(TM) and EQUANOX(TM) during cardiac surgery: a pilot study. *Heart, lung and vessels*. 2014;6(3):197-203.
11. Chung J, Ji S-H, Jang Y-E, et al. Evaluation of Different Near-Infrared Spectroscopy Devices for Assessing Tissue Oxygenation with a Vascular Occlusion Test in Healthy Volunteers. Article. *J Vasc Res*. 2020;57(6):341-347. doi:10.1159/000510072
12. Bickler PE, Feiner JR, Rollins MD. Factors affecting the performance of 5 cerebral oximeters during hypoxia in healthy volunteers. *Anesth Analg*. 2013;117(4):813-823. doi:10.1213/ANE.0b013e318297d763
13. Hyttel-Sorensen S, Hessel TW, Greisen G. Peripheral tissue oximetry: comparing three commercial near-infrared spectroscopy oximeters on the forearm. Original Paper. *Journal of Clinical Monitoring and Computing: Including a Specialty Section on Surgical Neuromonitoring*. 04/01/2014;28(2):149-155. doi:10.1007/s10877-013-9507-9
14. Schmidt C, Heringlake M, Kellner P, et al. The effects of systemic oxygenation on cerebral oxygen saturation and its relationship to mixed venous oxygen saturation: A prospective observational study comparison of the INVOS and ForeSight Elite cerebral oximeters. Original Paper. *Effets de l'oxygénation systémique sur la saturation cérébrale en oxygène et ses rapports avec la saturation veineuse mixte en oxygène : étude observationnelle prospective comparant les oxymètres cérébraux INVOS et ForeSight Elite*. 07/01/2018;65(7):766-775. doi:10.1007/s12630-018-1093-3
15. Thiele RH, Shaw AD, Bartels K, et al. American Society for Enhanced Recovery and Perioperative Quality Initiative Joint Consensus Statement on the Role of Neuromonitoring in Perioperative Outcomes: Cerebral Near-Infrared Spectroscopy. *Anesth Analg*. 2020;131(5):1444-1455. doi:10.1213/ANE.0000000000005081
16. Borg U. A Faster and Deeper Response to Hypoxia. Medtronic White Paper. 2017;12/2017-16-PM-0268
17. Moerman A, Vandenplas G, Bové T, Wouters PF, De Hert SG. Relation between mixed venous oxygen saturation and cerebral oxygen saturation measured by absolute and relative near-infrared spectroscopy during off-pump coronary artery bypass grafting. Article. *Br J Anaesth*. 02/01/February 2013 2013;110(2):258-265. doi:10.1093/bja/aes375

©2022 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 10/2022 - US-PM-2200583 - [WF#7830736]

6135 Gunbarrel Avenue Boulder, CO 80301 800.635.5267
[medtronic.com/covidien](https://www.medtronic.com/covidien)

Medtronic