

We take our

patients

pulse

performance

portfolio

partnership

to heart

Nellcor™ pulse oximetry

Never miss a beat

From neonate to the elderly, and virtually every patient in between, count on Nellcor™ pulse oximetry to provide quick and reliable information. We look at every heartbeat to ensure that readings are sensitive and timely, even in the most challenging monitoring conditions.†

†Oxygen saturation accuracy can be affected by certain environmental, equipment, and patient physiologic conditions (as discussed in the operator's manual for the monitor) that influence readings of SpO₂. Please consult the IFU and manual for full safety information.




We take our
patients
pulse
performance
portfolio
partnership
to heart

More than 40 years of innovation

Nellcor™ pulse oximetry has consistently delivered respiratory monitoring solutions that clinicians can count on. We continually build upon our technology with meaningful innovations that make a difference in accuracy, clinician workflow, and assessing the next steps in patient care.

We start with the heart



1983



First to market with **pulse oximetry**

1999



First to market with SatSeconds™ **alarm management** that dynamically assesses intensity of SpO₂ drop before alarming

2002



First to embed an OxiMax™ **digital memory chip** to provide sensor-specific calibration for low saturation accuracy

2003



First AI-enabled deployment of a neural network

2009



First to market with Nellcor™ **saturation pattern detection alert (SPD)** for a real-time indication of repetitive reductions in airflow (RRiA)

2012



First to market with Nellcor™ **respiration rate** software and finger sensor

2013



First company to receive FDA clearance (K123581) for a **motion tolerant** pulse oximeter that is also compliant with ISO 80601-2-61:2011

2022



First to market with a **silicone adhesive sensor OxySoftN** for gentleness and repositionability

We take our

patients

pulse

performance

portfolio

partnership

to heart

Patients are at the heart of all we do

You want the best for patients under your care. That's why when you use pulse oximetry monitoring, you want quick, reliable data no matter the situation. Discover how we met the performance needs of these patients during critical moments.



Meet Sean



Meet Michael



Meet Ann

We take our

patients

pulse

performance

portfolio

partnership

to heart

Patients are at the heart of all we do

Meet Sean

Meet Michael

Meet Ann

Meet Sean. Born prematurely with underdeveloped lungs, he needed to be monitored right after delivery for respiratory insufficiency. The labor and delivery team ensured that the Nellcor™ PM10N portable SpO₂ patient monitoring system was ready to begin surveillance. After his birth, a gentle, accurate Nellcor™ OxySoft™ SpO₂ sensor was promptly placed on his right hand. The Nellcor™ system quickly detected continued low pulse rate and oxygen saturation, so Sean was resuscitated and moved into the NICU for respiratory support.

Nellcor™ pulse oximetry has been shown to post on average

12 seconds faster than Masimo SET™* technology.¹

Nellcor™ pulse oximetry sensors were

60% more accurate than Masimo SET™* sensors at detecting pulse rate when subjected to motion.^{2,3,†}

†Internal head-to-head studies evaluating adults in motion. Calculations based on a Medtronic analysis of the data in the referenced studies.

1. Khoury R, Klinger G, Shir Y, Osovsky M, Bromiker R. Monitoring oxygen saturation and heart rate during neonatal transition: comparison between two different pulse oximeters and electrocardiography. *J Perinatol.* 2021 Apr;41(4):885-890. doi: 10.1038/s41372-020-00881-y. Epub 2020 Nov 30. PMID: 33250516. Study objective: Compare efficacy and reliability of two pulse-oximeters (POx) (Masimo Radical-7 and Nellcor™ Oximax Bedside). Study design: Prospective observational monocentric comparative clinical study. 60 newborns included in total. Funding & conflict of interest: Authors declare no conflict of interest, Nellcor™ and Masimo provided sensors free of charge.
2. Addison PS, Mannheimer PD, Ochs J. Pulse rate performance of two pulse oximeters during challenging monitoring conditions. 2013 [White Paper].
3. Batchelder K, Sethi R, Eng B, Pinto YJ. Pulse rate performance of two pulse oximeters in the NICU. 2015 [White Paper].

Learn more about our performance with

Speed to post

Accuracy

Motion



We take our

patients

pulse

performance

portfolio

partnership

to heart

Patients are at the heart of all we do

Meet Sean

Meet Michael

Meet Ann

Meet Michael, a 56-year-old male recovering from post cardiac surgery involving an aortic valve replacement and ascending aortic aneurysm repair. Michael was admitted to the cardiac ICU post surgery. He was intubated and on a ventilator with a core temperature of 33°C. His hands were very cold indicating low perfusion to the fingers. A Nellcor™ OxySoft™ SpO₂ sensor was placed on one of his fingers to monitor oxygenation. His care team chose the Nellcor™ OxySoft™ SpO₂ sensor with its brighter LED and thoughtful cord placement to overcome the limitations that impact those readings.¹

Nellcor™
technology
missed

74%

less hypoxemic events
in subjects with dark skin and low
perfusion than Masimo technology.^{2,†}

†Oxygen saturation accuracy can be affected by certain environmental, equipment, and patient physiologic conditions that influence readings of SpO₂. Please consult the instructions for use and manual for full safety information.

1. RE00368468, RevB - Expanded Claims Bench Test Report.

2. Gudelas MK, Lipnick M, Hendrickson C, et al. Low Perfusion and Missed Diagnosis of Hypoxemia by Pulse Oximetry in Darkly Pigmented Skin: A Prospective Study. *Anesth Analg*. 2024;138(3):552-561. doi:10.1213/ANE.0000000000006755.

Learn more about our performance with

Low perfusion



We take our

patients

pulse

performance

portfolio

partnership

to heart

Patients are at the heart of all we do

Meet Sean

Meet Michael

Meet Ann

Meet Ann, a 62-year-old female with COPD. She came into the ER complaining of shortness of breath and frequent coughing. She had a low-grade temperature and showed signs of dyspnea. Chest X-rays showed flattened diaphragm, bronchiectasis, and emphysema. Ann had been taking Ipratropium and Advair for COPD and lisinopril for hypertension. Since Ann had been on oral and inhaled steroids for some time, her skin became increasingly fragile over the years, prone to bruising and tears at the slightest touch. While monitoring her pulse oximetry, her care team decided to use a Nellcor™ OxySoft™ SpO₂ sensor with silicone adhesive to help protect her fragile skin.

Nellcor™ OxySoft™ sensors are the first with a gentle, silicone adhesive that


removes 75% less skin cells

than Masimo LNCS Neo™ sensors^{1,2} – so you can remove without pulling fragile skin.

1. Based on internal report MDT20006OXYVMT, Rev 4, SpO₂ Accuracy Validation of OxySoft during motion and nonmotion. April 2021.
2. Based on internal report RE00357465, revA, Marketing validation report from a blinded hands-on evaluation conducted with 12 clinicians (RNs, RTs). April 2021.

Learn more about our performance with

[Skin integrity](#)



We take our

patients

pulse

performance

portfolio


partnership

to heart

We start with the heart

Perfecting our technology for the **most challenging 5% of patients**

Our unique approach uses the **pulse as our foundation**. This results in moment-to-moment readings that help you keep up with changes in your patient's condition. The rigor we take in perfecting our technology for more challenging pulse oximetry situations - newborns, the elderly, and patients with darker skin tones - helps ensure that all patients wearing our sensors benefit from quick and reliable readings.



We take our

patients

pulse

performance

portfolio

partnership

to heart

We start with the heart

Key differentiator #1
Identify the pulse


Key differentiator #2
Eliminate noise

Key differentiator #3
Calibrate the signal

Bring it together

First, we **identify the pulse** with our pattern matching neural network (AI)

Our advanced algorithms have been refined using decades of data from multiple patient profiles and characteristics. We lock onto the pulse so we can track moment-to-moment changes in SpO₂ or Pulse Rate.



We take our

patients

pulse

performance

portfolio

partnership

to heart

We start with the heart

Key differentiator #1
Identify the pulse


Key differentiator #2
Eliminate noise

Key differentiator #3
Calibrate the signal

Bring it together

Next, we **eliminate noise** by scrubbing the signal

Noise can be generated from motion and low perfusion. We use **Cardiac Gated Averaging (CGA)** to “scrub” the signal to clearly read the pulse through the noise.



We take our

patients

pulse

performance

portfolio

partnership

to heart

We start with the heart

Key differentiator #1
Identify the pulse



Key differentiator #2
Eliminate noise

Key differentiator #3
Calibrate the signal

Bring it together

Then, we **calibrate the signal** with OxiMax™ technology

OxiMax™ technology is part of our OxiMax™ system – where we embed a digital memory chip in every sensor matched with individualized calibration curves to help enable higher accuracy and improved performance.



Digital memory chip embedded in each sensor

We take our

patients

pulse

performance

portfolio

partnership

to heart

We start with the heart

Key differentiator #1
Identify the pulse


Key differentiator #2
Eliminate noise

Key differentiator #3
Calibrate the signal

Bring it together

Lastly, we **bring it together** with the Nellcor™ OxiMax™ system

Our entire OxiMax™ system – from the boards in monitors to cables and sensors – are designed to work together to accurately capture, transmit, and convert each heartbeat into an SpO₂ and PR reading. The entire system is needed to deliver accurate moment-to-moment readings that are confidently grounded in arterial signals versus surrounding tissue.



We take our

patients

pulse

performance

portfolio

partnership

to heart

Designed with your everyday challenges in mind

Speed to post

Accuracy

Motion

Low perfusion

Alarms & workflow

Skin integrity

Nellcor™ pulse oximetry is **designed to integrate into your daily workflow.**

Our technologies help solve the sensor performance and alarm management issues that can distract you from patient care.

With Nellcor™ pulse oximetry, you can:

- Get the readings you need
- Believe the numbers
- Trust your alarms



We take our

patients

pulse

performance

portfolio

partnership

to heart

In action

See how our speed to post performance impacts patient care.

Meet Sean

Designed with your everyday challenges in mind

Speed to post

Accuracy

Motion

Low perfusion


Alarms & workflow

Skin integrity

Seconds count, especially in neonatal care decisions.¹
Don't lose them waiting for an accurate vital signs reading.


Nellcor™ pulse oximetry has been shown to post on average **12 seconds faster** than Masimo SET™* technology.¹

1. Wyckoff MH, Aziz K, Escobedo MB, Kapadia VS, Kattwinkel J, Perlman JM, Simon WM, Weiner GM, Zaichkin JG. Part 13: Neonatal Resuscitation: 2015 American Heart Association Guidelines Update for Cardiopulmonary Resuscitation and Emergency Cardiovascular Care. *Circulation*. 2015 Nov 3;132(18 Suppl 2):S543-60. doi: 10.1161/CIR.000000000000267. PMID: 26473001..
2. Khoury R, Klinger G, Shir Y, Osovsky M, Bromiker R. Monitoring oxygen saturation and heart rate during neonatal transition: comparison between two different pulse oximeters and electrocardiography. *J Perinatol*. 2021 Apr;41(4):885-890. doi: 10.1038/s41372-020-00881-y. Epub 2020 Nov 30. PMID: 33250516. Study objective: Compare efficacy and reliability of two pulse-oximeters (POx) (Masimo Radical-7 and Nellcor™ Oximax Bedside). Study design: Prospective observational monocentric comparative clinical study. 60 newborns included in total. Funding & conflict of interest: Authors declare no conflict of interest, Nellcor™ and Masimo provided sensors free of charge.
3. Medtronic Reports: BRPMS Essential Requirements Matrix (Doc no. 10099165); Clinical, Motion, Connery OEM Module (Doc no. 10096395); and Clinical Report-COVMOPR0384, Motion, LAMP-C (internal Doc no. 10099560).
4. Addison PS, Mannheimer PD, Ochs J (2013). *Pulse rate performance of two pulse oximeters during challenging monitoring conditions* [White Paper]. Medtronic.
5. Saraswat A, Simionato L, Dawson JA, Thio M, Kamlin CO, Owen L, Schmölder G, Davis PG. Determining the best method of Nellcor pulse oximeter sensor application in neonates. *Acta Paediatr*. 2012 May;101(5):484-7. doi: 10.1111/j.1651-2227.2011.02571.x. Epub 2012 Jan 23. PMID: 22181562.



Nellcor™
OxiMax™
Pulse Oximetry

vs.



Masimo
Radical-7™*
Pulse CO-Oximeter

100%

Stable signal obtained²
(% of patients)

92%

15_{sec}

Average time to stable signal²

27_{sec}

17_{sec}

Method 1³⁻⁵
(sensor to oximeter to patient)

23_{sec}

14_{sec}

Method 2³⁻⁵
(sensor to oximeter to investigator then patient)

18_{sec}

11_{sec}

Method 3³⁻⁵
(sensor to patient then oximeter)

13_{sec}

We take our

patients

pulse

performance

portfolio

partnership

to heart

In action

See how our accuracy performance impacts patient care.

Meet Sean

Designed with your everyday challenges in mind

Speed to post

Accuracy

Motion

Low perfusion

Alarms & workflow

Skin integrity


Inaccurate pulse rate readings may guide clinicians to inappropriate or unnecessary interventions.¹ Nellcor™ pulse oximetry has demonstrated **accuracy at saturation rates as low as 60%**.^{2,†} This is especially important for neonates whose saturation rates can be as low as 66% in the first minutes of life,^{3,4} and for patients with conditions causing low perfusion.

± 2% accuracy at 70% SpO₂^{2,†}

± 3% accuracy at 60% SpO₂^{2,†}

†Oxygen saturation accuracy can be affected by certain environmental, equipment, and patient physiologic conditions (as discussed in the operator's manual for the monitor) that influence readings of SpO₂. Please consult the IFU and manual for full safety information.
‡Range Applicability: Ranges apply to Nellcor™ pulse oximetry OXIMAX, MAX-A, MAX-AL, MAX-N, MAX-I, MAX-P sensors; see sensor IFUs for complete information.

- Wyckoff MH, Aziz K, Escobedo MB, et al. Part 13: neonatal resuscitation: 2015 American Heart Association guidelines update for cardiopulmonary resuscitation and emergency cardiovascular care. *Circulation*. 2015;132(suppl 2):S543-S560.
- PMB10N (LP100):10143258 - Clinical Report, OxiCable (LP100 embedded) Abbreviated Sensor Line including the validation at 70%-100% saturation range with MaxA, MaxFAST, SCA, OxicleqA, DS100A, DYS-E-RE00105732, SpO₂ Accuracy Clinical Report of USB Pulse Oximetry Monitor Interface Cable with FlexMax Reusable SpO₂ Sensor via Reference CO-OximetryRE00085099 Verification Test Report, Sensor Accuracy, Oxicable including simulated patient data -NELL1: 10011350 (N600x) - Clinical Summary Report for N600x including the following validation studies: 1) N600x performance with MaxA, MaxFast, SC-A, DS-100A, OxiCliq A, D-YSE at 70-100% saturation range and 2) N600x comparison to N595/MaxA. -10028895 - MP100_O6 (Nell 1 equivalent) - Clinical Summary Report including validation of accuracy of MP100-O6 with Max AL, DS100A, D-YSE, OxiCliq, SC-A and MaxFast.
- Rabi Y, Dawson JA. Oxygen therapy and oximetry in the delivery room. *Semin Fetal Neonatal Med*. 2013 Dec;18(6):330-5. doi: 10.1016/j.siny.2013.08.007. Epub 2013 Sep 10. PMID: 24035476.
- Dawson JA, Kamlin CO, Vento M, et al. Defining the reference range for oxygen saturation for infants after birth. *Pediatrics*. 2010;125(6):e1340-e1347. doi: 10.1542/peds.2009-1510.
- Khoury R, Klinger G, Shir Y, Osovsky M, Bromiker R. Monitoring oxygen saturation and heart rate during neonatal transition. comparison between two different pulse oximeters and electrocardiography. *J Perinatol*. 2021 Apr;41(4):885-890. doi: 10.1038/s41372-020-00881-y. Epub 2020 Nov 30. PMID: 33250516.




Nellcor™
OxiMax™
Pulse Oximetry

r=0.894
STRONGLY
correlated

0%

0%

vs.



Masimo
Radical-7™**
Pulse CO-Oximeter

Heart rate compared to ECG⁵
(correlation coefficient)

r=0.235
WEAKLY
correlated

Mismatch ≥ 40 bpm compared to ECG⁵
(% of patients)

31%

False bradycardia⁵
(% of patients)

35%

We take our

patients

pulse

performance

portfolio

partnership

to heart

In action

See how our motion performance impacts patient care.

[Meet Sean](#)

Designed with your everyday challenges in mind

Speed to post

Accuracy

Motion

Low perfusion

Alarms & workflow

Skin integrity

You need accurate readings, even when patients are simply getting up, eating, or moving. Nellcor™ pulse oximetry sensors with OxiMax™ technology was the **first motion-tolerant technology**¹ to comply with International Standard ISO 80601-2-61.

Nellcor™ pulse oximetry sensors were **60%** **more accurate** than Masimo SET™* sensors at detecting pulse rate when subjected to motion.^{2,3,†}

Detects hypoxia with **95%+** **specificity** for patients in motion.⁴



†Internal head-to-head studies evaluating adults in motion. Calculations based on a Medtronic analysis of the data in the referenced studies.
1. Medtronic Reports: BRPMS Essential Requirements Matrix (Doc no. 10099165); Clinical, Motion, Connery OEM Module (Doc no. 10096395); and Clinical Report-COVMOPR0384, Motion, LAMP-C (Doc no. 10099560).
2. Addison PS, Mannheim PD, Ochs J. White paper: Pulse rate performance of two pulse oximeters during challenging monitoring conditions. 2013.
3. Batchelder K, Sethi R, Eng B, Pinto YJ. White paper: Pulse rate performance of two pulse oximeters in the NICU. 2015.
4. Louie A, Feiner JR, Bickler PE, Rhodes L, Bernstein M, Lucero J. Four Types of Pulse Oximeters Accurately Detect Hypoxia during Low Perfusion and Motion. *Anesthesiology*. 2018 Mar;128(3):520-530. doi: 10.1097/ALN.0000000000002002. PMID: 29200008.

We take our

patients

pulse

performance

portfolio

partnership

to heart

In action

See how our low perfusion performance impacts patient care.

Meet Michael

Designed with your everyday challenges in mind

Speed to post

Accuracy

Motion

Low perfusion

Alarms & workflow

Skin integrity

When low perfusion limits a peripheral reading, use **Nellcor™ SpO₂ forehead sensor** with OxiMax™ technology. It allows you to **monitor even your most challenging patients**,¹ including those with intense vasoconstriction, hypothermia, low cardiac index, septic shock, and severe peripheral vascular diseases.

Reflects changes to SpO₂ values **1-2 minutes earlier** than digit sensors for patients with weak pulses.¹

Nellcor™ technology **missed 74%** less hypoxemic events in subjects with dark skin and low perfusion than Masimo technology.^{2,†}



†Oxygen saturation accuracy can be affected by certain environmental, equipment, and patient physiologic conditions that influence readings of SpO₂. Please consult the instructions for use and manual for full safety information.
1. Bebout DE, Mannheimer PD, Wun C-C. Site-dependent differences in the time to detect changes in saturation during low perfusion. *Crit Care Med.* 2001;29(12):A115. Study objective: To test the hypothesis that during low perfusion, forehead sensors will detect saturation changes substantially sooner than sensors placed on fingers. Study design: 10 healthy adults.
2. Gudelunas MK, Lipnick M, Hendrickson C, et al. Low Perfusion and Missed Diagnosis of Hypoxemia by Pulse Oximetry in Darkly Pigmented Skin: A Prospective Study. *Anesth Analg.* 2024;138(3):552-561. doi:10.1213/ANE.0000000000006755.

We take our

patients

pulse

performance

portfolio

partnership

to heart

Designed with your everyday challenges in mind

Speed to post

Accuracy

Motion

Low perfusion

Alarms & workflow

Skin integrity

Pulse oximetry is the most common alarm in the hospital.¹ We know alarm fatigue is real. And it can negatively impact your workflow, your patient's experience, and your ability to provide the best possible care.


85-99% of alarms **don't require clinical intervention.**^{2,3}

Nellcor™ SatSeconds alarm management **may reduce nuisance alarms** by up to **40%**^{4,5}

1. Graham KC, Cvach M. Monitor alarm fatigue: standardizing use of physiological monitors and decreasing nuisance alarms. *Am J Crit Care.* 2010;19(1):28-35.
 2. Ver Hage A. Alarm fatigue can endanger patients. *OR Today.* May 1, 2015. Accessed 8/23/2019.
 3. The Joint Commission. Medical device alarm safety in hospitals. Sentinel Event Alert. April 8, 2013; issue 50.
 4. Brostowicz, Heather & Rais-Bahrami, K.. (2010). Oxygen saturation monitoring in the neonatal intensive care unit (NICU): Evaluation of a new alarm management. *Journal of Neonatal-Perinatal Medicine.* 3. 201-205. 10.3233/NPM-2010-0116.
 5. Stefanescu BM, O'Shea TM, Haury F, Carlo WA, Slaughter JC. Improved Filtering of Pulse Oximeter Monitoring Alarms in the Neonatal ICU: Bedside Significance. *Respir Care.* 2016 Jan;61(1):85-9. doi: 10.4187/respcare.04177. Epub 2015 Oct 27. PMID: 26508772.

How it works

The proprietary **Nellcor™ SatSeconds** algorithm is based off of decades of data. It works on a basic formula – generating alarms based on the severity of a drop in SpO₂ and the length of the drop, and filters out a very short dip in SpO₂. If SpO₂ drops below the alarm threshold for a longer period, the alarm will sound.



We take our

patients

pulse

performance

portfolio

partnership

to heart

In action

See how our skin integrity performance impacts patient care.

Meet Ann

Designed with your everyday challenges in mind

Speed to post

Accuracy

Motion

Low perfusion

Alarms & workflow

Skin integrity

You want to protect the fragile skin of your tiniest patients. But oftentimes, necessary sensors can damage that delicate skin – leading to unnecessary stress and inflammation. The **Nellcor™ OxySoft™ SpO₂ sensor** is the first pulse oximetry sensor to use a silicone adhesive to protect fragile skin and improve repositionability.

Nellcor™ OxySoft™ sensors **remove**

75%

less skin cells than Masimo LNCS Neo™ sensors – so you can remove without pulling fragile skin.^{1,2}

Nellcor™ OxySoft™ sensors don't stick to themselves and **stay adhered**

2x

better than Masimo LNCS Neo™ sensors.³

Nellcor™ OxySoft™ sensors **maintain**

85%

adhesiveness after 18 repositions⁴ compared to Masimo LNCS Neo™ sensors which lose 50% of adhesiveness after only 2 repositions.⁵

1. Based on internal report MDT20006OXYVMT, Rev 4, SpO2 Accuracy Validation of OxySoft during motion and nonmotion. April 2021.
2. Based on internal report RE00357465, revA, Marketing validation report from a blinded hands-on evaluation conducted with 12 clinicians (RNs, RTs). April 2021.
3. Based on research report CyberDERM S21-16 post market adhesive comparison of pulse oximeter sensors on adults with fragile skin commissioned by Medtronic. April 2022.
4. Based on internal report RE00357465, RevA - Marketing Validation Report from a blinded hands-on evaluation conducted with 12 clinicians (RNs, RTs).
5. Based on research report CyberDERM S21-16 post market adhesive comparison of pulse oximeter sensors on adults with fragile skin commissioned by Medtronic. April 2022.



We take our

patients

pulse

performance

portfolio

partnership

to heart

Solutions for every area of care

MPM solutions

Monitors

Sensors

From our thoughtfully designed monitors and sensors to meet all patient types and situations, to our partnerships with Multiple Parameter Monitor (MPM) providers, we can help you ensure that your hospital has a **full suite of pulse oximetry products** to fit within your hospital's current workflow.



We take our

patients

pulse

performance

portfolio

partnership

to heart

Solutions for every area of care

MPM solutions

Monitors

Sensors

We value collaboration

Our strategic partnerships offer innovative solutions to our customers on their **multiple parameter monitors (MPM)**.

Dräger

Hill-Rom
WelchAllyn

PHILIPS

**FUKUDA
DENSHI**

mindray


SPACE LABS
HEALTHCARE

 **GE Healthcare**

 **NIHON KOHDEN**

We take our

patients

pulse

performance

portfolio

partnership

to heart

Solutions for every area of care

MPM solutions

Monitors

Sensors

Bedside and portable monitors to meet your needs

Nellcor™ pulse oximetry monitors read every single heartbeat. And by posting timely,¹ sensitive answers and reducing clinically insignificant alarms,^{2,3} it puts you closer to each patient.



**RespArray™
Patient Monitor**



**Nellcor™ Bedside
Respiratory Patient
Monitoring System,
PM1000N**



**Nellcor™ Portable SpO₂
Patient Monitoring
System, PM10N**



**Nellcor™ Bedside SpO₂
Patient Monitoring System**



**Nellcor™ Bedside SpO₂ Patient
Monitoring System PM100N**

Nellcor™ monitors

1. Khoury R, Klinger G, Shir Y, Osovsky M, Bromiker R. Monitoring oxygen saturation and heart rate during neonatal transition. comparison between two different pulse oximeters and electrocardiography. *J Perinatol.* 2021 Apr;41(4):885-890. doi: 10.1038/s41372-020-00881-y. Epub 2020 Nov 30. PMID: 33250516.
2. Brostowicz HM, Khodayar Rais-Bahrami K. Oxygen saturation monitoring in the neonatal intensive care unit: evaluation of a new alarm management. Presented at: American Academy of Pediatrics National Conference & Exhibition; October 17-20, 2009; Washington D.C.
3. Stefanescu BM, O'Shea TM, Haury F, Carlo WA, Slaughter JC. Improved Filtering of Pulse Oximeter Monitoring Alarms in the Neonatal ICU: Bedside Significance. *Respir Care.* 2016 Jan;61(1):85-9. doi: 10.4187/respcare.04177. Epub 2015 Oct 27. PMID: 26508772.

We take our

patients

pulse

performance

portfolio

partnership

to heart

Solutions for every area of care

MPM solutions

Monitors

Sensors

The right sensor for each patient

The **Nellcor™ pulse oximetry sensor portfolio** offers a range of options so you can choose the right sensor to meet the needs of the patient.



Silicone adhesive

Nellcor™ OxySoft™ SpO ₂ sensor	OXYSOFTN
---	----------

Adhesive

Nellcor™ SpO ₂ adhesive sensor, adult	MAXA/MAXAL
--	------------

Nellcor™ SpO ₂ adhesive sensor, pediatric	MAXP
--	------

Nellcor™ SpO ₂ adhesive sensor, infant	MAXI
---	------

Nellcor™ SpO ₂ adhesive sensor, neonatal/adult	MAXN
---	------

Nellcor™ SpO ₂ adhesive sensor, adult nasal	MAXR
--	------

Non-adhesive

Nellcor™ SpO ₂ nonadhesive sensor, adult	SC-A
---	------

Nellcor™ SpO ₂ nonadhesive sensor, neonatal	SC-NEO
--	--------

Nellcor™ SpO ₂ nonadhesive sensor, preterm infant	SC-PR
--	-------



Specialty: forehead

Nellcor™ SpO ₂ forehead sensor	MAXFAST
---	---------

Specialty: respiration rate

Nellcor™ respiration rate sensor	10068119
----------------------------------	----------

Reusable

Nellcor™ reusable SpO ₂ sensor	DS100A-1
---	----------

Nellcor™ reusable multisite SpO ₂ sensors, multisite	D-YS
---	------

Nellcor™ reusable multisite SpO ₂ sensors, ear clip	D-YSE
--	-------

Nellcor™ reusable multisite SpO ₂ sensors, pediatric clip	D-YSPD
--	--------

Nellcor™ flexible SpO ₂ sensor, large	FLEXMAX
--	---------

Nellcor™ flexible SpO ₂ sensor, small	FLEXMAX-P
--	-----------

2-piece reusable

Nellcor™ reusable SpO ₂ sensor	OXI-A/N
---	---------

Nellcor™ reusable SpO ₂ sensor	OXI-P/I
---	---------



Nellcor™ sensors

We take our

patients

pulse

performance

portfolio

partnership

to heart

Value beyond our products

Education and training

Clinical support

Total cost of ownership

Society partnerships

We provide you with more than just the right tools. We're there for you with the **education, service, and support** you need throughout our partnership with you. Our products and pricing are designed to help you reduce spending and know your true costs – without hidden fees that can disrupt your budget.



We take our

patients

pulse

performance

portfolio

partnership

to heart

Value beyond our products

Education and training

Clinical support

Total cost of ownership

Society partnerships

Your partner in education

We are committed to offering complimentary learning opportunities that meet your educational and training needs, especially in this era of high turnover. Our experienced Clinical Specialists offer **in-person and virtual training** courses that will broaden your clinical expertise and enhance your technical knowledge and skills.

[Explore virtual courses](#)

Explore the MedEd Learning Experience

A podcast series where experts share safe and effective use of therapies in patient monitoring and respiratory interventions.

[Subscribe to the series](#)

We take our

patients

pulse

performance

portfolio

partnership

to heart

Value beyond our products

Education and training

Clinical support

Total cost of ownership

Society partnerships

We're there for you

Experienced field-based clinicians, nurses, and respiratory therapists relate to your daily challenges while providing **training, device support, and best practices** for patient safety.

Online support services

- Product manuals
- Sensor application guides
- Hardware user guides
- FAQs and education links
- Add-on software

U.S.-based tech support

- Get the help you need by email or phone (**1-800-NELLCOR**)
- Connect with representatives cross-trained on all patient monitoring products
- Request replacement monitors

Online support services



We take our

patients

pulse

performance

portfolio

partnership

to heart

Value beyond our products

Education and training

Clinical support

Total cost of ownership

Society partnerships

Helping you lower your cost of ownership

We believe in **true partnership**, providing you quality products without hidden costs, including:

**No charge
for installation
or training**



**Standard
3-year monitor
warranty**



**Long lasting
cables supporting
8x more sensors^{1,†}**



Save with Nellcor™ pulse oximetry

[†]The analysis compared cable utilization at four different hospitals ranging in bedsize of 290-550 that switched from Nellcor to Masimo.

• Each hospital had at least 3 years of data from both Masimo and Nellcor, within the time period of 2012-2020.

• The main endpoint was to evaluate how many cables were needed per 100,000 sensors used, while controlling for the varying number of sensors used over the timeframe with the different technologies.

1. Calculations based on internal Medtronic data and analysis on file, October, 2021.

We take our

patients

pulse

performance

portfolio

partnership

to heart

Value beyond our products

Education and training

Clinical support

Total cost of ownership

Society partnerships

We're committed to patient safety research and education

We work with clinical and patient safety societies to **fund research** that puts patient and clinician safety first.



Over
\$7M

invested in patient safety organizations[†]

[†]Five-year period.