


Medtronic

Engineering the extraordinary

Permacol™ surgical implant
The biological choice for chest wall reconstruction

Apply the (bio)logic



The biological
choice for chest
wall reconstruction

Resistant and durable^{1-3,†}

Permacol™ surgical implant
The biological choice for chest wall reconstruction

Permacol™ is a well-proven biologic implant that has been used globally in more than 350,000 patients over 12 years.¹⁻⁷

However, in cases of chest wall resection and reconstruction caused by pulmonary disease, infection, trauma, oncologic treatment or tumours, there are still risks for patients.

What are the risks of thoracic surgery?

Even if a patient is undergoing minimally invasive thoracic surgery, there are several complications that could potentially occur including:



Infection, including pneumonia.⁸



Pain.⁸



Bleeding.⁸



Air leakage from lungs.⁸

Which patients are most at risk?

Each patient will be reviewed for thoracic surgery based on their own individual risks. However, some main risk factors include:⁹

- **Pre-existing lung diseases or cancer.**
- **Pre-existing heart diseases or conditions.**
- **Age.**
- **Male sex.**
- **Being overweight.**
- **Being a smoker.**

Permacol™ surgical implant (HDMI cross-linked acellular porcine dermis) has been associated with reliable clinical outcomes in contaminated or infected fields compared to non-crosslinked biologic meshes currently available in the market.^{1,2,4,5,†}

†Use of Permacol™ surgical implant in a contaminated or infected field may lead to a weakening or breakdown of the implant. Treat any existing or suspected infection according to accepted medical practice before implanting the device.



Care for the thoracic wall

Permacol™ surgical implant
The biological choice for chest wall reconstruction

Whether it's as a result of trauma, oncological treatment, infection or disease, thoracic surgery requires making insertions into the chest wall. Even in the case of minimally invasive surgery, some incisions will need to be made.

Typically, incisions of less than 5 cm do not require any form of reconstruction,¹⁰ although this might not be the case depending on the patient and location or specificity of surgery required. In the case of larger incisions, mesh may be required to support the soft tissue in the chest wall.

Meshes present the advantages of easy manipulation and handling and comply with the characteristics of ideal prosthetic material as determined by Le Roux and Sherma:¹⁰



- **rigidity** to abolish paradoxical movement
- **inertness** to allow in-growth of fibrous tissue and decrease the likelihood of infection
- **malleability** to fashion to the appropriate shape at the time of operation
- **radiolucency** to create an anatomic reference to do a better follow up and identify a possible local neoplastic relapse

Using synthetic meshes in thoracic surgery

Permacol™ surgical implant
The biological choice for chest wall reconstruction

Typically, synthetic mesh has been the material of choice for thoracic surgery because it's simple to use and usually well-tolerated.¹⁰

However, biological mesh has been shown to hold less risks:



† Use of Permacol™ surgical implant in a contaminated or infected field may lead to a weakening or breakdown of the implant. Treat any existing or suspected infection according to accepted medical practice before implanting the device.

Using bioprosthetic meshes in thoracic surgery

Permacol™ surgical implant
The biological choice for chest wall reconstruction

Bioprosthetic meshes have been shown to be beneficial in chest wall resection and reconstruction surgery.^{11-13,†}

Permacol™ surgical implants in a challenging environment:

Although bioprosthetic meshes can be more expensive and difficult to handle, they do have significant benefits when used in thoracic surgery, including:




The ability to revascularize and integrate quickly into native tissues.¹⁻³



Infection resistance.^{10,†}



Effective wound healing and long-term stability, even in large surgeries.^{1-3,13}



A crosslinked structure that provides strength, durability and flexibility thanks to their microfibrils made of chemically connected collagen molecules.¹⁻³



In short, choosing bioprosthetic meshes rather than synthetic meshes for chest wall reconstruction, particularly in large or contaminated cases, can improve patient outcomes.^{1,2,6,14}

† Use of Permacol™ surgical implant in a contaminated or infected field may lead to a weakening or breakdown of the implant. Treat any existing or suspected infection according to accepted medical practice before implanting the device.

Make the (bio)logical choice

Permacol™ surgical implant
The biological choice for chest wall reconstruction

Choose Permacol™ surgical implant

A Permacol™ surgical implant is a tough and flexible fibrous flat sheet of crosslinked acellular porcine dermal collagen and constituent elastin fibres.†

An alternative to synthetic meshes, Permacol™ is a well-proven biologic implant that has been used globally in a range of surgical procedures, including for hernia and abdominal wall repair.^{1,2}

Available in a range of sizes and with no rehydration or refrigeration required, Permacol™ surgical implant is ready to use when needed and will provide dimensional stability, whatever the size of the defect.

Thanks to the unique crosslinking technology, which is formed in a chemical process using HMDI (Hexamethylene Diisocyanate), Permacol™ surgical implant offers advantages over other collagen-based implants which are not crosslinked for hernia and abdominal wall repair.^{2,3,15}



There are three main advantages that Permacol™ provides over biosynthetic implants:




Biocompatibility

Studies have shown that Permacol™ has good fibroblastic and neovascular infiltration, excellent biocompatibility and resistance to degradation in potentially contaminated sites.^{1-3,‡}



Strength and durability

Permacol™ has shown greater tensile strength compared to non-crosslinked implants. While human cadaveric grafts and other non-crosslinked grafts may initially be successful, many lose tensile strength and have increased tissue laxity over time.¹⁶⁻²⁰



Dimensional stability

Optimal cross-linking and gentle processing methods means the Permacol™ surgical implant offers long-lasting dimensional stability. This ensures the integrity of the collagen graft through the body's wound healing process.¹⁻³

† See IFU_PT00099456 Permacol

‡ Use of Permacol™ surgical implant in a contaminated or infected field may lead to a weakening or breakdown of the implant. Treat any existing or suspected infection according to accepted medical practice before implanting the device.

To find out more about the benefits Permacol™ surgical implant can offer you and your patients, head to www.medtronic.com/covidien/en-gb/products/hernia-repair/permacol-surgical-implant.html



Order code	Description	Thickness	Dimensional area	Quantity
5220-100	Permacol™	1.0 mm	2 × 20 cm	1
5033-100	Permacol™	1.0 mm	3 × 3 cm	1
5416-100	Permacol™	1.0 mm	4 × 16 cm	1
5418-100	Permacol™	1.0 mm	4 × 18 cm	1
5000-100	Permacol™	1.0 mm	5 × 5 cm	1
5001-100	Permacol™	1.0 mm	5 × 10 cm	1
5616-100	Permacol™	1.0 mm	6 × 16 cm	1
5816-100	Permacol™	1.0 mm	8 × 16 cm	1
5110-100	Permacol™	1.0 mm	10 × 10 cm	1
5115-100	Permacol™	1.0 mm	10 × 15 cm	1
5210-100	Permacol™	1.0 mm	10 × 20 cm	1
5152-100	Permacol™	1.0 mm	15 × 20 cm	1
5120-100	Permacol™	1.0 mm	18 × 28 cm	1
5230-100	Permacol™	1.0 mm	20 × 30 cm	1
5033-150	Permacol™	1.5 mm	3 × 3 cm	1
5063-150	Permacol™	1.5 mm	3 × 6 cm	1
5000-150	Permacol™	1.5 mm	5 × 5 cm	1
5001-150	Permacol™	1.5 mm	5 × 10 cm	1
5110-150	Permacol™	1.5 mm	10 × 10 cm	1
5115-150	Permacol™	1.5 mm	10 × 15 cm	1
5152-150	Permacol™	1.5 mm	15 × 20 cm	1
5120-150	Permacol™	1.5 mm	18 × 28 cm	1
5230-150	Permacol™	1.5 mm	20 × 30 cm	1
5240-150	Permacol™	1.5 mm	20 × 40 cm	1
5250-150	Permacol™	1.5 mm	20 × 50 cm	1
5284-150	Permacol™	1.5 mm	28 × 40 cm	1

References

1. N. J. Smart, N. Bryan, J. A. Hunt. A scientific evidence for the efficacy of biologic implants for soft tissue reconstruction. First published: 08 November 2012; Special Issue: Biologic Meshes in Colorectal Surgery. Volume 14, Issue Supplement s3, pages 1-6, December 2012. | <https://doi.org/10.1111/codi.12042>.
2. Smart NJ, Bloor S. Durability of Biologic Implants for Use in Hernia Repair: A Review. *Surg Innov*, First Published December 4, 2011; 221 -229; <https://doi.org/10.1177/1553350611429027>.
3. De Castro Brás LE, Proffitt JL, Bloor S, Sibbons PD. Effect of crosslinking on the performance of a collagen-derived biomaterial as an implant for soft tissue repair: a rodent model. *Biomed Mater Res B Appl Biomater*. 2010; 95:239-249*. Published online 28 September 2010 in Wiley Online Library (wileyonlinelibrary.com). DOI: 10.1002/jbm.b.31704.
4. Chand B, Indeck M, Needleman B, Finnegan M, Van Sickle KR, Ystgaard B, Gossetti F, Pullan RD, Giordano P, McKinley A. A retrospective study evaluating the use of Permacol™ surgical implant in incisional and ventral hernia repair. *Int J Surg*. 2014;12(4):296-303.
5. Neil J. Smart, Morwena Marshall, Ian R. Daniels. Biological meshes: A review of their use in abdominal wall hernia repairs. *The Surgeon*, Available online March 2012;159-171
6. Hiles M. et al. Are Biologic Grafts Effective for Hernia Repair? A Systematic Review of the Literature. *Surgical Innovation / Vol. 16, No. 1, March 2009; 26-37* © 2009 Sage Publications 10.1177/1553350609331397.
7. Complaints & Vigilances Review - Permacol™ Surgical Implant and ENDURAGen™ Collagen Implant. CVR020j. February 2021
8. Video Assisted Thoracoscopic Surgery- (VATS) Procedure. University Hospitals Coventry and Warwickshire NHS Trust. Pge3.pdf. HIC/LFT/2407/20 February 2020.
9. Irons JF, and Martinez G. Complex, high-risk thoracic surgery—does risk always outweigh the benefit or can we manage it safely? *Video-Assisted Thoracic Surgery*. Vol 2, no.9,2017. <https://vats.amegroups.com/article/view/4059>.
10. Sanna S, Brandolini J, Pardolesi A, et al. Materials and techniques in chest wall reconstruction: a review. *J Vis Surg*. 2017;3:95. Published 2017 Jul 26. doi:10.21037/jovs.2017.06.10.
11. Barua A, Catton JA, Socci L, Raurell A, Malik M, Internullo E, Martin-Ucar AE. Initial experience with the use of biological implants for soft tissue and chest wall reconstruction in thoracic surgery. *Ann Thorac Surg*. 2012 Nov;94(5):1701-5. doi: 10.1016/j.athoracsur.2012.07.001. Epub 2012 Sep 7. PMID: 22959570.
12. Shankaran V, Weber DJ, Reed RL, 2nd, Luchette FA. A review of available prosthetics for ventral hernia repair. *Ann Surg*. Jan 2011;253(1):16-26. doi:10.1097/SLA.0b013e3181f9b6e6
13. D'Amico G, Manfredi R, Nita G, et al. Reconstruction of the Thoracic Wall With Biologic Mesh After Resection for Chest Wall Tumors: A Presentation of a Case Series and Original Technique. *Surgical Innovation*. 2018;25(1):28-36. doi:10.1177/1553350617745954.
14. Gaertner, W.G., et al. Experimental Evaluation of Four Biologic Prostheses for Ventral Hernia Repair. *J. Gastrointest Surg*. 2007.
15. Annor AH et al. Effect of enzymatic degradation on the mechanical properties of biological scaffold materials. *Surg Endosc* (2012) 26:2767-2778.
16. Zheng, F, et al. Host Response after reconstruction of abdominal wall defects with porcine dermal collagen in a rat model. *American Journal of Obstetrics and Gynecology*(2004)191,1961-70.
17. Jin, J, et al. Abdominal Hernia Repair with Bridging Acellular Dermal Matrix - An Expensive Hernia Sac. *American Hernia Society* March 2007.
18. Jin, J, et al. Use of Acellular Dermal Matrix for Complicated Ventral Hernia Repair: Does Technique Affect Outcomes? *Journal of the American College of Surgeons* Volume 205, Issue 5, November 2007, Pages 654-660.
19. Rodriguez, Eduardo, et al. Abdominal Wall Reconstruction following Severe Loss of Domain: The R Adams Cowley Shock Trauma Center Algorithm. *PRSJournal* September 2007 DIO:10.1097/01/prs0000270303.44219.76.
20. Newcomb, W.; Heniford, T, et al . Long-term Outcomes After Abdominal Wall Reconstruction With Acellular Human Dermis in Patients With Abdominal Injuries. 2007 AAST Session IV, Poster #10.

Important: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2023 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. EMEA-HR-2200050-permacol-in-thoracic-brochure-en-emea-9879772