

MEDTRONIC DEEP BRAIN STIMULATION (DBS) FOR EPILEPSY

SUMMARY OF CLASS I EVIDENCE FOR NEUROMODULATION THERAPIES

Therapies	MEDTRONIC DEEP BRAIN STIMULATION (DBS)	LIVANOVA VAGUS NERVE STIMULATION (VNS)
	Implanted brain stimulation system for epilepsy, Parkinson's disease, essential tremor, dystonia	Implanted autonomic nerve stimulator for epilepsy and Treatment-Resistant Depression

Results are from different studies and are shown for illustrative purposes only. Results may differ in a head-to-head study. Refer to product labeling for specific information regarding indications, contraindications, warnings, precautions and safety.

	DBS	VNS		
Regulatory Status for Epilepsy	CE Mark (2010) FDA Approval (2018)	CE Mark (1994) FDA Approval (1997)		
Indication	Bilateral stimulation of the anterior nucleus of the thalamus (ANT) as an adjunctive therapy for reducing the frequency of seizures in individuals 18 years of age or older diagnosed with epilepsy characterized by partial-onset seizures, with or without secondary generalization, that are refractory to antiepileptic medications. The Medtronic DBS System for Epilepsy has demonstrated safety and effectiveness for patients who average six or more seizures per month over the three most recent months prior to implant of the DBS system (with no more than 30 days between seizures). The Medtronic DBS System for Epilepsy has not been evaluated in patients with less frequent seizures. ¹	An adjunctive therapy in reducing the frequency of seizures in patients 4 years of age and older with partial onset seizures that are refractory to antiepileptic medications. ²		
Published Clinical Studies	SANTE Pivotal Study - RCT <ul style="list-style-type: none"> ■ Blinded Phase Results³ ■ Long-term Results^{4,5} 	E03 / E05 Studies - RCTs <ul style="list-style-type: none"> ■ Blinded Phase Results^{6,7} ■ Long-term Results⁹ 		
Stimulation / Target Procedure	Stereotactic placement of bilateral leads in the Anterior Nucleus of Thalamus (ANT)	Electrodes wrapped around cranial nerve accessed via the neck		
MRI	Conditional; full-body 1.5T Scan ⁹	Conditional: device set to 0mA; transmit RF body coil would require system explant ²		
Implant Location of the Neurostimulator	Pectoral or Abdominal	Pectoral		
Clinical Programming / Therapy Maintenance	Programming for optimal therapy benefit	Programming for optimal therapy benefit		
Blinded Phase Results - Median % Change	Active: -35.0% Control: -21.1% ¹⁰	<table border="1" style="width: 100%; text-align: center;"> <tr> <td>E03² High: -23% Low: -6%</td> <td>E05² High: -23% Low: -21%</td> </tr> </table>	E03² High: -23% Low: -6%	E05² High: -23% Low: -21%
E03² High: -23% Low: -6%	E05² High: -23% Low: -21%			
Blinded Phase Results - Responder Rate (% of subjects with ≥50% reduction in seizures)	Active: 29.6% Control: 25.9% ¹⁰	<table border="1" style="width: 100%; text-align: center;"> <tr> <td>E03² High: 30% Low: 14%</td> <td>E05² High: 23% Low: 16%</td> </tr> </table>	E03² High: 30% Low: 14%	E05² High: 23% Low: 16%
E03² High: 30% Low: 14%	E05² High: 23% Low: 16%			
Seizure Reduction 3 Years - Median % Change	57% Last Observation Carried Forward ¹⁰	44% Last Value Carried Forward ⁸		
Seizure Reduction 6 Years - Median % Change	70% Last Observation Carried Forward ¹⁰	Data Not Available		
Quality of Life Results Blinded Phase - Responder Rate (5 point or greater improvement on the QOLIE Assessment)	Active: 48.1% Control: 32.1% ¹⁰	Not Assessed in E03/E05		
Quality of Life Results: Long Term - Responder Rate (5 point or greater improvement on the QOLIE Assessment)	1 Year: 46.1% 2 Year: 38.8% 5 Year: 47.5% 7 Year: 43.3% ¹⁰	Improvement on various QOL metrics subjectively rated by physicians Provider Survey Data ¹¹		
Epilepsy Related Injuries : Blinded Phase	Active Group exhibited significantly fewer epilepsy related injuries (Active: 7.4% vs Control: 25.5%) ¹⁰	NA		

Healthcare Resource Utilization (Epilepsy and Device Related Visits Normalized to 1 Year)	Active Group had fewer hospital visits but difference not statistically significant (Active: 4.2 vs Control: 20.2) Blinded Phase Results ¹⁰	Reduction in seizure-related hospitalizations. RTI Health Solutions ¹²
Mood and Memory Outcomes	Subjective reports of depression and memory impairment during the blinded phase. Neuropsych scores are stable or improved over the long term (7 years). ⁵	No clear cognitive changes during the blinded phase (E05). Long-term results not reported. ¹³
Most Common Serious Device Related Adverse Event	Implant Site Infection (10.9%) 713 Device Years ¹⁰	Surgically Related (4.1%)* 591 Device Years in E03/E05 ²
Sudden Unexplained Death in Epilepsy (SUDEP) Rate (95% Confidence Interval)	2.5 deaths / 1000 Years (0.31–9.16) ¹⁰	4.1 deaths / 1000 Years (CI not available) ¹⁴

* Included infection, nerve paralysis, hypesthesia, facial paresis, left vocal cord paralysis, left facial paralysis, left hemidiaphragm paralysis, left recurrent laryngeal nerve injury, urinary retention, and low-grade fever.

REFERENCES

1. Medtronic. *Physician Indication Insert for the DBS indication*. M933213A002.
2. LivaNova. VNS Therapy® System Epilepsy Physician's Manual. April 2019. US Version. 76-0000-5400/3 (US).
3. Fisher R, Salanova V, Witt T, et al. SANTE Study Group. Electrical stimulation of the anterior nucleus of thalamus for treatment of refractory epilepsy. *Epilepsia* 2010;51:899-908.
4. Salanova V, Witt T, Worth R, et al. Long-term efficacy and safety of thalamic stimulation for drug-resistant partial epilepsy. *Neurology* 2015;84:1017-1025.
5. Tröster AI, Meador KJ, Irwin CP, Fisher RS; SANTE Study Group. Memory and mood outcomes after anterior thalamic stimulation for refractory partial epilepsy. *Seizure*. 2017 Feb;45:133-141.
6. The Vagus Nerve Stimulation Study Group. A randomized controlled trial of chronic vagus nerve stimulation for treatment of medically intractable seizures. *Neurology* 1995;45:224-230.
7. Handforth A, DeGiorgio CM, Schachter SC, et al. Vagus nerve stimulation therapy for partial-onset seizures: a randomized active-control trial. *Neurology* 1998;51:48-55.
8. Morris GL, Mueller WM, and The Vagus Nerve Stimulation Study Group E01-E05. Long-term treatment with vagus nerve stimulation in patients with refractory epilepsy. *Neurology* 1999;53:1731-1735.
9. Medtronic MRI Guidelines for deep brain stimulation systems. Medtronic Inc. 2015. M929535A043 REV.A
10. Medtronic. *Medtronic DBS Therapy for Epilepsy Clinical Summary*. 2018-10-15. M959283A001 Rev C.
11. Englot DJ, Hassnain KH, Rolston JD, Harward SC, Sinha SR, Haglund MM. Quality-of-life metrics with vagus nerve stimulation for epilepsy from provider survey data. *Epilepsy Behav*. 2017 Jan;66:4-9.
12. Purser MF, Mladsi DM, Beckman A, Barion F, Forsey J. Expected Budget Impact and Health Outcomes of Expanded Use of Vagus Nerve Stimulation Therapy for Drug-Resistant Epilepsy. *Adv Ther*. 2018 Oct;35(10):1686-1696.
13. Dodrill CB, Morris GL. Effects of Vagal Nerve Stimulation on Cognition and Quality of Life in Epilepsy. *Epilepsy Behav*. 2001 Feb;2(1):46-53.
14. Annegers JF, Coan SP, Hauser WA, Leestma J. Epilepsy, vagal nerve stimulation by the NCP system, all-cause mortality, and sudden, unexpected, unexplained death. *Epilepsia*. 2000 May;41(5):549-53.

Medtronic

Europe
 Medtronic International Trading Sarl
 Route du Molliiau 31
 Case postale
 1131 Tolochenaz
 Switzerland
 Tel: +41 (0) 21 802 70 00
 Fax: +41 (0) 21 802 79 00

Brief Statement:

See the device manual for detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events. If using an MRI SureScan® device, see the MRI SureScan® technical manual before performing an MRI. For further information, contact your local Medtronic representative and/or consult the Medtronic website at medtronic.eu.

UC202004637 EE © Medtronic 2023.
 All rights reserved.