

Medtronic

Engineering the extraordinary

Animal Health

Surgical devices, respiratory and monitoring

Product catalog



800.962.9888

Table of Contents

Wound closure and hernia mesh products	2
Skin, open, and laparoscopic stapling	25
Electrosurgery and vessel sealing	34
Smoke evacuation	50
Patient monitoring	52
Access and trocar instruments	54
Laparoscopic hand instruments	66



Wound closure and hernia mesh products



Family of wound closure products

Absorbable



Velosorb™ Fast Suture **Polysorb™ Suture** **Caprosyn™ Suture** **Biosyn™ Suture** **Maxon™ Suture**

Structure	Braided	Braided	Monofilament	Monofilament	Monofilament
Suture type	Short-Term Absorbable	Mid-Term Absorbable	Short-Term Absorbable	Mid-Term Absorbable	Long-Term Absorbable
Composition	Synthetic polyester comprised of glycolide and lactide	LACTOMER™ 9-1 Glycolide and Lactide	POLYGLYTON™ 6211 Glycolide, Caprolactone, Trimethylene Carbonate, and Lactide	GLYCOMER™ 631 Glycolide, Trimethylene Carbonate and Dioxanone	POLYGYCONATE Copolymer of glycolic acid and trimethylene carbonate
Coating	Glycolide and lactide copolymer and calcium stearate	Glycolide, Caprolactone and Calcium Stearoyl Lactylate	Uncoated	Uncoated	Uncoated
Color of material	Violet Undyed	Violet Undyed	Violet Undyed	Violet Undyed	Green Clear
Tensile strength retention in-vivo	7 days	3 Weeks	10 Days	3 Weeks	6 Weeks
Absorption profile	40-50 Days	56-70 Days	56 Days	90-110 Days	180 Days
Sizes	Sizes 6-0 to 1	Sizes 8-0 to 2	Sizes 6-0 to 1	Sizes 6-0 to 1	Sizes 7-0 to 1

Non-absorbable



Surgipro™ and Surgipro™ II Suture **Vasculfil™ Suture** **Ti-Cron™ / Surgidac™ Suture** **Monosof™ / Dermalon™ Suture** **Surgilon™ Suture**

Structure	Monofilament	Monofilament	Braided	Monofilament	Braided
Suture type	Non-Absorbable	Non-Absorbable	Non-Absorbable	Non-Absorbable	Non-Absorbable
Composition	Polypropylene	Polybutester	Polyester	Nylon	Nylon
Coating	Uncoated	Polytribolate	Ticron: Coated with silicon and uncoated Surgidac: Uncoated and Coated (Polybutylene Adipate)	Uncoated	Silicone
Color of material	Blue Clear	Blue	Blue or white (Ticron™) Green or white (Surgidac™)	Black/Clear (Monosof™) Blue (Dermalon™)	Black White
Tensile strength retention in-vivo	Permanent	Permanent	Permanent	Gradual Loss Over Time	Gradual Loss Over Time
Absorption profile	Permanent	Permanent	Permanent	Permanent	Permanent
Sizes	Sizes 10-0 to 2 (Surgipro™) Sizes 8-0 to 3-0 (Surgipro™ II)	Sizes 7-0 to 2-0	Sizes 5- 8/0 (Ticron™) Sizes 2-6-0 (Surgidac™)	Sizes 11-0 to 2 (Monosof™) Sizes 2-0 to 6-0 (Dermalon™)	Sizes 5-0 - 2

Family of wound closure products

Absorbable



	V-Loc™ 90 Device	V-Loc™ 180 Device
Structure	Barbed Monofilament	Barbed Monofilament
Suture Type	Mid-Term Absorbable	Long-Term Absorbable
Composition	GLYCOMER™ 631 Glycolide, dioxanone and trimethylene carbonate	POLYGYCONATE Copolymer of glycolic acid and trimethylene carbonate
Coating	Uncoated	Uncoated
Color of Material	Violet Undyed	Green Clear
Tensile Strength Retention In-Vivo	2 Weeks	3 Weeks
Absorption Profile	90-110 days	180 Days
Sizes	Sizes 4-0 to 2-0	Sizes 4-0 to 0

Natural (Absorbable)



	Chromic Gut	Plain Gut
Structure	Natural	Natural
Suture Type	Short/Mid-term Absorbable	Short/Mid-term Absorbable
Composition	Purified connective tissue (mostly collagen) derived from the serosal layer of beef (bovine) intestines	Purified connective tissue (mostly collagen) derived from the serosal layer of beef (bovine) intestines
Coating	Uncoated	Uncoated
Color of Material	Undyed	Undyed
Tensile Strength Retention In-Vivo	21-28 days	7-10 days
Absorption Profile	90 days	70 days
Sizes	7-0 to 2	6-0 to 1

Non-absorbable



	Novafil™ Suture	Flexon™ Temporary Cardiac Pacing Lead	Steel Suture	V-Loc™ PBT Device
Structure	Monofilament	Multifilament	Monofilament	Barbed Monofilament
Suture Type	Non-Absorbable	Non-Absorbable	Non-Absorbable	Non-Absorbable
Composition	Polybutester	Stainless Steel	Stainless Steel	Copolymer of butylene terephthalate and polytetramethylene ether glycol
Coating	Uncoated	FEP Polymer or PTFE Polymer	Uncoated	Uncoated
Color of Material	Blue Clear	Clear Orange White	Steel	Blue
Tensile Strength Retention In-Vivo	Permanent	Permanent	Permanent	Permanent
Absorption Profile	Permanent	Permanent	Permanent	Permanent
Sizes	Sizes 7-0 to 2	Sizes 3-0 to 0	Sizes 3-0 to 7	Sizes 3-0 to 1

Natural (Non-absorbable)




	Sofsilk™ Suture
Structure	Natural
Suture Type	Non-Absorbable
Composition	Silk
Coating	Silicone/Wax
Color of Material	Black White
Tensile Strength Retention In-Vivo	Gradual Loss Over Time
Absorption Profile	Permanent
Sizes	Sizes 5-0 and 2-0 through 8-0

Sutures

Needle chart

Reverse cutting
3/8 circle



Caprosyn™ Monofilament Absorbable Suture


Provides short-term tensile strength with the benefits of rapid absorption (56 days).

Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
C-13	reverse cutting, 3/8 Circle, 19 mm	3-0	30"	Undyed	SC692	36
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	30"	Undyed	SC691	36
P-13	premium reverse cutting, 3/8 Circle, 13 mm	4-0	18"	Violet	SC5590G	12
P-13	premium reverse cutting, 3/8 Circle, 13 mm	5-0	18"	Violet	SC5587G	12
P-12	premium reverse cutting, 3/8 Circle, 19 mm	4-0	18"	Violet	SC5627G	12
CV-23	taper point, 1/2 Circle, 17 mm	3-0	30"	Undyed	UC204	36
CV-23	taper point, 1/2 Circle, 17 mm	4-0	30"	Undyed	UC203	36
CV-23	taper point, 1/2 Circle, 17 mm	5-0	30"	Undyed	UC202	36


Biosyn™ Monofilament Absorbable Suture

Provides mid-term tensile strength.

Reverse Cutting
1/2 circle








Taper Point
1/2 circle



Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
C-13	reverse cutting, 3/8 Circle, 19 mm	3-0	30"	Undyed	SM822	36
C-13	reverse cutting, 3/8 Circle, 19 mm	3-0	30"	Violet	SM922	36
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	30"	Undyed	SM691	36
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	30"	Violet	SM791	36
C-13	reverse cutting, 3/8 Circle, 19 mm	5-0	30"	Undyed	SM690	36
C-14	reverse cutting, 3/8 Circle, 24 mm	2-0	30"	Undyed	SM823	36
C-14	reverse cutting, 3/8 Circle, 24 mm	2-0	30"	Violet	SM923	36
C-14	reverse cutting, 3/8 Circle, 24 mm	3-0	30"	Undyed	SM693	36
C-23	reverse cutting, 1/2 circle, 24 mm	3-0	30"	Violet	SM632	36
GS-10	reverse cutting, 1/2 circle, 26 mm	2-0	30"	Violet	CM879	36
GS-11	reverse cutting, 1/2 circle, 37 mm	0	30"	Violet	CM817	36
GS-11	reverse cutting, 1/2 circle, 37 mm	2-0	30"	Violet	CM816	36
CV-23	taper point, 1/2 Circle, 17 mm	3-0	30"	Violet	UM204	36
CV-23	taper point, 1/2 Circle, 17 mm	4-0	30"	Violet	UM203	36
CV-23	taper point, 1/2 Circle, 17 mm	4-0	30"	Violet	UM403	36
CV-23	taper point, 1/2 Circle, 17 mm	5-0	30"	Violet	UM202	36
V-20	taper point, 1/2 circle, 26 mm	2-0	30"	Violet	GM123	36
V-20	taper point, 1/2 circle, 26 mm	3-0	30"	Violet	GM122	36
V-20	taper point, 1/2 circle, 26 mm	4-0	30"	Violet	GM121	36
GS-22	taper point, 1/2 circle, 27 mm	0	30"	Violet	CM884	36
GS-22	taper point, 1/2 circle, 27 mm	2-0	30"	Violet	CM883	36
GS-22	taper point, 1/2 circle, 27 mm	3-0	30"	Violet	CM882	36

Sutures

Needle chart

 Reverse cutting
3/8 circle

 Reverse Cutting
1/2 circle

 Taper Point
1/2 circle


Biosyn™ Monofilament Absorbable Suture continued

Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
GS-22	taper point, 1/2 circle, 30 mm	0	36"	Violet	GM327	36
GS-21	taper point, 1/2 Circle, 37 mm	0	30"	Violet	CM812	36
GS-21	taper point, 1/2 Circle, 37 mm	0	36"	Violet	CM924	36
GS-21	taper point, 1/2 Circle, 37 mm	1	30"	Violet	CM813	36
GS-21	taper point, 1/2 Circle, 37 mm	2-0	30"	Violet	CM811	36
GS-21	taper point, 1/2 Circle, 37 mm	2-0	36"	Violet	CM923	36
GS-21	taper point, 1/2 Circle, 37 mm	3-0	30"	Violet	CM810	36
GS-25	taper point, 1/2 Circle, 48 mm	1	36"	Violet	CM905	36

Maxon™ Monofilament Absorbable Suture

Provides long-term tensile strength and security for about six weeks.

Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
C-13	reverse cutting, 3/8 Circle, 19 mm	3-0	30"	Green	8886663941	36
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	30"	Clear	8886661231	36
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	30"	Green	8886663931	36
C-13	reverse cutting, 3/8 Circle, 19 mm	5-0	18"	Clear	8886660821	36
C-14	reverse cutting, 3/8 Circle, 24 mm	2-0	30"	Clear	8886661851	36
C-14	reverse cutting, 3/8 Circle, 24 mm	2-0	30"	Green	8886664151	36
C-14	reverse cutting, 3/8 Circle, 24 mm	3-0	30"	Clear	8886661841	36
C-14	reverse cutting, 3/8 Circle, 24 mm	3-0	30"	Green	8886664141	36
GS-10	reverse cutting, 1/2 circle, 26 mm	2-0	30"	Green	8886674551	36
GS-11	reverse cutting, 1/2 circle, 37 mm	0	30"	Green	8886675561	36
GS-11	reverse cutting, 1/2 circle, 37 mm	1	30"	Green	8886675571	36
GS-11	reverse cutting, 1/2 circle, 37 mm	2-0	30"	Green	8886675551	36
CV-23	taper point, 1/2 Circle, 17 mm	4-0	30"	Green	8886620731	36
CV-23	taper point, 1/2 Circle, 17 mm	5-0	30"	Green	8886620721	36
CV-24	taper point, 1/2 circle, 20 mm	3-0	30"	Green	8886621741	36
CV-24	taper point, 1/2 circle, 20 mm	4-0	30"	Green	8886621731	36
V-20	taper point, 1/2 circle, 26 mm	2-0	30"	Green	8886623351	36
V-20	taper point, 1/2 circle, 26 mm	3-0	30"	Green	8886623341	36
V-20	taper point, 1/2 circle, 26 mm	3-0	36"	Green	8886623541	36
V-20	taper point, 1/2 circle, 26 mm	4-0	42"	Green	8886614331	36
V-20	taper point, 1/2 circle, 26 mm	4-0	30"	Green	8886623331	36
GS-22	taper point, 1/2 circle, 27 mm	0	30"	Green	8886624561	36



Sutures

Needle chart


Reverse cutting
3/8 circle



Reverse Cutting
1/2 circle



Taper Point
1/2 circle



Maxon™ Monofilament Absorbable Suture continued

Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
GS-22	taper point, 1/2 circle, 27 mm	2-0	30"	Green	8886624551	36
GS-22	taper point, 1/2 circle, 27 mm	3-0	30"	Green	8886624541	36
GS-21	taper point, 1/2 Circle, 37 mm	0	30"	Green	8886626761	36
GS-21	taper point, 1/2 Circle, 37 mm	0	36"	Green	8886626961	36
GS-21	taper point, 1/2 Circle, 37 mm	1	30"	Green	8886626771	36
GS-21	taper point, 1/2 Circle, 37 mm	1	36"	Green	8886626971	36
GS-21	taper point, 1/2 Circle, 37 mm	2-0	30"	Green	8886626751	36
GS-21	taper point, 1/2 Circle, 37 mm	3-0	30"	Green	8886626741	36
GS-25	taper point, 1/2 Circle, 48 mm	1	60"	Green	8886630972	24
GS-26	taper point, 1/2 Circle, 65 mm	1	96"	Green	GMM344L	12
GS-26	taper point, 1/2 Circle, 65 mm	1	60"	Green	8886631972	24

Velosorb™ Braided Fast Absorbable Suture

Provides soft tissue approximation when only 7-10 days of wound support is required.

Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
C-13	reverse cutting, 3/8 Circle, 19 mm	3-0	30"	Undyed	SV2298	36
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	30"	Undyed	SV2297	36
C-14	reverse cutting, 3/8 Circle, 24 mm	2-0	30"	Undyed	SV2253	36
C-14	reverse cutting, 3/8 Circle, 24 mm	3-0	30"	Undyed	SV2252	36
C-14	reverse cutting, 3/8 Circle, 24 mm	4-0	30"	Undyed	SV2251	36
CV-23	taper point, 1/2 Circle, 17 mm	3-0	30"	Undyed	CV1215	36
CV-23	taper point, 1/2 Circle, 17 mm	5-0	30"	Undyed	CV1213	36

Polysorb™ Braided Absorbable Suture


Provides mid-term tensile strength with exclusive braiding process.

Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
P-10	premium reverse cutting, 3/8 Circle, 11 mm	6-0	18"	Undyed	SL1625	36
P-13	premium reverse cutting, 3/8 Circle, 13 mm	4-0	18"	Undyed	SL5690	36
P-13	premium reverse cutting, 3/8 Circle, 13 mm	5-0	18"	Undyed	SL5687G	12
PC-13	conventional cutting, 3/8 circle, 11 mm	5-0	18"	Undyed	SL1659	36
C-12	reverse cutting, 3/8 Circle, 12 mm	5-0	30"	Violet	SL630	36
C-1	reverse cutting, 3/8 Circle, 16 mm	5-0	18"	Undyed	SL172	36
C-13	reverse cutting, 3/8 Circle, 19 mm	3-0	30"	Violet	SL636	36


Sutures

Needle chart


**Reverse cutting
3/8 circle**



**Reverse Cutting
1/2 circle**



**Taper Point
1/2 circle**




Polysorb™ Braided Absorbable Suture continued

Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
C-13	reverse cutting, 3/8 Circle, 19 mm	3-0	30"	Undyed	SL822	36
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	30"	Violet	SL635	36
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	30"	Undyed	SL691	36
C-13	reverse cutting, 3/8 Circle, 19 mm	5-0	30"	Violet	SL634	36
C-13	reverse cutting, 3/8 Circle, 19 mm	5-0	30"	Undyed	SL690	36
C-14	reverse cutting, 3/8 Circle, 24 mm	2-0	30"	Violet	CL823	36
C-14	reverse cutting, 3/8 Circle, 24 mm	2-0	30"	Undyed	SL823	36
C-14	reverse cutting, 3/8 Circle, 24 mm	3-0	30"	Undyed	SL693	36
C-14	reverse cutting, 3/8 Circle, 24 mm	3-0	30"	Violet	SL694	36
C-14	reverse cutting, 3/8 Circle, 24 mm	4-0	30"	Undyed	SL692	36
C-16	reverse cutting, 3/8 Circle, 30 mm	2-0	60"	Undyed	CL977	36
C-16	reverse cutting, 3/8 Circle, 30 mm	2-0	30"	Undyed	SL611	36
C-16	reverse cutting, 3/8 Circle, 30 mm	3-0	30"	Undyed	SL696	36
C-17	reverse cutting, 3/8 Circle, 39 mm	2	36"	Undyed	SL675	36
C-23	reverse cutting, 1/2 circle, 24 mm	3-0	30"	Violet	SL632	36
GS-10	reverse cutting, 1/2 circle, 26 mm	3-0	30"	Undyed	CL868	36
GS-11	reverse cutting, 1/2 circle, 37 mm	0	30"	Undyed	CL544	36
GS-11	reverse cutting, 1/2 circle, 37 mm	0	30"	Violet	CL817	36
GS-11	reverse cutting, 1/2 circle, 37 mm	1	30"	Undyed	CL545	36
GS-11	reverse cutting, 1/2 circle, 37 mm	1	30"	Violet	CL818	36
GS-11	reverse cutting, 1/2 circle, 37 mm	2	36"	Undyed	CL533	36
GS-11	reverse cutting, 1/2 circle, 37 mm	2	36"	Violet	CL919	36
GS-11	reverse cutting, 1/2 circle, 37 mm	2-0	30"	Undyed	CL543	36
GS-11	reverse cutting, 1/2 circle, 37 mm	2-0	30"	Violet	CL816	36
CV-11	taper point, 3/8 circle, 13mm	5-0	30"	Violet	GL890	36
CVF-21	taper point, 1/2 circle, 12 mm	5-0	18"	Undyed	SL432	36
CV-23	taper point, 1/2 Circle, 17 mm	2-0	30"	Undyed	UL216	36
CV-23	taper point, 1/2 Circle, 17 mm	3-0	30"	Violet	UL204	36
CV-23	taper point, 1/2 Circle, 17 mm	3-0	30"	Undyed	UL215	36
CVF-23	taper point, 1/2 circle, 17 mm	4-0	30"	Violet	UL203	36
CV-23	taper point, 1/2 Circle, 17 mm	5-0	30"	Violet	UL202	36
CV-23	taper point, 1/2 Circle, 17 mm	5-0	30"	Undyed	UL213	36
CV-24	taper point, 1/2 circle, 20 mm	3-0	30"	Violet	CL461	36
CV-24	taper point, 1/2 circle, 20 mm	4-0	30"	Violet	CL460	36
CV-25	taper point, 1/2 circle, 22 mm	3-0	30"	Violet	GL182	36
CV-25	taper point, 1/2 circle, 22 mm	4-0	30"	Violet	GL181	36


Sutures

Needle chart


**Reverse cutting
3/8 circle**



**Reverse Cutting
1/2 circle**



**Taper Point
1/2 circle**




Polysorb™ Braided Absorbable Suture continued

Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
CV-25	taper point, 1/2 circle, 22 mm	4-0	30"	Violet	GL181	36
GS-23	taper point, 1/2 circle, 22 mm	3-0	30"	Violet	CL892	36
V-20	taper point, 1/2 circle, 26 mm	0	30"	Violet	GL124	36
V-20	taper point, 1/2 circle, 26 mm	2-0	36"	Violet	GL125	36
V-20	taper point, 1/2 circle, 26 mm	2-0	30"	Violet	GLS123	36
V-20	taper point, 1/2 circle, 26 mm	3-0	36"	Violet	GL226	36
V-20	taper point, 1/2 circle, 26 mm	3-0	36"	Violet	GL332	36
V-20	taper point, 1/2 circle, 26 mm	3-0	30"	Violet	GL67M	12
V-20	taper point, 1/2 circle, 26 mm	3-0	30"	Violet	GLS122	36
V-20	taper point, 1/2 circle, 26 mm	3-0	30"	Undyed	GLS322	36
V-20	taper point, 1/2 circle, 26 mm	4-0	30"	Violet	GLS121	36
V-20	taper point, 1/2 circle, 26 mm	4-0	30"	Undyed	GLS321	36
GS-22	taper point, 1/2 circle, 27 mm	0	30"	Violet	CL884	36
GS-22	taper point, 1/2 circle, 27 mm	2-0	30"	Violet	CL883	36
GS-22	taper point, 1/2 circle, 27 mm	3-0	30"	Violet	CL882	36
GS-21	taper point, 1/2 Circle, 37 mm	0	30"	Violet	CL812	36
GS-21	taper point, 1/2 Circle, 37 mm	0	36"	Violet	CL924	36
GS-21	taper point, 1/2 Circle, 37 mm	1	30"	Violet	CL813	36
GS-21	taper point, 1/2 Circle, 37 mm	1	36"	Violet	CL925	36
GS-21	taper point, 1/2 Circle, 37 mm	2	30"	Violet	CL814	36
GS-21	taper point, 1/2 Circle, 37 mm	2-0	30"	Violet	CL811	36
GS-21	taper point, 1/2 Circle, 37 mm	2-0	36"	Violet	CL923	36
GS-21	taper point, 1/2 Circle, 37 mm	3-0	36"	Violet	CL922	36
V-26	taper point, 1/2 circle, 37 mm	0	30"	Violet	GL128	36
V-26	taper point, 1/2 circle, 37 mm	3-0	30"	Violet	GL126	36
GS-26	taper point, 1/2 Circle, 65 mm	2	60"	Violet	CL825	24
GS-26	taper point, 1/2 Circle, 65 mm	2	30"	Violet	CL659M	6
	Non Needled Reel	2-0	60"	Violet	LL113	24
	Non Needled Reel	3-0	18"	Violet	L72	36


Sutures

Needle chart


Reverse cutting
3/8 circle



Reverse Cutting
1/2 circle



Taper Point
1/2 circle



Plain Gut Absorbable Suture

Natural gut fiber with short to mid term absorption with variable tensile strength retention and absorption time.

Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
C-13	reverse cutting, 3/8 Circle, 19 mm	3-0	30"	None	SG822	36
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	18"	None	SG752	36
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	30"	None	SG821	36
C-13	reverse cutting, 3/8 Circle, 19 mm	5-0	18"	None	SG751	36
C-13	reverse cutting, 3/8 Circle, 19 mm	5-0	30"	None	SG820	36
C-14	reverse cutting, 3/8 Circle, 24 mm	2-0	30"	None	SG823	36
C-14	reverse cutting, 3/8 Circle, 24 mm	3-0	30"	None	SG810	36
GS-21	taper point, 1/2 Circle, 37 mm	2-0	30"	None	CG843	36
GS-21	taper point, 1/2 Circle, 37 mm	3-0	30"	None	CG842	36






Chromic Gut Absorbable Suture

Natural gut fiber treated with chromic salt solution which offer short to mid term absorption.



Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
C-13	reverse cutting, 3/8 Circle, 19 mm	3-0	30"	None	SG636	36
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	30"	None	SG635	36
C-13	reverse cutting, 3/8 Circle, 19 mm	5-0	30"	None	SG634	36
C-14	reverse cutting, 3/8 Circle, 24 mm	2-0	30"	None	SG637	36
GS-10	reverse cutting, 1/2 circle, 26 mm	3-0	30"	None	CG415	36
GS-11	reverse cutting, 1/2 circle, 37 mm	0	30"	None	CG817	36
GS-11	reverse cutting, 1/2 circle, 37 mm	2-0	30"	None	CG816	36
CV-23	taper point, 1/2 Circle, 17 mm	3-0	30"	None	UG204	36
CV-23	taper point, 1/2 Circle, 17 mm	4-0	30"	None	UG203	36
CV-25	taper point, 1/2 circle, 22 mm	3-0	30"	None	GG182	36
V-20	taper point, 1/2 circle, 26 mm	0	30"	None	GG124	36
V-20	taper point, 1/2 circle, 26 mm	2-0	30"	None	GG123	36
V-20	taper point, 1/2 circle, 26 mm	3-0	30"	None	GG122	36
V-20	taper point, 1/2 circle, 26 mm	4-0	30"	None	GG121	36
GS-22	taper point, 1/2 circle, 27 mm	0	30"	None	CG884	36
GS-22	taper point, 1/2 circle, 27 mm	2-0	30"	None	CG883	36
GS-21	taper point, 1/2 Circle, 37 mm	0	30"	None	CG812	36
GS-21	taper point, 1/2 Circle, 37 mm	0	36"	None	CG924	36
GS-21	taper point, 1/2 Circle, 37 mm	1	36"	None	CG925	36
GS-21	taper point, 1/2 Circle, 37 mm	2-0	30"	None	CG811	36
GS-21	taper point, 1/2 Circle, 37 mm	2-0	36"	None	CG923	36

Needle chart





**Reverse cutting
3/8 circle**

- C-13  19 mm
- C-14  24 mm
- C-15  26 mm
- C-16  30 mm
- C-17  39 mm

**Reverse Cutting
1/2 circle**

- GS-10  26 mm
- GS-11  37 mm

**Taper Point
1/2 circle**

- CV-23  17 mm
- V-20  26 mm
- GS-22  27 mm
- GS-21  37 mm

Chromic Gut Absorbable Suture continued

Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
GS-21	taper point, 1/2 Circle, 37 mm	3-0	36"	None	CG922	36
GS-24	taper point, 1/2 circle, 40 mm	1	30"	None	CG803	36
GS-26	taper point, 1/2 Circle, 65 mm	2	60"	None	CG825L	24

Surgipro™ Monofilament Non-Absorbable Suture

Permanent polypropylene sutures indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological surgery.

Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
P-12	premium reverse cutting, 3/8 Circle, 19 mm	3-0	18"	Blue	SP5687	36
PC-12	conventional cutting, 3/8 circle, 19 mm	3-0	18"	Blue	SP1643	36
C-13	reverse cutting, 3/8 Circle, 19 mm	3-0	18"	Blue	SP665	36
C-13	reverse cutting, 3/8 Circle, 19 mm	3-0	18"	Blue	SP665G	12
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	18"	Blue	SP683	36
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	18"	Blue	SP683G	12
C-13	reverse cutting, 3/8 Circle, 19 mm	5-0	18"	Blue	SP661	36
C-13	reverse cutting, 3/8 Circle, 19 mm	5-0	18"	Blue	SP661G	12
C-14	reverse cutting, 3/8 Circle, 24 mm	3-0	18"	Blue	SP684G	12
C-14	reverse cutting, 3/8 Circle, 24 mm	4-0	18"	Blue	SP629	36
C-15	reverse cutting, 3/8 Circle, 26 mm	2-0	30"	Blue	SP685	36
C-15	reverse cutting, 3/8 Circle, 26 mm	2-0	18"	Blue	SP868	36
C-15	reverse cutting, 3/8 Circle, 26 mm	2-0	18"	Blue	SP868G	12
C-16	reverse cutting, 3/8 Circle, 30 mm	3-0	30"	Blue	SP675	36
C-17	reverse cutting, 3/8 Circle, 39 mm	0	30"	Blue	SP690	36
GS-18	reverse cutting, 3/8 Circle, 77 mm	2	30"	Blue	CP460	24
KV-34	tapercutting, 1/2 circle, 37 mm	0	30"	Blue	CP844	36
GS-21	taper point, 1/2 Circle, 37 mm	0	30"	Blue	CP424	36
GS-21	taper point, 1/2 Circle, 37 mm	1	30"	Blue	CP425	36
GS-21	taper point, 1/2 Circle, 37 mm	2-0	30"	Blue	CP423	36

Needle chart

**Reverse cutting
3/8 circle**



**Reverse Cutting
1/2 circle**



**Taper Point
1/2 circle**



Surgipro II™ Monofilament Non-Absorbable Suture






Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
P-13	premium reverse cutting, 3/8 Circle, 13 mm	4-0	18"	Blue	SP5699	36
P-13	premium reverse cutting, 3/8 Circle, 13 mm	5-0	18"	Blue	SP5698G	12
P-12	premium reverse cutting, 3/8 Circle, 19 mm	3-0	18"	Blue	SP5687G	12
P-12	premium reverse cutting, 3/8 Circle, 19 mm	4-0	18"	Blue	SP5682	36
CV-23	taper point, 1/2 Circle, 17 mm	3-0	30"	Blue	VP872X	36
CV-23	taper point, 1/2 Circle, 17 mm	5-0	30"	Blue	VP870X	36

Monosof™ Monofilament Nylon Non-Absorbable Suture



Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
PC-13	conventional cutting, 3/8 circle, 11 mm	6-0	18"	Black	SN1956	36
SCE-4	conventional cutting, 3/8 circle, 19 mm	4-0	18"	Black	SN568	36
C-1	reverse cutting, 3/8 Circle, 16 mm	5-0	18"	Black	SN871	36
C-13	reverse cutting, 3/8 Circle, 19 mm	2-0	36"	Black	SN666	36
C-13	reverse cutting, 3/8 Circle, 19 mm	3-0	30"	Black	SN641	36
C-13	reverse cutting, 3/8 Circle, 19 mm	3-0	18"	Black	SN653	36
C-13	reverse cutting, 3/8 Circle, 19 mm	3-0	36"	Black	SN665	36
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	36"	Black	SN644	36
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	18"	Black	SN662	36
C-13	reverse cutting, 3/8 Circle, 19 mm	5-0	18"	Black	SN661	36
C-14	reverse cutting, 3/8 Circle, 24 mm	2-0	18"	Black	SN758	36
C-14	reverse cutting, 3/8 Circle, 24 mm	3-0	36"	Black	SN648	36
C-14	reverse cutting, 3/8 Circle, 24 mm	3-0	18"	Black	SN663	36
C-14	reverse cutting, 3/8 Circle, 24 mm	3-0	30"	Black	SN669	36
C-14	reverse cutting, 3/8 Circle, 24 mm	4-0	18"	Black	SN658	36
C-15	reverse cutting, 3/8 Circle, 26 mm	2-0	30"	Black	SN664	36
C-15	reverse cutting, 3/8 Circle, 26 mm	2-0	18"	Black	SN764	36
C-16	reverse cutting, 3/8 Circle, 30 mm	2-0	30"	Black	SN688	36
C-16	reverse cutting, 3/8 Circle, 30 mm	3-0	30"	Black	SN671	36
C-17	reverse cutting, 3/8 Circle, 39 mm	2-0	30"	Black	SN674	36
SC-1	cutting, straight, 51 mm	3-0	36"	Black	SN616	36
SC-2	cutting, straight, 60 mm	2-0	30"	Black	SN628	36
SC-2	cutting, straight, 60 mm	3-0	30"	Black	SN627	36
SS-24	premium spatula, 1/4 circle, 95 degrees	5-0	18"	Black	N2781K	12

Needle chart





**Reverse cutting
3/8 circle**

- C-13  19 mm
- C-14  24 mm
- C-15  26 mm
- C-16  30 mm
- C-17  39 mm

**Reverse Cutting
1/2 circle**

- GS-10  26 mm
- GS-11  37 mm

**Taper Point
1/2 circle**




- CV-23  17 mm
- V-20  26 mm
- GS-22  27 mm
- GS-21  37 mm


Novafil™ Monofilament Non-Absorbable Suture

Unique dynamic suture made of polybutester. It stretches with edema and returns to its original shape when the swelling subsides.

Needle type	Needle description	Suture size	Suture length	Suture color	Item code	Box quantity
P-13	premium reverse cutting, 3/8 Circle, 13 mm	5-0	18"	Blue	SPB5223G	12
C-1	reverse cutting, 3/8 Circle, 16 mm	5-0	18"	Blue	8886441923	12
C-13	reverse cutting, 3/8 Circle, 19 mm	3-0	18"	Blue	8886442043	12
C-13	reverse cutting, 3/8 Circle, 19 mm	4-0	18"	Blue	8886442033	12
C-13	reverse cutting, 3/8 Circle, 19 mm	5-0	18"	Blue	8886442023	12
C-14	reverse cutting, 3/8 Circle, 24 mm	2-0	18"	Blue	8886442253	12
C-14	reverse cutting, 3/8 Circle, 24 mm	3-0	18"	Blue	8886442243	12
C-14	reverse cutting, 3/8 Circle, 24 mm	4-0	18"	Blue	8886442233	12
C-16	reverse cutting, 3/8 Circle, 30 mm	2-0	30"	Blue	8886442451	36
C-16	reverse cutting, 3/8 Circle, 30 mm	3-0	30"	Blue	8886442441	36
GS-11	reverse cutting, 1/2 circle, 37 mm	1	30"	Blue	8886443371	36
GS-22	taper point, 1/2 circle, 27 mm	0	30"	Blue	8886445961	36
GS-21	taper point, 1/2 Circle, 37 mm	0	30"	Blue	8886445561	36







V-Loc™ wound closure device

	V-Loc™ 90 device	V-Loc™ 180 device	V-Loc™ PBT nonabsorbable device
Tensile Strength	14 days	21 days	Permanent
Absorption Profile	90-110 days	180 days	Permanent
Procedural Applications	Soft tissue approximation where support is required consistent with the absorption profile	Soft tissue approximation where support is required consistent with the absorption profile	Soft tissue approximation
Color	Undyed, violet 	Undyed, green 	Blue 
Composition	Glycolide, dioxanone and trimethylene carbonate	Copolymer of glycolic acid and trimethylene carbonate	Copolymer of butylene terephthalate and polytetramethylene ether glycol
Indications	V-Loc™ 90 and V-Loc™ 180 absorbable wound closure devices are indicated for soft tissue approximation where use of an absorbable suture is appropriate.		V-Loc™ PBT nonabsorbable wound closure devices are indicated for soft tissue approximation.



V-Loc™ wound closure device






V-Loc™ 90 absorbable wound closure device

Premium reverse cutting 	Suture			Units per box	4-0	3-0	2-0	0
	Length (cm)	Length (in)	Color					
3/8 Circle 16 mm  P-11 5 mm Trocar	30	12	Undyed	12	VLOCM1413 ❖			
	45	18	Undyed	12	VLOCM1423 ❖			
3/8 Circle 19 mm  P-12 5 mm Trocar	15	6	Undyed	12	VLOCM0003 ❖	VLOCM0004 ❖	VLOCM0005	
	30	12	Undyed	12	VLOCM0013 ❖	VLOCM0014 ❖	VLOCM0015	
	45	18	Undyed	12	VLOCM0023 ❖	VLOCM0024 ❖	VLOCM0025	
	58	23	Undyed	12	VLOCM0033 ❖	VLOCM0034 ❖	VLOCM0035	
3/8 Circle 24 mm  P-14 8 mm Trocar	23	9	Violet	12			VLOCM3245	
	30	12	Undyed	12	VLOCM0113 ❖	VLOCM0114 ❖	VLOCM0115	
	45	18	Undyed	12	VLOCM0123 ❖	VLOCM0124 ❖	VLOCM0125	
	58	23	Undyed	12	VLOCM0133 ❖	VLOCM0134 ❖	VLOCM0135	
3/8 Circle 30 mm  P-17 11 mm Trocar	45	18	Undyed	12		VLOCM1824		
Taper point	Suture			Units per box	4-0	3-0	2-0	0
	Length (cm)	Length (in)	Color					
3/8 Circle 17 mm  CV-15 5 mm Trocar	15	6	Violet	12	VLOCM1203 ❖	VLOCM1204 ❖		

❖ Surgalloy™ Needle

V-Loc™ wound closure device

V-Loc™ 90 absorbable wound closure device

Premium reverse cutting 	Suture			Units per box	4-0	3-0	2-0	0
	Length (cm)	Length (in)	Color					
1/2 Circle 17 mm  CV-23 5 mm Trocar‡	15	6	Violet	12		VLOCM0804 ❖		
	15	6	Undyed	12		VLOCM1904 ❖		
	23	9	Violet	12	VLOCM0843 ❖	VLOCM0844 ❖		
	23	9	Undyed	12		VLOCM1944 ❖		
	30	12	Violet	12	VLOCM0813 ❖	VLOCM0814 ❖		
	45	18	Violet	12		VLOCM0824 ❖		
1/2 Circle 26 mm  V-20 8 mm Trocar‡	15	6	Violet	12	VLOCM0603 ❖	VLOCM0604 ❖		
	15	6	Undyed	12		VLOCM2004 ❖		
	23	9	Violet	12		VLOCM0644 ❖		
	23	9	Undyed	12		VLOCM2044 ❖		
	30	12	Violet	12		VLOCM0614 ❖		
	45	18	Violet	12	VLOCM0623 ❖	VLOCM0624 ❖	VLOCM0625	
	58	23	Violet	12			VLOCM0635	
1/2 Circle 30 mm  V-30 11 mm Trocar	23	9	Violet	12		VLOCM1544	VLOCM1545	
1/2 Circle 27 mm  GS-22 8 mm Trocar‡	15	6	Violet	12			VLOCM2105	VLOCM2106 ^Ω
	15	6	Undyed	12			VLOCM2205	
	23	9	Violet	12			VLOCM2145	VLOCM2146 ^Ω
	23	9	Undyed	12			VLOCM2245	
	30	12	Violet	12			VLOCM2115	VLOCM2116 ^Ω
	45	18	Violet	12				VLOCM2126 ^Ω
	55	22	Violet	12				VLOCM2136 ^Ω









❖ Surgalloy™ Needle

‡ Results can vary; depending on the trocar brand and manufacturer due to variations in trocar's internal diameter.

Ω New V-Loc™ device code.

V-Loc™ wound closure device

V-Loc™ 90 absorbable wound closure device






Premium reverse cutting 	Suture			Units per box	4-0	3-0	2-0	0
	Length (cm)	Length (in)	Color					
1/2 Circle 37 mm  GS-21 11 mm Trocar‡	15	6	Violet	12			VLOCM0305	VLOCM0306 ^Ω
	45	18	Violet	12			VLOCM0325	
	23	9	Violet	12		VLOCM0344	VLOCM0345	VLOCM0346 ^Ω
	30	12	Violet	12			VLOCM0315	VLOCM0316 ^Ω
	45	18	Violet	12		VLOCM0324	VLOCM0325	VLOCM0326 ^Ω
	45	18	Undyed	12			VLOCM3225	VLOCM3326 ^Ω
	55	22	Violet	12				VLOCM0336 ^Ω
	55	22	Undyed	12				VLOCM3336 ^Ω
	58	23	Violet	12			VLOCM0335	
1/2 Circle 40 mm  GS-24 15 mm Trocar	45	18	Violet	12				VLOCM4726 ^Ω
1/2 Circle 48 mm  GS-25 15 mm Trocar	55	22	Violet	12				VLOCM0436 ^Ω
5/8 Circle 27 mm  U-46 11 mm Trocar	15	6	Violet	12		VLOCM1704		
	23	9	Violet	12		VLOCM1744		
Straight cutting 	Suture			Units per box	4-0	3-0	2-0	0
	Length (cm)	Length (in)	Color					
Straight 51 mm  SC-1	45	18	Undyed	12		VLOCM1624		
Straight 60 mm  SC-2	45	18	Undyed	12	VLOCM0223	VLOCM0224		

‡ Results can vary; depending on the trocar brand and manufacturer due to variations in trocar's internal diameter.

Ω New V-Loc™ device code.

V-Loc™ wound closure device

V-Loc™ 180 absorbable wound closure device








Premium reverse cutting <input checked="" type="checkbox"/>	Suture			Units per box	4-0	3-0	2-0	0
	Length (cm)	Length (in)	Color					
3/8 Circle 16 mm  P-11 5 mm Trocar	30	12	Undyed	12	VLOCL1413 ❖			
3/8 Circle 19 mm  P-12 5 mm Trocar	15	6	Undyed	12	VLOCL0003 ❖	VLOCL0004		
	30	12	Undyed	12	VLOCL0013 ❖	VLOCL0014	VLOCL0015	
	45	18	Undyed	12	VLOCL0023 ❖	VLOCL0024	VLOCL0025	
	58	23	Undyed	12	VLOCL0033 ❖			
	60	24	Undyed	12		VLOCL0034		
3/8 Circle 24 mm  P-14 8 mm Trocar	30	12	Undyed	12		VLOCL0114	VLOCL0115	
	30	12	Green	12				VLOCL3216
	45	18	Undyed	12	VLOCL0123 ❖	VLOCL0124	VLOCL0125	
	58	23	Undyed	12	VLOCL0133 ❖			
	60	24	Undyed	12		VLOCL0134	VLOCL0135	
Premium reverse cutting <input checked="" type="checkbox"/>	Suture			Units per box	4-0	3-0	2-0	0
Length (cm)	Length (in)	Color						
1/2 Circle 37 mm  GS-11 12 mm Trocar‡	30	12	Undyed	12			VLOCL2815	
	45	18	Undyed	12				VLOCL2826
	45	18	Green	12			VLOCL2825	
	60	24	Green	12				VLOCL2836
1/2 Circle 40 mm  HOS-12 12 mm Trocar	60	24	Green	12				VLOCL4636

❖ Surgalloy™ Needle

‡ Results can vary; depending on the trocar brand and manufacturer due to variations in trocar's internal diameter.

V-Loc™ wound closure device

V-Loc™ 180 absorbable wound closure device

Premium reverse cutting 	Suture			Units per box	4-0	3-0	2-0	0
	Length (cm)	Length (in)	Color					
1/2 Circle 17 mm  CV-23 5 mm Trocar†	15	6	Green	12	VLOCL0803 ❖	VLOCL0804		
	15	6	Undyed	12		VLOCL1904 ❖		
	23	9	Green	12		VLOCL0844 ❖		
	30	12	Green	12	VLOCL0813 ❖	VLOCL0814		
	45	18	Green	12		VLOCL0824		
1/2 Circle 26 mm  V-20 8 mm Trocar†	15	6	Green	12	VLOCL0603 ❖	VLOCL0604 ❖		
	23	9	Green	12	VLOCL0643 ❖	VLOCL0644 ❖		
	30	12	Green	12	VLOCL0613 ❖	VLOCL0614 ❖	VLOCL0615	
	45	18	Green	12		VLOCL0624 ❖	VLOCL0625	
	60	24	Green	12			VLOCL0635	
1/2 Circle 30 mm  V-30 11 mm Trocar	23	9	Green	12		VLOCL1544	VLOCL1545	
1/2 Circle 27 mm  GS-22 8 mm Trocar†	15	6	Green	12			VLOCL2105	
	23	9	Green	12			VLOCL2145	VLOCL2246
	30	12	Green	12				VLOCL2216
1/2 Circle 37 mm  GS-21 11 mm Trocar†	15	6	Green	12			VLOCL0305	VLOCL0306
	23	9	Green	12		VLOCL0344	VLOCL0345	VLOCL0346
	30	12	Green	12			VLOCL0315	VLOCL0316
	45	18	Green	12		VLOCL0324	VLOCL0325	VLOCL0326
	60	24	Green	12			VLOCL0335	VLOCL0336
1/2 Circle 48 mm  GS-25 15 mm Trocar	30	12	Green	12				VLOCL0416
	45	18	Green	12				VLOCL0426
	60	24	Green	12				VLOCL0436

❖ Surgalloy™ Needle

† Results can vary; depending on the trocar brand and manufacturer due to variations in trocar's internal diameter.







V-Loc™ wound closure device

V-Loc™ 180 absorbable wound closure device

Taper point 	Suture			Units per box	3-0	2-0	0
	Length (cm)	Length (in)	Color				
ski 23 mm  SK	23	9	Green	12	VLOCL2744		
Straight cutting 	Length (cm)	Suture Length (in)	Color	Units per box	3-0	2-0	0
Straight 60 mm  SC-2	45	18	Clear	12	VLOCL0224		

V-Loc™ wound closure device

V-Loc™ PBT nonabsorbable wound closure device

Premium reverse cutting 	Suture			Units per box	4-0	3-0	2-0	0	1
	Length (cm)	Length (in)	Color						
1/2 Circle 17 mm  CV-23 5 mm Trocar‡	15	6	Blue	12		VLOCN0804			
	30	12	Blue	12		VLOCN0814 ❖			
1/2 Circle 26 mm  V-20 8 mm Trocar‡	15	6	Blue	12		VLOCN0604 ❖	VLOCN0605		
	23	9	Blue	12		VLOCN0644 ❖			
	30	12	Blue	12		VLOCN0614 ❖	VLOCN0615		
1/2 Circle 27 mm  GS-22 8 mm Trocar‡	15	6	Blue	12			VLOCN2105	VLOCN2106	
	23	9	Blue	12			VLOCN2145	VLOCN2146	VLOCN2147
	30	12	Blue	12					VLOCN2117
	45	18	Blue	12			VLOCN1125	VLOCN1126	
1/2 Circle 37 mm  GS-21 11 mm Trocar‡	15	6	Blue	12			VLOCN0305	VLOCN0306	
	23	9	Blue	12			VLOCN0345	VLOCN0346	VLOCN0347
	30	12	Blue	12					VLOCN0317
	45	18	Blue	12			VLOCN0325	VLOCN0326	VLOCN0327
1/2 Circle 48 mm  GS-25 15 mm Trocar	60	24	Blue	12					VLOCN0437

❖ Surgalloy™ Needle

‡ Results can vary; depending on the trocar brand and manufacturer due to variations in trocar's internal diameter.

Endo Stitch™ suturing device



Endo Stitch™ Suturing Device 10 mm

Endoscopic suturing device

Code	Description	Box quantity
173016	Endo Stitch™ Suturing Device 10 mm, 15" / 34 cm Shaft Length	3



Endo Stitch™ Triple-Stitch Loading Unit

Single use loading units for use with Endo Stitch™ 10 mm suturing device.

Each loading unit contains three strands of suture, each attached to a 9 mm needle.



Code	Suture material	USP size	Metric size	Color	Length	Box quantity
170043	Surgidac™ Ctd, Brd Polyester, Non Absorbable	0	3.5	Green	48" / 120 cm	3
170040	Surgidac™ Ctd, Brd Polyester, Non Absorbable	0	3.5	Green	7" / 18 cm	3
170041	Surgidac™ Ctd, Brd Polyester, Non Absorbable	2-0	3	Green	7" / 18 cm	3
170044	Surgidac™ Ctd, Brd Polyester, Non Absorbable	2-0	3	Green	48" / 120 cm	3

Endo Stitch™ suturing device


Endo Stitch™ Single-Stitch Loading Unit

Single use loading units for use with Endo Stitch™ 10 mm suturing device.



Code	Suture material	USP size	Metric size	Color	Length	Box quantity
170001	Sofsilk™ Ctd, Brd Silk	0	3.5	Black	7" / 18 cm	12
170003	Sofsilk™ Ctd, Brd Silk	0	3.5	Black	48" / 120 cm	12
170002	Sofsilk™ Ctd, Brd Silk	2-0	3	Black	7" / 18 cm	12
170004	Sofsilk™ Ctd, Brd Silk	2-0	3	Black	48" / 120 cm	12
170010	Bralon™ Ctd, Brd Nylon	0	3.5	Black	7" / 18 cm	12
170012	Bralon™ Ctd, Brd Nylon	0	3.5	Black	48" / 120 cm	12
170011	Bralon™ Ctd, Brd Nylon	2-0	3	Black	7" / 18 cm	12
170013	Bralon™ Ctd, Brd Nylon	2-0	3	Black	48" / 120 cm	12
173026	Surgidac™ Ctd, Brd Polyester	0	3.5	Green	7" / 18 cm	12
173024	Surgidac™ Ctd, Brd Polyester	0	3.5	Green	48" / 120 cm	12
173021	Surgidac™ Ctd, Brd Polyester	2-0	3	Green	7" / 18 cm	12
173023	Surgidac™ Ctd, Brd Polyester	2-0	3	Green	48" / 120 cm	12
170050	Polysorb™ Ctd, Brd Syn. Abs	0	3.5	Violet	7" / 18 cm	12
170052	Polysorb™ Ctd, Brd Syn. Abs	0	3.5	Violet	48" / 120 cm	12
170054	Polysorb™ Ctd, Brd Syn. Abs	0	3.5	Undyed	7" / 18 cm	12
170051	Polysorb™ Ctd, Brd Syn. Abs	2-0	3	Violet	7" / 18 cm	12
170053	Polysorb™ Ctd, Brd Syn. Abs	2-0	3	Violet	48" / 120 cm	12
170055	Polysorb™ Ctd, Brd Syn. Abs	2-0	3	Undyed	7" / 18 cm	12
170056	Polysorb™ Ctd, Brd Syn. Abs	0	3.5	Undyed	48" / 120 cm	12
170057	Polysorb™ Ctd, Brd Syn. Abs	2-0	3	Undyed	48" / 120 cm	12
170070	Polysorb™ Ctd, Brd Syn. Abs	3-0	2	Violet	17" / 18 cm	12
170071	Polysorb™ Ctd, Brd Syn. Abs	3-0	2	Violet	48" / 120 cm	12
170072	Polysorb™ Ctd, Brd Syn. Abs	3-0	2	Undyed	7" / 18 cm	12
170073	Polysorb™ Ctd, Brd Syn. Abs	3-0	2	Undyed	48" / 120 cm	12
170092	Polysorb™ Ctd, Brd Syn. Abs	4-0	1.5	Violet	7" / 18 cm	12
170090	Polysorb™ Ctd, Brd Syn. Abs	4-0	1.5	Violet	48" / 120 cm	12
170094	Polysorb™ Ctd, Brd Syn. Abs	4-0	1.5	Undyed	7" / 18 cm	12
170096	Polysorb™ Ctd, Brd Syn. Abs	4-0	1.5	Undyed	48" / 120 cm	12

Parietene™ macroporous mesh†



Code	Description	Box quantity
PPM1106	Monofilament polypropylene mesh, 11 x 6 cm	1
PPM1106X3	Monofilament polypropylene mesh, 11 x 6 cm	3
PPM1106X6	Monofilament polypropylene mesh, 11 x 6 cm	6
PPM1508	Monofilament polypropylene mesh, 15 x 8 cm	1
PPM1508X3	Monofilament polypropylene mesh, 15 x 8 cm	3
PPM1508X6	Monofilament polypropylene mesh, 15 x 8 cm	6
PPM1510	Monofilament polypropylene mesh, 15 x 10 cm	1
PPM1510X3	Monofilament polypropylene mesh, 15 x 10 cm	3
PPM1510X6	Monofilament polypropylene mesh, 15 x 10 cm	6
PPM1515	Monofilament polypropylene mesh, 15 x 15 cm	1
PPM1515X3	Monofilament polypropylene mesh, 15 x 15 cm	3
PPM1515X6	Monofilament polypropylene mesh, 15 x 15 cm	6
PPM2020	Monofilament polypropylene mesh, 20 x 20 cm	1
PPM2020X3	Monofilament polypropylene mesh, 20 x 20 cm	3
PPM3030	Monofilament polypropylene mesh, 30 x 30 cm	1
PPM3030X3	Monofilament polypropylene mesh, 30 x 30 cm	3
PPM4530	Monofilament polypropylene mesh, 45 x 30 cm	1
PPM4530X3	Monofilament polypropylene mesh, 45 x 30 cm	3
PPMK1106	Monofilament polypropylene mesh with keyhole, 11 x 6 cm	1
PPMK1106X3	Monofilament polypropylene mesh with keyhole, 11 x 6 cm	3
PPMK1106X6	Monofilament polypropylene mesh with keyhole, 11 x 6 cm	6
PPMK1508	Monofilament polypropylene mesh with keyhole, 15 x 8 cm	1
PPMK1508X3	Monofilament polypropylene mesh with keyhole, 15 x 8 cm	3
PPMK1508X6	Monofilament polypropylene mesh with keyhole, 15 x 8 cm	6

† Contact your local representative for availability

Skin, open, and laparoscopic stapling



Skin Staplers




Auto Suture™ Skin Stapler

Single use skin staplers with stainless steel staples

Code	Description	Ship case
8886803512	Skin Stapler with 35 Regular Staples	12
8886803712	Skin Stapler with 35 Wide Staples	12

Clip applicators



Premium Surgiclip™ III Clip Applier

Single use automatic clip applicator with super titanium clips

Code	Description	Ship case
133650	Premium Surgiclip™ III S- 9.0 Clip Applier, 20 Small Titanium Clips	6



Premium Surgiclip™ II Clip Applier

Single use automatic clip applicators with Super Interlock™ titanium clips

Code	Description	Ship case
134051	Premium Surgiclip™ II M-9.75 Clip Applier, 20 Medium Titanium Clips	6
134053	Premium Surgiclip™ II M-11.5 Clip Applier, 30 Medium Titanium Clips	6



Premium Surgiclip™ Clip Applier

Single use automatic clip applicators with Super Interlock™ titanium clips

Code	Description	Ship case
134048	Premium Surgiclip™ L-13.0 Clip Applier, 15 Large Titanium Clips	6

Surgical staplers



TA™ 30 Stapler with DST Series™ Technology

Directional Stapling Technology (DST)

Reloadable stapler with titanium staples

Code	Description	Color	Ship case
TA30V3S	TA™ 30-V3 (Vascular) Single Use Reloadable Stapler	White	3
TA3035S	TA™ 30 mm-3.5 mm Single Use Reloadable Stapler	Blue	3
TA3048S	TA™ 30 mm-4.8 mm Single Use Reloadable Stapler	Green	3



TA™ 30 disposable loading units (DLU) with titanium staples

Code	Description	Color	Ship case
TA30V3L	TA™ 30 mm-V3 (Vascular) Single Use Loading Unit	White	6
TA3035L	TA™ Reload, Single Use Loading Unit, 30 mm-3.5 mm	Blue	6
TA3048L	TA™ Reload, Single Use Loading Unit, 30 mm-4.8 mm	Green	6



TA™ 45 Stapler with DST Series™ Technology

Directional Stapling Technology (DST)

Reloadable stapler with titanium staples

Code	Description	Color	Ship case
TA4535S	TA™ 45 mm-3.5 mm Single Use Reloadable Stapler	Blue	3
TA4548S	TA™ 45 mm-4.8 mm Single Use Reloadable Stapler	Green	3



TA™ 45 disposable loading units (DLU) with titanium staples

Code	Description	Color	Ship case
TA4535L	TA™ Reload, Single Use Loading Unit, 45 mm-3.5 mm	Blue	6
TA4548L	TA™ Reload, Single Use Loading Unit, 45 mm-4.8 mm	Green	6



TA™ 60 Stapler with DST Series™ Technology

Directional Stapling Technology (DST)

Reloadable stapler with titanium staples

Code	Description	Color	Ship case
TA6035S	TA™ 60 mm-3.5 mm Single Use Reloadable Stapler	Blue	3
TA6048S	TA™ 60 mm-4.8 mm Single Use Reloadable Stapler	Green	3

Surgical staplers



TA™ 60 disposable loading units (DLU) with titanium staples

Code	Description	Color	Ship case
TA6035L	TA™ Reload, Single Use Loading Unit, 60 mm-3.5 mm	Blue	6
TA6048L	TA™ Reload, Single Use Loading Unit, 60 mm-4.8 mm	Green	6



TA™ 90 Stapler with DST Series™ Technology

Directional Stapling Technology (DST)

Reloadable stapler with titanium staples

Code	Description	Color	Ship case
TA9035s	TA™ 90 mm-3.5 mm Single Use Reloadable Stapler	Blue	3
TA9048s	TA™ 90 mm-4.8 mm Single Use Reloadable Stapler	Green	3



TA™ 90 disposable loading units (DLU) with titanium staples

Code	Description	Color	Ship case
TA9035L	TA™ Reload, Single Use Loading Unit, 90 mm-3.5 mm	Blue	6
TA9048L	TA™ Reload, Single Use Loading Unit, 90 mm-4.8 mm	Green	6

Product specifications

Stapler	Product description	Cartridges available	Staple length before closure	Staple leg length closed	Color	Number of reloads/total number of firings	Rows/number of staples	Staple cross section
TA30V3S	TA™ Stapler with DST Series™ Technology 30 mm	TA30V3L	2.5 mm	1.0 mm	White	7/8	3/23	.18 X .27 mm
TA3035S	TA™ Stapler with DST Series™ Technology 30 mm	TA3035L	3.5 mm	1.5 mm	Blue	7/8	2/11	.19 X .30 mm
TA3048S	TA™ Stapler with DST Series™ Technology 30 mm	TA3048L	4.8 mm	2.0 mm	Green	7/8	2/11	.24 X .35 mm
TA4535S	TA™ Stapler with DST Series™ Technology 45 mm	TA4535L	3.5 mm	1.5 mm	Blue	7/8	2/15	.19 X .30 mm
TA4548S	TA™ Stapler with DST Series™ Technology 45 mm	TA4548L	4.8 mm	2.0 mm	Green	7/8	2/15	.24 X .35 mm
TA6035S	TA™ Stapler with DST Series™ Technology 60 mm	TA6035L	3.5 mm	1.5 mm	Blue	7/8	2/21	.19 X .30 mm
TA6048S	TA™ Stapler with DST Series™ Technology 60 mm	TA6048L	4.8 mm	2.0 mm	Green	7/8	2/21	.24 X .35 mm
TA9035S	TA™ Stapler with DST Series™ Technology 90 mm	TA9035L	3.5 mm	1.5 mm	Blue	7/8	2/33	.19 X .30 mm
TA9048S	TA™ Stapler with DST Series™ Technology 90 mm	TA9048L	4.8 mm	2.0 mm	Green	7/8	2/33	.24 X .35 mm

Gastrointestinal anastomosis

GIA™ 60 Stapler with DST Series™ Technology

Directional Stapling Technology (DST)

Reloadable stapler with titanium staples



Code	Description	Color	Ship case
GIA6025S	GIA™ 60 mm-2.5 mm, Single Use Reloadable Stapler	White	3
GIA6038S	GIA™ 60 mm-3.8 mm, Single Use Reloadable Stapler	Blue	3
GIA6048S	GIA™ 60 mm-4.8 mm, Single Use Reloadable Stapler	Green	3

GIA™ 60 Loading Unit with DST Series™ Technology

Places two double staggered rows of titanium staples with dividing knife in between

Interchangeability of reloads within same stapler length



Code	Description	Color	Ship case
GIA6025L	GIA™ Reload 60 mm-2.5 mm, Single Use Loading Unit	White	6
GIA6038L	GIA™ Reload 60 mm-3.8 mm, Single Use Loading Unit	Blue	6
GIA6048L	GIA™ Reload 60 mm-4.8 mm, Single Use Loading Unit	Green	6

GIA™ 80 Stapler with DST Series™ Technology

Directional Stapling Technology (DST)

Reloadable stapler with titanium staples



Code	Description	Color	Ship case
GIA8038S	GIA™ 80 mm-3.8 mm, Single Use Reloadable Stapler	Blue	3
GIA8048S	GIA™ 80 mm-4.8 mm, Single Use Reloadable Stapler	Green	3

GIA™ 80 Loading Unit with DST Series™ Technology

Places two double staggered rows of titanium staples with dividing knife in between

Interchangeability of reloads within same stapler length



Code	Description	Color	Ship case
GIA8038L	GIA™ Reload 80 mm-3.8 mm, Single Use Loading Unit	Blue	6
GIA8048L	GIA™ Reload 80 mm-4.8 mm, Single Use Loading Unit	Green	6

Gastrointestinal anastomosis

GIA™ 100 Stapler with DST Series™ Technology

Directional Stapling Technology (DST)

Reloadable stapler with titanium staples



Code	Description	Color	Ship case
GIA10038S	GIA™ 100 mm-3.8 mm, Single Use Reloadable Stapler	Blue	3
GIA10048S	GIA™ 100 mm-4.8 mm, Single Use Reloadable Stapler	Green	3

GIA™ 100 Loading Unit with DST Series™ Technology

Places two double staggered rows of titanium staples with dividing knife in between

Interchangeability of reloads within same stapler length



Code	Description	Color	Ship case
GIA10038L	GIA™ Reload 100 mm-3.8 mm, Single Use Loading Unit	Blue	6
GIA10048L	GIA™ Reload 100 mm-4.8 mm, Single Use Loading Unit	Green	6

Product specifications

Stapler	Product description	Cartridges available	Staple length before closure	Staple leg length closed	Color	Number of reloads/total number of firings	Rows/number of staples	Staple cross section	Length of Staple Line
GIA6025S	GIA™ Stapler with DST Series™ Technology 60 mm	GIA6025L	2.5 mm	1.0 mm	White	7/8	64	4	66 mm
GIA6038S	GIA™ Stapler with DST Series™ Technology 60 mm	GIA6038L	3.8 mm	1.5 mm	Blue	7/8	64	4	66 mm
GIA6048S	GIA™ Stapler with DST Series™ Technology 60 mm	GIA6048L	4.8 mm	2.0 mm	Green	7/8	64	4	66 mm
GIA8038S	GIA™ Stapler with DST Series™ Technology 80 mm	GIA8038L	3.8 mm	1.5 mm	Blue	7/8	84	4	86 mm
GIA8048S	GIA™ Stapler with DST Series™ Technology 80 mm	GIA8048L	4.8 mm	2.0 mm	Green	7/8	84	4	86 mm
GIA10038S	GIA™ Stapler with DST Series™ Technology 100 mm	GIA10038L	3.8 mm	1.5 mm	Blue	7/8	104	4	106 mm
GIA10048S	GIA™ Stapler with DST Series™ Technology 100 mm	GIA10048L	4.8 mm	2.0 mm	Green	7/8	104	4	106 mm

Staplers




Endo GIA™ Ultra Universal Staplers

Single use instrument

Places two triple staggered rows of titanium staples and simultaneously, the knife divides the tissue in between

Code	Description	Ship case
EGIAUSHORT	Endo GIA™ 12 mm Ultra Universal Short Single Use Instrument	3
EGIAUSTND	Endo GIA™ 12 mm Ultra Universal Single Use Instrument	3
EGIAUXL	Endo GIA™ 12 mm Ultra Universal XL Single Use Instrument	3

Varied-height staples and stepped cartridge face



† Products depicted are not actual size.

Ordering information

Reorder code	Description	Cartridge color	Staple size (inner to outer row)
EGIA45AV‡	Endo GIA™ 45 mm Articulating Vascular Reload	Gray	2, 2, 2 mm
SIG45CTAV‡	Endo GIA™ 45 mm Curved Tip Articulating Vascular Reload	Gray	2, 2, 2 mm
SIG30AVM	Endo GIA™ 30 mm Articulating Vascular/Medium Reload with Tri-Staple™ Technology	Tan	2, 2.5, 3 mm
SIG30CTAVM	Endo GIA™ 30 mm Curved Tip Articulating Vascular/Medium Reload with Tri-Staple™ Technology	Tan	2, 2.5, 3 mm
EGIA45AVM	Endo GIA™ 45 mm Articulating Vascular/Medium Reload with Tri-Staple™ Technology	Tan	2, 2.5, 3 mm
SIG45CTAVM	Endo GIA™ 45 mm Curved Tip Articulating Vascular/Medium Reload with Tri-Staple™ Technology	Tan	2, 2.5, 3 mm
EGIA60AVM	Endo GIA™ 60 mm Articulating Vascular/Medium Reload with Tri-Staple™ Technology	Tan	2, 2.5, 3 mm
SIG60CTAVM	Endo GIA™ 60 mm Curved Tip Articulating Vascular/Medium Reload with Tri-Staple™ Technology	Tan	2, 2.5, 3 mm
SIG30AMT	Endo GIA™ 30 mm Articulating Medium/Thick Reload with Tri-Staple™ Technology	Purple	3, 3.5, 4 mm
EGIA45AMT	Endo GIA™ 45 mm Articulating Medium/Thick Reload with Tri-Staple™ Technology	Purple	3, 3.5, 4 mm
EGIA60AMT	Endo GIA™ 60 mm Articulating Medium/Thick Reload with Tri-Staple™ Technology	Purple	3, 3.5, 4 mm
SIG60CTAMT	Endo GIA™ 60 mm Curved Tip Articulating Medium/Thick Reload with Tri-Staple™ Technology	Purple	3, 3.5, 4 mm
SIG45XT	Endo GIA™ 45 mm Articulating Extra-Thick Reload with Tri-Staple™ Technology	Black	4, 4.5, 5 mm
SIG60XT	Endo GIA™ 60 mm Articulating Extra-Thick Reload with Tri-Staple™ Technology	Black	4, 4.5, 5 mm

‡ Reorder codes EGIA45AV and SIG45CTAV, while offering a slimmer fixed anvil, are vascular reloads that do not come with Tri-Staple™ Technology.

Clip Appliers

Endo Clip™ II 10 mm Clip Applier

Single use clip applier with Super Interlock™ titanium clips

Designed for introduction and use through an appropriately sized trocar sleeve, or larger with the use of a converter



Code	Description	Ship case
176657	Endo Clip™ II ML 10 mm Pistol Grip	6



Endo Clip™ III 5 mm Clip Applier with Clip Logic™ Technology

Single use clip applier with Clip Logic™ technology clips

Designed for introduction and use through an appropriately sized trocar sleeve, or larger with the use of a converter

Code	Description	Ship case
176630	Endo Clip™ III ML 5 mm	6

Electrosurgery and vessel sealing



Generators



ForceTriad™ Energy Platform

Delivers enhanced monopolar and bipolar energy with the generation of LigaSure™ tissue fusion technology all in one unit

LigaSure™ tissue fusion

TissueFect™ sensing technology

Automatic instrument recognition

Utilizes most LigaSure™ instruments

Unique Valleylab™ mode for enhanced dissection with hemostasis

Code	Description	Ship case
FORCETRIAD	ForceTriad™ Energy Platform	1

Footswitches



Valleylab™ Monopolar Footswitch

Compatible with all Force™ series and SurgiStat™ II generators

4.6 m (15 ft) cord

Code	Description	Ship case
E6008B	Valleylab™ Monopolar Footswitch	1




Valleylab™ Dome Bipolar Footswitch

Compatible with Force™ series generators

4.6 m (15 ft) cord

Code	Description	Ship case
E6019	Valleylab™ Dome Bipolar Footswitch	1

Footswitches



Monopolar Footswitch

Compatible with all Force™ series and SurgiStat™ II generators.
4.6m (15ft) cord

Code	Description	Ship case
E6008	Monopolar Footswitch	1



LigaSure™ Single Pedal Footswitch

Compatible with ForceTriad™ energy platform only.
4.6m (15ft) cord

Code	Description	Ship case
LF0500	LigaSure™ Single Pedal Footswitch	1



LigaSure™ Single Pedal Footswitch

Compatible with LigaSure™ vessel sealing system
and ForceTriad™ energy platform.
4.6m (15ft) cord

Code	Description	Ship case
LS0300	LigaSure™ Single Pedal Footswitch	1

LigaSure™ tissue fusion instruments

Sterile, single use, for open procedures†



LigaSure Impact™ Instrument Sealer/Divider, Nano-Coated

Length: 18.0 cm (7.1 in) Seal Length: 36 mm Jaw Design: Curved, 14 Degrees
 Shaft: 12.6 mm x 7.6 mm Cut Length: 34 mm
 Shaft Rotation: 180 degrees Hand or Foot Switching

Code	Description	Ship case
LF4418	LigaSure Impact™ Tissue Fusion Open Instrument	6



LigaSure Atlas™ Tissue Fusion Open Sealer/Divider

Length: 20 cm (7.9 in) Seal Length: 22 mm Hand or Foot Switching
 Shaft Diameter: 10 mm Cut Length: 20 mm Shaft Rotation: < 359 Degrees
 Jaw Design: Straight, Blunt

Code	Description	Ship case
LS1020	LigaSure Atlas™ 20 cm – Sealer/divider	6



Sealer/Divider

Length: 18.8 cm (7.4 in) Jaw Design: 28 Degrees Hand or Foot Switching
 Cut Length: 15.2 mm Seal Length: 17.3 mm

Code	Description	Ship case
LF1212A	LigaSure™ Small Jaw Instrument – Sealer/divider	6



BiZact™ device

Length: 12 cm (4.7 in.) Jaw Design: 24 Degrees Hand or Foot Switching
 Cut Length: 4.3 mm Seal Length: 4.9 mm

Code	Description	Ship case
BZ4212A	The BiZact™ device uses advanced bipolar energy to permanently seal vessels up to and including 3 mm in diameter.	6



LigaSure™ Maryland Jaw Open Vessel Sealer/Divider, Nano-Coated

Length: 23 cm (9.1 in) Seal Length: 20.3 mm Jaw Design: Maryland-style
 Shaft Diameter: 5 mm Cut Length: 18.5 mm Hand or Foot Switching
 Shaft Rotation: 350 Degrees

Code	Description	Ship case
LF1923	LigaSure™ Maryland Jaw Open Vessel Sealer/divider	6



LigaSure™ 5 mm Blunt Tip Open Sealer/Divider, Nano-Coated

Length: 23 cm (9.1 in) Seal Length: 18.3 mm Jaw Design: Straight, Blunt
 Shaft Diameter: 5 mm Cut Length: 16.8 mm Hand or Foot Switching
 Shaft Rotation: 180 Degrees

Code	Description	Ship case
LF1823	LigaSure™ 5 mm Blunt Tip Open Sealer/Divider	6

† Compatible with ForceTriad™ energy platform only.

LigaSure™ tissue fusion instruments


Sterile, single use, for laparoscopic procedures



LigaSure Atlas™ Tissue Fusion Laparoscopic Sealer/Divider

Length: 37 cm (14.6 in) Cut Length: 20 mm Shaft Rotation: < 359 Degrees
 Shaft Diameter: 10 mm Jaw Design: Straight, Blunt
 Seal Length: 22 mm Hand or Foot Switching


Code	Description	Ship case
LS1037	LigaSure Atlas™ Laparoscopic Instrument - Sealer/divider	6



LigaSure™ Maryland Jaw Vessel Sealer/Divider, Nano-Coated

Length: 37, 44 cm (14.6, 17.3 in) Seal Length: 20.3 mm Hand or Foot Switching
 Shaft Diameter: 5 mm Cut Length: 18.5 mm
 Shaft Rotation: 350 Degrees Jaw Design: Maryland-style

Code	Description	Ship case
LF1937	LigaSure™ Maryland Jaw Laparoscopic Vessel Sealer/Divider	6
LF1944	LigaSure™ Maryland Jaw Laparoscopic Vessel Sealer/Divider	6



LigaSure™ 5 mm Blunt Tip Sealer/Divider, Nano-Coated

Length: 37 cm (14.6 in) Seal Length: 18.3 mm Jaw Design: Straight, Blunt
 Shaft Diameter: 5 mm Cut Length: 16.8 mm Hand or Foot Switching
 Shaft Rotation: 180 Degrees

Code	Description	Ship case
LF1837	LigaSure™ 5 mm Blunt Tip Laparoscopic Sealer/Divider	6
LF1844	LigaSure™ 5 mm Blunt Tip Laparoscopic Sealer/Divider	6

Electrosurgical pencils



Force TriVerse™ Electrosurgical Device, Single Use

Cord Length: 10' (3 m) With Edge™ coated hex- Compatible with the
 Single use locking blade electrode ForceTriad™ energy
 Access to Valleylab™ mode Sterile-field power control platform.

Code	Description	Quantity
FT3000	Force TriVerse™ Electrosurgical Device, Holster, 10' (3 m)	25/case
FT3000DB	Force TriVerse™ Electrosurgical Device, Holster, 15' (4.6 m)	25/case

Rocker switch pencils



Rocker Switch Pencils, Single Use

Code	Description	Quantity
E2350H	Edge™ Rocker Switch Pencil, Blade Electrode, Holster, 10' (3 m)	50/case
E2350HDB	Edge™ Rocker Switch Pencil, Blade Electrode, Holster, 15' (4.6 m)	50/case
E2515	Valleylab™ Rocker Switch Pencil, 10' (3 m)	50/case
E2515H	Valleylab™ Rocker Switch Pencil, Holster, 10' (3 m)	50/case
E2515HDA	Valleylab™ Rocker Switch Pencil, Needle Electrode, Holster, 10' (3 m)	50/case
E2515HDB	Valleylab™ Rocker Switch Pencil, Holster, 15' (4.6 m)	50/case

Button switch pencils



Button Switch Pencils, Single Use

Code	Description	Quantity
E2450H	Edge™ Button Switch Pencil, Blade Electrode, Holster, 10' (3 m)	50/case
E2516	Valleylab™ Button Switch Pencil, 10' (3 m)	50/case
E2516H	Valleylab™ Button Switch Pencil, Holster, 10' (3 m)	50/case
E2516HDA	Valleylab™ Button Switch Pencil, Needle Electrode, Holster, 10' (3 m)	50/case

Assorted pencils and accessories



Rocker Switch Pencil, Reusable

Code	Description	Quantity
E2100	Valleylab™ Rocker Switch Pencil, Reusable, 15' (4.6 m)	10/case

Nonsterile. Steam or EtO sterilizable.

Coated electrodes, single use



Blade Electrodes

Code	Description	Quantity
E1450G	Edge™ Blade Electrode, 3" (7.62 cm)	50/case
E1450X	Edge™ Hex-locking Blade Electrode, 2.5" (6.35 cm)	50/case
E14504	Edge™ Blade Electrode, 4" (10.16 cm)	50/case
E14506	Edge™ Blade Electrode, 6.5" (16.51 cm)	50/case
E1475X	Edge™ Hex-locking Blade Electrode, 3" (7.62 cm)	50/case



Needle Electrodes

Code	Description	Quantity
E1452	Edge™ Needle Electrode, 2.8" (7.2 cm)	50/case
E14526	Edge™ Needle Electrode, 6.5" (16.51 cm)	50/case

Insulated coated electrodes, single use



Blade Electrodes

Code	Description	Quantity
E1455	Edge™ Insulated Blade Electrode, 2.75" (7.0 cm)	50/case
E1455B	Edge™ Insulated Blade Electrode, Safety Sleeve™, 2.75" (7.0 cm)	25/case
E14554	Edge™ Insulated Blade Electrode, 4" (10.16 cm)	50/case
E1455B4	Edge™ Insulated Blade Electrode, Safety Sleeve™, 4" (10.16 cm)	25/case
E14556	Edge™ Insulated Blade Electrode, 6.5" (16.51 cm)	50/case



Needle Electrodes

Code	Description	Quantity
E1465	Edge™ Insulated Needle Electrode, 2.8" (7.2 cm)	50/case
E1465B	Edge™ Insulated Needle Electrode, Safety Sleeve™, 2.8" (7.2 cm)	25/case
E14654	Edge™ Insulated Needle Electrode, 4" (10.16 cm)	50/case
E1465B4	Edge™ Insulated Needle Electrode, Safety Sleeve™, 4" (10.16 cm)	25/case
E14656	Edge™ Insulated Needle Electrode, 6.5" (16.51 cm)	50/case

Stainless steel electrodes, single use

Blade Electrodes



Code	Description	Quantity
E1551G	Valleylab™ Blade Electrode, 2.4" (6.2 cm)	150/case
E1551X	Valleylab™ Hex-locking Blade Electrode, 2.4" (6.2 cm)	150/case
E15516	Valleylab™ Blade Electrode, 6.5" (16.51 cm)	50/case

Needle Electrodes



Code	Description	Quantity
E1552	Valleylab™ Needle Electrode, 2.8" (7.2 cm)	150/case
E15526	Valleylab™ Needle Electrode, 6.5" (16.51 cm)	50/case

Ball Electrodes



Code	Description	Quantity
E1550	Valleylab™ Ball Electrode, 5 mm	150/case
E1563	Valleylab™ Ball Electrode, 3 mm	10/case
E1564	Valleylab™ Ball Electrode, 5 mm	10/case

Stainless steel electrodes, reusable

Blade Electrodes



Code	Description	Quantity
E1020	Valleylab™ Hex-locking Blade Electrode, Reusable, 2.75" (7.0 cm)	1 each

Needle Electrodes



Code	Description	Quantity
E1003	Valleylab™ Fine Needle Electrode, Reusable, 2.6" (6.6 cm)	1 each

Ball Electrodes



Code	Description	Quantity
E1000	Valleylab™ Ball Electrode, Reusable, 3.2 mm	1 each
E1002	Valleylab™ Ball Electrode, Reusable, 5.6 mm	1 each

Loop Electrodes



Code	Description	Quantity
E1005	Valleylab™ Tungsten Loop Electrode, Reusable, 9.5 mm	1 each

Electrode extensions, reusable



Electrode Extensions

Code	Description	Quantity
E1502	Valleylab™ Straight Extension, 13 cm	1 each
E1504	Valleylab™ Straight Extension, 34.3 cm	1 each

Tungsten microsurgical needle electrodes, single use



Needle Electrodes

Code	Description	Quantity
E1650	Valleylab™ Straight Micro-Needle Electrode, 2 cm	10/case
E1651	Valleylab™ Straight Micro-Needle Electrode, 3 cm	10/case
E1652	Valleylab™ Angled Micro-Needle Electrode, 45°, 3 mm from tip, 3 cm	10/case
E1653	Valleylab™ Angled Micro-Needle Electrode, 45°, 1 cm from tip, 3 cm	10/case
E1654B	Valleylab™ Insulated Micro-Needle Electrode, 9 cm	10/case

Arthroscopic electrodes, single use



Arthroscopic Electrodes

Code	Description	Quantity
E1510	Valleylab™ Arthroscopic Hook Electrode, 90°	10/case
E1511	Valleylab™ Arthroscopic Hook Electrode, 45°	10/case
E1512	Valleylab™ Insulated Arthroscopic Blade Electrode	10/case

Tungsten loop electrodes, single use



Round Loops

Code	Description	Quantity
E1559	Valleylab™ Tungsten Loop Electrode, 10 mm x 10 mm	10/case
E1560	Valleylab™ Tungsten Loop Electrode, 15 mm x 12 mm	10/case
E1561	Valleylab™ Tungsten Loop Electrode, 20 mm x 12 mm	10/case
E1562	Valleylab™ Tungsten Loop Electrode, 20 mm x 15 mm	10/case

Ball Electrodes



Code	Description	Quantity
E1563	Valleylab™ Ball Electrode, 3 mm	10/case
E1564	Valleylab™ Ball Electrode, 5 mm	10/case

Square Loops



Code	Description	Quantity
1p6	Valleylab™ Tungsten Square Loop Electrode, 5 mm x 5 mm	10/case
E1567	Valleylab™ Tungsten Square Loop Electrode, 10 mm x 8 mm	10/case

Patient return electrodes



Adult REM™ Patient Return Electrodes, Single Use

For patients > 30 lbs (> 13.6 kg)

For use with RECQM generators

Code	Description	Quantity
E7507	Valleylab™ REM Polyhesive™ Adult Patient Return Electrode, 9' (2.7 m)	50/case



REM™ Infant Patient Return Electrodes, Single Use

For patients 6-30 lbs (2.72 - 13.6 kg)

For use with RECQM generators

Code	Description	Quantity
E751025	Valleylab™ REM Polyhesive™ Infant Patient Return Electrode, 9' (2.7 m)	25/case



REM™ Neonatal Patient Return Electrode, Single Use

For patients 1-6 lbs (0.45 kg - 2.72 kg)

For use with RECQM generators

Code	Description	Quantity
E7512	Valleylab™ REM Polyhesive™ Neonatal Patient Return Electrode, 9' (2.7 m)	12/case

Forceps, reusable



Code	Description	Quantity
E4051CT	Valleylab™ Bayonet Forceps, Scoville-Greenwood, 19.7 cm-1.5 mm	1 each



Code	Description	Quantity
E4052CT	Valleylab™ Forceps, Cushing, 15.2 cm-0.7 mm	1 each



Code	Description	Quantity
E4053CT	Valleylab™ Forceps, Cushing, 17.8 cm-1.5 mm	1 each



Code	Description	Quantity
E4054CT	Valleylab™ Bayonet Forceps, Cushing, 19.1 cm-2.0 mm	1 each



Code	Description	Quantity
E4055CT	Valleylab™ Straight Tip Forceps, Jewelers, 10.2 cm-0.4 mm	1 each



Code	Description	Quantity
E4057CT	Valleylab™ Bayonet Forceps, Gerald, 19.1 cm-0.7 mm	1 each



Code	Description	Quantity
E4058CT	Valleylab™ Bayonet Forceps, Cushing, 19.1 cm-0.7 mm	1 each



Code	Description	Quantity
E4059CT	Valleylab™ Bayonet Forceps with Stops, Hardy, 20.9 cm-0.5 mm	1 each



Code	Description	Quantity
E4060CT	Valleylab™ Forceps with Stops, Semkin, 14.0 cm-0.5 mm	1 each



Code	Description	Quantity
E4062CT	Valleylab™ Curved Forceps, Iris, 8.9 cm-0.5 mm	1 each



Code	Description	Quantity
E4073CT	Valleylab™ Bayonet Forceps, Angled Titanium, 22.2 cm-1.0 mm	1 each

Cords, single use



Code	Description	Quantity
E0503	Valleylab™ Monopolar Cord, 10' (3 m)	50/case



Code	Description	Quantity
E0510	Valleylab™ Laparoscopic Monopolar Cord, 10' (3 m)	50/case



Code	Description	Quantity
E0512	Valleylab™ Bipolar Forceps Cord, 12' (3.6 m)	50/case



Code	Description	Quantity
FT0510 ⁺	ForceTriad™ Monopolar Cord, 15' (4.6 m)	25/case

Cords, reusable



Code	Description	Quantity
E0020V	Valleylab™ Bipolar Forceps Cord, Reusable, 15' (4.6 m)	1 each



Code	Description	Quantity
E0021S	Valleylab™ Bipolar Forceps Cord, Storz, Reusable, 15' (4.6 m)	1 each



Code	Description	Quantity
E0022W	Valleylab™ Bipolar Forceps Cord, Wolf, Reusable, 15' (4.6 m)	1 each



Code	Description	Quantity
E0560	Valleylab™ REM™ Patient Return Electrode Cord, Reusable, 15' (4.6 m)	1 each



Code	Description	Quantity
FT0021S	ForceTriad™ Bipolar Resection Cord, Storz, Reusable, 15' (4.6 m)	1 each



Code	Description	Quantity
FT0022W	ForceTriad™ Bipolar Resection Cord, Wolf, Reusable, 15' (4.6 m)	1 each

Active and specialty adapters



Code	Description	Quantity
------	-------------	----------

E0017	Valleylab™ Universal Monopolar Adapter	1 each
-------	--	--------




E0021B	Valleylab™ Macrobipolar Adapter	1 each
--------	---------------------------------	--------



E05021	Valleylab™ Monopolar Adapter	1 each
--------	------------------------------	--------

REM adapters



Code	Description	Quantity
------	-------------	----------

E0600	Valleylab™ REM Patient Return Electrode Adapter	1 each
-------	---	--------

Non-REM adapters



Code	Description	Quantity
------	-------------	----------

E0507B	Valleylab™ Multiple Return/S Cord Adapter	1 each
--------	---	--------

E05042	Valleylab™ Patient Return Electrode Adapter, Non-REM	1 each
--------	--	--------

Solid PTFE coated laparoscopic electrodes, single use

For use with standard electrosurgical pencils



Code	Description	Quantity
E3770-36C	CleanCoat™ Laparoscopic Straight Spatula Electrode, 36 cm, coated	10/case



Code	Description	Quantity
E3770-36C	CleanCoat™ Laparoscopic Straight Spatula Electrode, 36 cm, coated	10/case



Code	Description	Quantity
E3772-36C	CleanCoat™ Laparoscopic Wire J-Hook Electrode, 36 cm, coated	10/case



Code	Description	Quantity
E3773-36C	CleanCoat™ Laparoscopic Wire L-Hook Electrode, 36 cm, coated	10/case

E3773-45C	CleanCoat™ Laparoscopic Wire L-Hook Electrode, 45 cm, coated	10/case
-----------	--	---------



Code	Description	Quantity
E3774-36C	CleanCoat™ Laparoscopic Flat L-Hook Electrode, 36 cm, coated	10/case

E3774-45C	CleanCoat™ Laparoscopic Flat L-Hook Electrode, 45 cm, coated	10/case
-----------	--	---------

Solid stainless steel laparoscopic electrodes, single use

For use with standard electrosurgical pencils



Code	Description	Quantity
E3770-36	Valleylab™ Laparoscopic Straight Spatula Electrode, Solid, 36 cm	10/case



Code	Description	Quantity
E3771-36	Valleylab™ Laparoscopic Curved Spatula Electrode, Solid, 36 cm	10/case

E3771-45	Valleylab™ Laparoscopic Curved Spatula Electrode, Solid, 45 cm	10/case
----------	--	---------



Code	Description	Quantity
E3772-36	Valleylab™ Laparoscopic Wire J-Hook Electrode, Solid, 36 cm	10/case



Code	Description	Quantity
E3773-36	Valleylab™ Laparoscopic Wire L-Hook Electrode, Solid, 36 cm	10/case

E3773-45	Valleylab™ Laparoscopic Wire L-Hook Electrode, Solid, 45 cm	10/case
----------	---	---------



Code	Description	Quantity
E3774-36	Valleylab™ Laparoscopic Flat L-Hook Electrode, Solid, 36 cm	10/case

Hollow laparoscopic electrodes, single use

Retractable

For suction and irrigation performance with the four-function handset



Code	Description	Quantity
E2781R28ASP	Valleylab™ Laparoscopic Curved Spatula Electrode, Retractable, Aspiration Holes, 28 cm	10/case



Code	Description	Quantity
E2782R28ASP	Valleylab™ Laparoscopic Wire J-Hook Electrode, Retractable, Aspiration Holes, 28 cm	10/case



Code	Description	Quantity
E2783R28ASP	Valleylab™ Laparoscopic Wire L-Hook Electrode, Retractable, Aspiration Holes, 28 cm	10/case



Code	Description	Quantity
E2783R36ASP	Valleylab™ Laparoscopic Wire L-Hook Electrode, Retractable, Aspiration Holes, 36 cm	10/case

Code	Description	Quantity
E2784R28ASP	Valleylab™ Laparoscopic Flat L-Hook Electrode, Retractable, Aspiration Holes, 28 cm	10/case

Hollow laparoscopic electrodes, single use

Non-Retractable

Standard electrodes for use with the four-function handset



Code	Description	Quantity
E278028	Valleylab™ Laparoscopic Straight Spatula Electrode, Hollow, 28 cm	10/case



Code	Description	Quantity
E278128	Valleylab™ Laparoscopic Curved Spatula Electrode, Hollow, 28 cm	10/case



Code	Description	Quantity
E278228	Valleylab™ Laparoscopic Wire J-Hook Electrode, Hollow, 28 cm	10/case



Code	Description	Quantity
E278328	Valleylab™ Laparoscopic Wire L-Hook Electrode, Hollow	10/case



Code	Description	Quantity
E278628	Valleylab™ Laparoscopic Cylindrical Tip Electrode, Hollow, 28 cm	10/case

Smoke evacuation





RapidVac™ Smoke Evacuator

Code	Description	Quantity
SE3690	RapidVac™ Smoke Evacuator, 110 V	1 each



Valleylab™ Laparoscopic Smoke Evacuation System



Code	Description	Quantity
SEL7000A	Valleylab™ Laparoscopic Smoke Evacuation System	10/case

RapidVac™ Smoke Evacuator Accessories

Code	Description	Quantity
SEA3700	RapidVac™ ULPA Filter	2/case
SEA3730	RapidVac™ Interlink Cable	1 each
SEA3741	RapidVac™ Remote Switch Activator	1 each
SEA3745	RapidVac™ Footswitch	1 each
RVFT10	RapidVac™ Fluid Trap	1 each
E3810	Valleylab™ Smoke Attachment 10' smoke tubing	20/case



Patient monitoring



Nellcor™ Portable SpO₂ Patient Monitoring System, PM10N

Code	Description	Quantity
V-SAT	Veterinary Tongue Sensor with Small and Large Clip	1 each
PMAC10N-P	Portable SpO ₂ Protective Cover (Pink)	1 each
PMAC10N-B	Portable SpO ₂ Protective Cover (Blue)	1 each
PMAC10N-N	Portable SpO ₂ Protective Cover (Navy)	1 each
PMAC10N-G	Portable SpO ₂ Protective Cover (Green)	1 each
PMAC10N-T	Portable SpO ₂ Transport Protective Cover	1 each
PMAC10N-CC	Portable SpO ₂ Carrying Case	1 each


Nellcor™ Bedside SpO₂ Patient Monitoring System, PM100N

Code	Description	Quantity
PM100N	Nellcor™ Bedside SpO ₂ Patient Monitoring System	1 each



Capnostream™ 35 Portable Respiratory Monitor

Code	Description	Quantity
PM35MN	Capnostream™ 35 Portable Respiratory Monitor	1 each



RespArray Nellcor™ Kit


Code	Description	Quantity
RESPDOC10KITAHA	RespArray Nellcor Kit includes RespArray patient monitor, Quick Start Kit, DOC-10 cable	1 each

Access and trocar instruments




Optical trocars

VersaOne™ Optical Trocars



Code	Description	Quantity
ONB5SHF	VersaOne™ Optical Trocar with Fixation Cannula, Short Length 5 mm	6
ONB5STF	VersaOne™ Optical Trocar with Fixation Cannula, Standard Length 5 mm	6
ONB5STF2C	VersaOne™ Optical Trocar with an Extra Standard Fixation Cannula, Standard Length 5 mm	6
ONB5LGF	VersaOne™ Optical Trocar with Fixation Cannula, Long Length 5 mm	6
ONB11STF	VersaOne™ Optical Trocar with Fixation Cannula, Standard Length 11 mm	6
ONB11LGF	VersaOne™ Optical Trocar with Fixation Cannula, Long Length 11 mm	6
ONB12SHF	VersaOne™ Optical Trocar with Fixation Cannula, Short Length 12 mm	6
ONB12STF	VersaOne™ Optical Trocar with Fixation Cannula, Standard Length 12 mm	6
ONB12STS	VersaOne™ Optical Trocar with Smooth Cannula, Standard Length 12 mm	6
ONB12LGF	VersaOne™ Optical Trocar with Fixation Cannula, Long Length 12 mm	6
ONB15STF	VersaOne™ Optical Trocar with Fixation Cannula, Standard Length 15 mm	6

VersaOne™ Universal Cannula



Code	Description	Quantity
UNVCA5SHF	VersaOne™ Universal Cannula, Fixation, Short Length 5 mm	6
UNVCA5STF	VersaOne™ Universal Cannula, Fixation, Standard Length 5 mm	6
UNVCA5LGF	VersaOne™ Universal Cannula, Fixation, Long Length 5 mm	6
UNVCA11STF	VersaOne™ Universal Cannula, Fixation, Standard Length 11 mm	6
UNVCA12STF	VersaOne™ Universal Cannula, Fixation, Standard Length 12 mm	6
UNVCA12STS	VersaOne™ Universal Cannula, Smooth, Standard Length 12 mm	6
UNVCA12LGF	VersaOne™ Universal Cannula, Fixation, Long Length 12 mm	6

Optical trocars



Visiport™ Plus RPF Optical Trocars

Code	Description	Quantity
176673P	Visiport™ Plus RPF 5 mm-11 mm Standard Optical Trocar with Smooth Cannula	3
176674PF	Visiport™ Plus RPF 5 mm-12 mm Standard Optical Trocar with Fixation Cannula	3
176674P	Visiport™ Plus RPF 5 mm-12 mm Standard Optical Trocar with Smooth Cannula	3


Bladed trocars



VersaOne™ Bladed Trocar

Single Use Trocar System with Universal Seal

Code	Description	Quantity
B5SHF	VersaOne™ Bladed Trocar with Fixation Cannula, Short Length 5 mm	6
B5SHS	VersaOne™ Bladed Trocar with Smooth Cannula, Short Length 5 mm	6
B5STF	VersaOne™ Bladed Trocar with Fixation Cannula, Standard Length 5 mm	6
B5STS	VersaOne™ Bladed Trocar with Smooth Cannula, Standard Length 5 mm	6
B11STF	VersaOne™ Bladed Trocar with Fixation Cannula, Standard Length 11 mm	6
B11STS	VersaOne™ Bladed Trocar with Smooth Cannula, Standard Length 11 mm	6
B12STF	VersaOne™ Bladed Trocar with Fixation Cannula, Standard Length 12 mm	6
B12STS	VersaOne™ Bladed Trocar with Smooth Cannula, Standard Length 12 mm	6
B12LGF	VersaOne™ Bladed Trocar with Fixation Cannula, Long Length 12 mm	6
B12LGS	VersaOne™ Bladed Trocar with Smooth Cannula, Long Length 12 mm	6



VersaOne™ Universal Cannula

Code	Description	Quantity
UNVCA5SHF	VersaOne™ Universal Cannula, Fixation, Short Length 5 mm	6
UNVCA5STF	VersaOne™ Universal Cannula, Fixation, Standard Length 5 mm	6
UNVCA5LGF	VersaOne™ Universal Cannula, Fixation, Long Length 5 mm	6
UNVCA11STF	VersaOne™ Universal Cannula, Fixation, Standard Length 11 mm	6
UNVCA12STF	VersaOne™ Universal Cannula, Fixation, Standard Length 12 mm	6
UNVCA12STS	VersaOne™ Universal Cannula, Smooth, Standard Length 12 mm	6
UNVCA12LGF	VersaOne™ Universal Cannula, Fixation, Long Length 12 mm	6



Autosuture™ Locking Trocar

Single use trocars

Code	Description	Quantity
240-55	Autosuture™ Locking Trocar	5

Bladeless trocars


VersaOne™ Bladeless Trocars

Port fixation and security



Code	Description	Quantity
NONB5SHF	VersaOne™ Bladeless Trocar with Fixation Cannula, Short Length 5 mm	6
NONB5STF	VersaOne™ Bladeless Trocar with Fixation Cannula, Standard Length 5 mm	6
NONB5LGF	VersaOne™ Bladeless Trocar with Fixation Cannula, Long Length 5 mm	6
NONB8STF	VersaOne™ Bladeless Trocar with Fixation Cannula, Standard Length 8 mm	6
NONB11STF	VersaOne™ Bladeless Trocar with Fixation Cannula, Standard Length 11 mm	6
NONB12SHF	VersaOne™ Bladeless Trocar with Fixation Cannula, Short Length 12 mm	6
NONB12STF	VersaOne™ Bladeless Trocar with Fixation Cannula, Standard Length 12 mm	6
NONB12STS	VersaOne™ Bladeless Trocar with Smooth Cannula, Standard Length 12 mm	6
NONB12LGF	VersaOne™ Bladeless Trocar with Fixation Cannula, Long Length 12 mm	6
NONB15STF	VersaOne™ Bladeless Trocar with Fixation Cannula, Standard Length 15 mm	6
NONB15LGF	VersaOne™ Bladeless Trocar with Fixation Cannula, Long Length 15 mm	6


VersaOne™ Universal Cannula



Code	Description	Quantity
UNVCA5SHF	VersaOne™ Universal Cannula, Fixation, Short Length 5 mm	6
UNVCA5STF	VersaOne™ Universal Cannula, Fixation, Standard Length 5 mm	6
UNVCA5LGF	VersaOne™ Universal Cannula, Fixation, Long Length 5 mm	6
UNVCA11STF	VersaOne™ Universal Cannula, Fixation, Standard Length 11 mm	6
UNVCA12STF	VersaOne™ Universal Cannula, Fixation, Standard Length 12 mm	6
UNVCA12STS	VersaOne™ Universal Cannula, Smooth, Standard Length 12 mm	6
UNVCA12LGF	VersaOne™ Universal Cannula, Fixation, Long Length 12 mm	6

Blunt tip trocars

Blunt tip trocars with fixation balloon, inflated balloon 25 cc
Designed for surgeons that use a Hasson or open technique



Code	Description	Quantity
OMST10BTNL	10 mm Single Use Sterile Blunt Tip Trocar Built-in Converters, Sizes 5 mm and 7 mm/8 mm Latex-free, Standard Length	5
OMST10BTSNL	10 mm Single Use Sterile Blunt Tip Trocar Built-in Converters, Sizes 5 mm and 7 mm/8 mm Latex-free, Short Length	5
OMST12BTNL	12 mm Single Use Sterile Blunt Tip Trocar Built-in Converters, Sizes 5 mm and 7 mm/8 mm Latex-free, Standard Length	5
OMST10BTS	10 mm Single Use Sterile Blunt Tip Trocar Built-in Converters, Sizes 5 mm and 7 mm/8 mm Short Length	5
OMST10BT	10 mm Single Use Sterile Blunt Tip Trocar Built-in Converters, Sizes 5 mm and 7 mm/8 mm Standard Length	5
OMST12BT	12 mm Single Use Sterile Blunt Tip Trocar Built-in Converters, Sizes 5 mm and 10 mm Standard Length	5



VersaOne™ Blunt Trocars


Single use converterless 5 mm to 12 mm trocar system

Code	Description	Quantity
BPT12STS	VersaOne™ Blunt Trocar 5 mm - 12 mm With Threaded Anchoring Device	3

Trocars

VersaStep™ Plus Single Use Radially Expandable Access System

Includes Sleeves, Seals, and Needles



Code	Description	Quantity
VS070000	VersaStep™ Short Insufflation/Access Needle, 14 Gauge	12
S100000	Step™ Insufflation/Access Needle, 14 Gauge	10
VS150000	VersaStep™ Long Insufflation/Access Needle, 14 Gauge	12
VS100700	VersaStep™ Short Radially Expandable Sleeve	6
VS101000	VersaStep™ Radially Expandable Sleeve	6
VS101500	VersaStep™ Long Radially Expandable Sleeve	6
VSR100005	VersaStep™ Reposable 5 mm VersaSeal™ Seal with Radially Expandable Sleeve	6
VSR100812P	VersaStep™ Plus Reposable 5 mm-12 mm VersaSeal™ Plus Seal with Radially Expandable Sleeve	6

Trocars

VersaStep™ Plus Single Use Radially Expandable Access System

Includes Cannula and Dilator with Radially Expandable Sleeve




Code	Description	Quantity
VS100705	VersaStep™ Short 5 mm Cannula and Dilator with Radially Expandable Sleeve	3
VS101005	VersaStep™ Standard 5 mm Cannula and Dilator with Radially Expandable Sleeve	3
VS101505	VersaStep™ Long 5 mm Cannula and Dilator with Radially Expandable Sleeve	3
VS100711P	VersaStep™ Plus Short 11 mm Cannula and Dilator with Radially Expandable Sleeve	3
VS101011P	VersaStep™ Plus Standard 11 mm Cannula and Dilator with Radially Expandable Sleeve	3
VS100712P	VersaStep™ Plus Short 12 mm Cannula and Dilator with Radially Expandable Sleeve	3
VS101012P	VersaStep™ Plus Standard 12 mm Cannula and Dilator with Radially Expandable Sleeve	3
VS101512P	VersaStep™ Plus Long 12 mm Cannula and Dilator with Radially Expandable Sleeve	3
VS101015P	VersaStep™ Plus Standard 15 mm Cannula and Dilator with Radially Expandable Sleeve	3

Step™ Trocars



Step™ Single Use Radially Expandable Access System

Code	Description	Quantity
S100705	Step™ Short 5 mm Insufflation / Access Needle	10
S101005	Step™ 5 mm Insufflation / Access Needle	10
S100710	Step™ Short 10 mm Insufflation / Access Needle	10
S101005	Step™ 10 mm Insufflation / Access Needle	10
S100712	Step™ Short 12 mm Insufflation / Access Needle	10
S101012	Step™ 12 mm Insufflation / Access Needle	10



Mini Step™ Single Use Radially Expandable Access System

Code	Description	Quantity
MS100703	Mini Step™ Short 2 mm/3 mm Cannula and Dilator with Radially Expandable Sleeve	10
MS101003	Mini Step™ 2 mm/3 mm Cannula and Dilator with Radially Expandable Sleeve	10
MS100705	Mini Step™ Short 5 mm Cannula and Dilator with Radially Expandable Sleeve and 5 mm-2 mm/3 mm Reducer	10
MS101005	Mini Step™ 5 mm Cannula and Dilator with Radially Expandable Sleeve and 5 mm-2 mm/3 mm Reducer	10
MS100708	Mini Step™ Short 7 mm/8 mm Cannula and Dilator with Radially Expandable Sleeve and 7 mm/8 mm-5 mm Reducer	10
MS101008	Mini Step™ 7 mm/8 mm Cannula and Dilator with Radially Expandable Sleeve and 7 mm/8 mm-5 mm Reducer	10

Insufflation Access Needle and Sleeves, Compatible with both Step™ and Mini Step™ Systems



Code	Description	Quantity
S110000	Step™ Short Insufflation/Access Needle, 14 Gauge	10
S100000	Step™ Insufflation/Access Needle, 14 Gauge	10
S100700	Step™ Short Radially Expandable Sleeve	10
S101000	Step™ Radially Expandable Sleeve	10

Thoracic access products

Thoracoport™

Single use thoracic trocars with non-conductive sleeve



Code	Description	Quantity
179305	Thoracoport™ single use trocar, 5.5 mm	12
179301	Thoracoport™ single use trocar, 10.5 mm	12
179303	Thoracoport™ single use trocar, 11.5 mm	12
179307	Thoracoport™ single use trocar, 15 mm	12
179308	Thoracoport™ single use trocar, Soft 5 mm	12
179309	Thoracoport™ single use trocar, Soft 12 mm	12
179310	Thoracoport™ single use trocar, Soft 15 mm	12

SILS™ port

SILS™ Port



Code	Description	Quantity
SILSPT5TA	(1) SILS™ 5 mm Flexible Port With insufflation / smoke evacuation tubing and stopcock	6
	(3) 5 mm Cannula Suitable for use with 5 mm instrumentation only	
	(1) 5 mm Blunt-Tipped Obturator	
SILSPT12TA	(1) SILS™ 12 mm Flexible Port With insufflation / smoke evacuation tubing and stopcock	6
	(3) 5 mm Cannula	
	(1) 5 mm Blunt-Tipped Obturator	
	(1) 5 mm-12 mm Cannula with VersaSeal™ Plus Self Adjusting Seal Suitable for use with instrumentation ranging from 5 mm - 12 mm	
	(1) 12 mm Blunt-Tipped Obturator	

Note: The SILS™ Port has been designed for use only with the SILS™ port components

Fascial closure system

VersaOne™ Fascial Closure System

Multifunctional tool that serves as a trocar and a fascial closure device



Code	Description	Quantity
ONB12STFCS	VersaOne™ Fascial Closure System Optical Trocar, 12 mm	6
B12STFCS	VersaOne™ Fascial Closure System Bladed Trocar, 12 mm	6
FCSPASS	VersaOne™ Fascial Closure System Reusable Suture Passer	1

Visualization System



Clearify™ Visualization System

Code	Description	Quantity
21-345	Clearify™ Visualization System Warming hub, microfiber cloths, and trocar wipe	10

Anti-Fog and Lubricant Solutions




Code	Description	Quantity
220-50	Fred™ Anti-Fog Kit	20
220-90	Endo-Lube™ Seal and Instrument Lubrication Kit	20

Accessories

Surgineedle™ Insufflation Needle

Single use pneumoperitoneum needles



Code	Description	Quantity
172015	Surgineedle™ 120 mm Insufflation Needle	12
172016	Surgineedle™ 150 mm Long Insufflation Needle	12

Surgisleeve™ wound protector

SurgiSleeve™ Wound Protector

Maximize wound exposure and help to protect the wound site from contamination, while providing a strong film barrier.




Order Code	Description	Box Quantity
WPXSM24	Wound Protector Extra Small Incision Size 2 - 4 cm	5
WPSM256	Wound Protector Small Incision Size 2.5 - 6 cm	5
WPLG914	Wound Protector Large Incision Size 9 - 14 cm	5
WPLGR914	Wound Protector Large Incision Size 9 - 14 cm with Retraction Ring	5
WPXLGR1117†	Wound Protector Extra Large Incision Size 11 - 17 cm with Retraction Ring†	5

Laparoscopic hand instruments



Accessories



Laparoscopic Grasping Instruments

Single use

Code	Description	Quantity
173030	Endo Grasp™ 5 mm Single Use Instrument with On/Off Ratchet Switch - 31 cm Length	6
174317	Endo Clinch™ II 5 mm Single Use Instrument with On/Off Ratchet Switch - 31 cm Length	6
174001	Endo Babcock™ 10 mm Single Use Clamp - 31 cm Length	6
173051	Endo Lung™ 12 mm Single Use Clamp - 19 cm Length	6
174213	Endo Mini Dissect™ 5 mm Dissector with Roticulator™ Technology	6
176645	Endo Dissect™ 5 mm Dissector	6

Retractors

Laparoscopic Retractors

Single use



Code	Description	Quantity
176613	Endo Retract™ 10 mm Single Use Instrument - 31 cm Length, 3 Finger Articulating	6
176647	Endo Retract™ II 10 mm Single Use Instrument - 32 cm Length, 5 Finger Articulating	6
174311	Endo Retract™ Maxi 10 mm Single Use Instrument - 32 cm Length	6
173046	Endo Paddle Retract™ 12 mm Single Use Instrument - 46 cm Length	3
174209	Endo Mini-Retract™ 5 mm Single Use Instrument - 31 cm Length	6
173019	Endo Peanut™ 5 mm Single Use Device - 45 cm Length (3 units per pouch - 4 pouches per box)	3

Specimen retrieval products



Single Use Specimen Retrieval Products

Code	Description	Quantity
173050G	Endo Catch™ Gold 10 mm 2.5" x 6" Specimen Pouch - 34.5 cm	6
173049	Endo Catch™ II 15 mm 5" x 9" Specimen Pouch - 34.5 cm	3



©2023 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.™ Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 10/2023 - US-WC-2300044 - [WF# 630375]



[medtronic.com/animalhealth](https://www.medtronic.com/animalhealth)

Medtronic