

Medtronic

VersaOne™ optical trocar with fixation balloon cannula

Trocar stability – from access to closure.^{1,†,‡}

Every surgery has its challenges – but a trocar that doesn't stay in place shouldn't be one of them.

Introducing the VersaOne™ optical trocar with fixation balloon cannula – a new access solution with a number of benefits:


- Consistent, reliable access^{1,2,†,‡,††}
- Reduced migration risk^{1,‡,Ω}
- Adaptability when you want it^{3,‡,††}

Compatible with 5 mm standard length VersaOne™ optical obturators: passive scope lock design to hold laparoscope in place while in the obturator

Low profile: provides greater flexibility when placing multiple ports

Antimigration collar: reduces the risk of internal trocar movement^{1,‡,Ω}

Stabilizer balloon: provides superior abdominal fixation^{1,‡,Ω}



Insufflation port

Balloon inflation port: where the syringe connects to inflate balloon

Expanding the
VersaOne™
access family

†Based on in vivo fixation testing 4 hours after penetration with an in vivo porcine model. ‡Preclinical results may not correlate with clinical performance in humans. §VersaOne™ fixation balloon cannula portfolio includes 5 mm short, standard, and long devices. ΩCompared to competitive Endopath XCEL™ ribbed, Kii™ Z thread trocars and Medtronic VersaOne™ optical trocar with fixation cannula. ††17 out of 17 surgeons surveyed after use agreed.

Your instruments need to move, your trocar doesn't.

Supports
instrument
exchanges^{1,2,†,‡}

Improved
fixation^{1,‡,§}

Minimal depth
penetration^{1,2,‡,Ω}

Ordering information

For more information on the VersaOne™ optical trocar with fixation balloon cannula, call 800-772-8772 or visit us at [medtronic.com/versaone-balloon](https://www.medtronic.com/versaone-balloon)



Order code	Diameter	Length
ONB5SHB	5 mm	70 mm short optical trocar with fixation balloon cannula
ONB5STB	5 mm	100 mm standard optical trocar with fixation balloon cannula
ONB5LGB	5 mm	150 mm long optical trocar with fixation balloon cannula
UNVCA5STB	5 mm	100 mm standard optical trocar with fixation balloon cannula only

Quantity per box: 6



©2024 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 03/2024 - US-AI-2400008 - [WF#10893172]

[medtronic.com/covidien](https://www.medtronic.com/covidien)

1. Based on internal report #RE00453137, VersaOne™ optical trocar with fixation balloon cannula 5 mm in vivo claims evaluation lab report. June 26, 2023.
2. Based on internal report #RE00471794, VersaOne™ optical trocar with fixation balloon cannula 5 mm surgeon use questionnaire. July 26, 2023.
3. Based on internal report #RE00322029, 5 mm VersaOne™ advanced fixation trocar (Project Benri phase 1) validation/summative report. Sept. 21, 2021.

Medtronic