

Medtronic

Evolut™ FX+ TAVR System

Product specifications

Bioprosthesis




Evolut FX+ 23 mm Valve

Model number:
EVFXPLUS-23

Size: 23 mm

Aortic annulus
diameter:
17†/18-20 mm




Evolut FX+ 26 mm Valve

Model number:
EVFXPLUS-26

Size: 26 mm

Aortic annulus
diameter:
20-23 mm



Evolut FX+ 29 mm Valve

Model number:
EVFXPLUS-29

Size: 29 mm

Aortic annulus
diameter:
23-26 mm



Evolut FX+ 34 mm Valve

Model number:
EVFXPLUS-34

Size: 34 mm

Aortic annulus
diameter:
26-30 mm


Loading system



Evolut FX+ Loading System

Model number:
L-EVOLUTFX-2329

23, 26, and 29 mm TAV



Evolut FX+ Loading System

Model number:
L-EVOLUTFX-34

34 mm TAV


†Measurement for TAV-in-SAV only.

All items above are single-use. Refer to Instructions for Use for a complete list of warnings, precautions, and contraindications.

Delivery catheter system


Evolut FX+ delivery catheter system with InLine sheath 14 Fr equivalent (6.0 mm capsule outer diameter)

Model number:
D-EVOLUTFX-2329
23, 26, and 29 mm TAV



Evolut FX+ delivery catheter system with InLine sheath 18 Fr equivalent (7.33 mm capsule outer diameter)

Model number:
D-EVOLUTFX-34
34 mm TAV



Medtronic

710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA

Toll-free: 800.633.8766
Tel: +1.763.514.4000

[medtronic.com](https://www.medtronic.com)

LifeLine
CardioVascular Technical Support
Toll-free: 877.526.7890
Tel: +1.763.526.7890
rs.structuralheart@medtronic.com

©2024 Medtronic. All rights reserved.
Medtronic, Medtronic logo, and Engineering
the extraordinary are trademarks of Medtronic.
All other brands are trademarks of a Medtronic
company.

UC202402378 EN
03/2024