








Medtronic

Capnostream™ 20p bedside patient monitor

Power, connectivity
and mounting
accessories



Power and connectivity accessories

Product		Part number
	U.S. power cord (3 meters)	014256
	Battery pack	BP01000
	Nurse call system cable Supplied unterminated	011149
	Power cable strain relief	CS08807
	Fuse, 3.15 amp/250 volt (2/package)	010543
	Paper (6 rolls/package)	010516

Mounting accessories

Product	Part number	
<p>Wireless bridge adapter (Bernoulli-compatible)</p>	011892	Supports all Nellcor™ OxiNet™ technology bridges: <ul style="list-style-type: none"> - Moxa TCC-80 (DB 9) RS-232 to 422 Adapter - Moxa NPort RS 232 to wired ethernet adapter - Nellcor™ brand version of wireless bridge Reuses Capnostream™ 20p patient monitor feet for table-top use. Handle designed for monitor transport and use with GCX 5" Slide-in plunger mounting plate. Hook and loop straps offer storage of the SpO ₂ cabling and power cords.
<p>Quick release pole clamp kit</p>	011782	Clamp fits pole mount (0.75"-1.5"/19 mm-38 mm) or rail mount (0.39"/10 mm x 0.98"/25 mm). Note: Never use a temporary mount when a falling Capnostream™ 20p monitor could cause patient injury.
<p>VESA standard mounting adapter plate</p>	010713	For mounting Capnostream™ 20p monitor on any VESA standard 75 x 75 mounting system as required by various industry mounting solutions.
<p>Roll stand, handle and basket</p>	10002385	Features tilt/swivel adjustments. Post height is 38"/96.5 cm, overall height with casters, floor to mounting plate is 43"/109.2 cm. Mounting plate sold separately (010713).
<p>CPU channel mount for Nuvon VEGA™ system</p>	GXC Part #RS-0009-40 [†]	When the Capnostream™ 20p patient monitor is connected to the Nuvon VEGA™ System, the IDM-MG 3000 network appliance can be mounted on the post of a GCX Corporation roll stand.
<p>Wall channel, vertical</p>	10003482	Standard channel. Channel permits monitor to be set anywhere in the channel.
<p>Wall mount pivot arm</p>	10003484	12" pivot arm. Pivots at the slide for lateral adjustability. Wall Channel sold separately (10003482).

[†]Please visit gcx.com for product details and availability. Capnostream™ 20p bedside monitor accessories are designed to accommodate several mounting options. Additional mounting options are available through GCX Corporation.