



Endo Stitch™ Single-Stitch Loading Units

PRODUCT INFORMATION

Code	Description	Qty Per Case
170001	Sofsilk™ Ctd, Brd Silk	12
170002	Sofsilk™ Ctd, Brd Silk	12
170003	Sofsilk™ Ctd, Brd Silk	12
170004	Sofsilk™ Ctd, Brd Silk	12
170010	Bralon™ Ctd, Brd Nylon	12
170011	Bralon™ Ctd, Brd Nylon	12
170012	Bralon™ Ctd, Brd Nylon	12
170013	Bralon™ Ctd, Brd Nylon	12
170040	Surgidac™ Ctd, Brd Polyester non absorbable	3
170041	Surgidac™ Ctd, Brd Polyester non absorbable	3
170043	Surgidac™ Ctd, Brd Polyester non absorbable	3
170044	Surgidac™ Ctd, Brd Polyester non absorbable	3
170050	Polysorb™ Ctd, Brd Syn. Abs.	12
170051	Polysorb™ Ctd, Brd Syn. Abs.	12
170052	Polysorb™ Ctd, Brd Syn. Abs.	12
170053	Polysorb™ Ctd, Brd Syn. Abs.	12
170054	Polysorb™ Ctd, Brd Syn. Abs.	12
170055	Polysorb™ Ctd, Brd Syn. Abs.	12
170070	Polysorb™ Ctd, Brd Syn. Abs.	12
170071	Polysorb™ Ctd, Brd Syn. Abs.	12
170072	Polysorb™ Ctd, Brd Syn. Abs.	12
170073	Polysorb™ Ctd, Brd Syn. Abs.	12
170090	Polysorb™ Ctd, Brd Syn. Abs.	12
170092	Polysorb™ Ctd, Brd Syn. Abs.	12
170094	Polysorb™ Ctd, Brd Syn. Abs.	12
173021	Surgidac™ Ctd, Brd Polyester	12
173023	Surgidac™ Ctd, Brd Polyester	12
173024	Surgidac™ Ctd, Brd Polyester	12
173026	Surgidac™ Ctd, Brd Polyester	12

MATERIALS OF MANUFACTURE

170001	ABS GE White, 303 Stainless Steel, SilkBlack 1/0, Nylon
170002	ABS GE White, 303 Stainless Steel, SilkBlack 1/0, Nylon
170003	ABS GE White, 303 Stainless Steel, SilkBlack 1/0, Nylon
170004	ABS GE White, 303 Stainless Steel, SilkBlack 1/0, Nylon
170010	ABS GE White, 303 Stainless Steel, SilkBlack 1/0, Nylon braided Black size 0, 2/0
170011	ABS GE White, 303 Stainless Steel, SilkBlack 1/0, Nylon braided Black size 2/0
170012	ABS GE White, 303 Stainless Steel, SilkBlack 1/0, Nylon braided Black size 0, 2/0
170013	ABS GE White, 303 Stainless Steel, SilkBlack 1/0, Nylon braided Black size 2/0
170040	ABS GE White, 303 Stainless Steel, Brd cut Green size 0, 2/0 uncoated
170041	ABS GE White, 303 Stainless Steel, Brd cut Green size 2/0 uncoated
170043	ABS GE White, 303 Stainless Steel, Brd cut Green size 0, 2/0 uncoated
170044	ABS GE White, 303 Stainless Steel, Brd cut Green size 2/0 uncoated
170050	ABS GE, 303 Stainless Steel, BSA Violet L23V Size 0 52/48 coated, L6 Drawn Yarn Violet, BSA Braider Bobbin Core Bobbin, L-Lactide Reagent grade L6 Lactomer, Nylon LNP Lubricomp, Monadnock 12PT Suture Board
170051	ABS GE, 303 Stainless Steel, BSA Violet L23V Size 2/0 52/48 coated, L6 Drawn Yarn Violet, BSA Braider Bobbin Core Bobbin, L-Lactide Reagent grade L6 Lactomer, Nylon LNP Lubricomp, Monadnock 12PT Suture Board
170052	ABS GE, 303 Stainless Steel, BSA Violet L23V Size 0 52/48 coated, L6 Drawn Yarn Violet, BSA Braider Bobbin Core Bobbin, L-Lactide Reagent grade L6 Lactomer, Nylon LNP Lubricomp, Monadnock 12PT Suture Board
170053	ABS GE, 303 Stainless Steel, BSA Violet L23V Size 2/0 52/48 coated, L6 Drawn Yarn Violet, BSA Braider Bobbin Core Bobbin, L-Lactide Reagent grade L6 Lactomer, Nylon LNP Lubricomp, Monadnock 12PT Suture Board
170054	ABS GE, 303 Stainless Steel, BSA Undyed L23V Size 0, 52/48 coated, L6 Drawn Yarn Natural, BSA Braider Bobbin Undyed Core Bobbin, L-Lactide Reagent grade L6 Lactomer, Nylon LNP Lubricomp, Monadnock 12PT Suture Board
170055	ABS GE, 303 Stainless Steel, BSA Undyed L23V Size 2/0, 52/48 coated, L6 Drawn Yarn Natural, BSA Braider Bobbin Undyed Core Bobbin, L-Lactide Reagent grade L6 Lactomer, Nylon LNP Lubricomp, Monadnock 12PT Suture Board
170070	ABS GE, 303 Stainless Steel, BSA Violet L23V Size 2/0 52/48 coated, L6 Drawn Yarn Violet, BSA Braider Bobbin Core Bobbin, L-Lactide Reagent grade L6 Lactomer, Nylon LNP Lubricomp, Monadnock 12PT Suture Board
170071	ABS GE, 303 Stainless Steel, BSA Violet L23V Size 2/0 52/48 coated, L6 Drawn Yarn Violet, BSA Braider Bobbin Core Bobbin, L-Lactide Reagent grade L6 Lactomer, Nylon LNP Lubricomp, Monadnock 12PT Suture Board
170072	ABS GE, 303 Stainless Steel, BSA Undyed L23V Size 3/0, 52/48 coated, L6 Drawn Yarn Natural, BSA Braider Bobbin Undyed Core Bobbin, L-Lactide Reagent grade L6 Lactomer, Nylon LNP Lubricomp, Monadnock 12PT Suture Board
170073	ABS GE, 303 Stainless Steel, BSA Undyed L23V Size 3/0, 52/48 coated, L6 Drawn Yarn Natural, BSA Braider Bobbin Undyed Core Bobbin, L-Lactide Reagent grade L6 Lactomer, Nylon LNP Lubricomp, Monadnock 12PT Suture Board
170090	ABS GE, 303 Stainless Steel, BSA Violet L23V Size 4/0 52/48 coated, L6 Drawn Yarn Violet, BSA Braider Bobbin, L-Lactide Reagent grade L6 Lactomer, Nylon LNP Lubricomp, Monadnock 12PT Suture Board
170092	ABS GE, 303 Stainless Steel, BSA Violet L23V Size 4/0 52/48 coated, L6 Drawn Yarn Violet, BSA Braider Bobbin L-Lactide Reagent grade L6 Lactomer, Nylon LNP Lubricomp, Monadnock 12PT Suture Board
170094	ABS GE, 303 Stainless Steel, BSA Undyed L23V Size 4/0, 52/48 coated, L6 Drawn Yarn Natural, BSA Braider Bobbin, L-Lactide Reagent grade L6 Lactomer, Nylon LNP Lubricomp, Monadnock 12PT Suture Board
173021	ABS GE, 303 Stainless Steel, Green Size 2/0 uncoated, Nylon LNP Lubricomp
173023	ABS GE, 303 Stainless Steel, Green Size 2/0 uncoated, Nylon LNP Lubricomp
173024	ABS GE, 303 Stainless Steel, Green Size 2/0 uncoated, BRD cut Green Size 0, Nylon LNP Lubricomp
173026	ABS GE, 303 Stainless Steel, Green size 2/0 uncoated, BRD Cut Green Size 0, Nylon LNP Lubricomp

SHIPPING INFORMATION

	Inner Box Size (mm)
170001	L 120 × W 59 × H 145
170002	L 120 × W 59 × H 145
170003	L 120 × W 59 × H 145
170004	L 120 × W 59 × H 145
170010	L 120 × W 59 × H 145
170011	L 120 × W 59 × H 145
170012	L 120 × W 59 × H 145
170013	L 120 × W 59 × H 145
170040	L 121 × W 69 × H 144
170041	L 121 × W 69 × H 144
170043	L 121 × W 69 × H 144
170044	L 121 × W 69 × H 144
170050	L 120 × W 59 × H 145
170051	L 120 × W 59 × H 145
170052	L 120 × W 59 × H 145
170053	L 120 × W 59 × H 145
170054	L 120 × W 59 × H 145
170055	L 120 × W 59 × H 145
170070	L 120 × W 59 × H 145
170071	L 120 × W 59 × H 145
170072	L 120 × W 59 × H 145
170073	L 120 × W 59 × H 145
170090	L 120 × W 59 × H 145
170092	L 120 × W 59 × H 145
170094	L 120 × W 59 × H 145
173021	L 120 × W 59 × H 145
173023	L 120 × W 59 × H 145
173024	L 120 × W 59 × H 145
173026	L 120 × W 59 × H 145

STERILIZATION METHOD

	Method
170001	Ethylene Oxide/Radiation
170002	Ethylene Oxide/Radiation
170003	Ethylene Oxide/Radiation
170004	Ethylene Oxide/Radiation
170010	Ethylene Oxide/Radiation
170011	Ethylene Oxide/Radiation
170012	Ethylene Oxide/Radiation
170013	Ethylene Oxide/Radiation
170040	Ethylene Oxide/Radiation
170041	Ethylene Oxide/Radiation
170043	Ethylene Oxide/Radiation
170044	Ethylene Oxide/Radiation
170050	Ethylene Oxide
170051	Ethylene Oxide
170052	Ethylene Oxide
170053	Ethylene Oxide
170054	Ethylene Oxide
170055	Ethylene Oxide
170070	Ethylene Oxide
170071	Ethylene Oxide
170072	Ethylene Oxide
170073	Ethylene Oxide
170090	Ethylene Oxide
170092	Ethylene Oxide
170094	Ethylene Oxide
173021	Ethylene Oxide/Radiation
173023	Ethylene Oxide/Radiation
173024	Ethylene Oxide/Radiation
173026	Ethylene Oxide/Radiation

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Shelf Life	5 years
CE Marked	Class III According to Directive 93/42/EEC on Medical Devices (MDD), Annex II



IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 17-weu-loading-units-tech-sheet-1515729