



Surgitie™

single use ligating loops with delivery system

PRODUCT INFORMATION

Code	Description	Qty Per Case
ED-20-L	Surgitie™ Ligating Loop with Surgidac™ Coated, Braided, Polyester Suture	12
ED-21-L	Surgitie™ Ligating Loop with Surgidac™ Coated, Braided, Polyester Suture	12
EL-20-L	Surgitie™ Ligating Loop with Delivery System with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	6
EL-21-L	Surgitie™ Ligating Loop with Delivery System with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	6
EL-22-LN	Surgitie™ Ligating Loop with Delivery System with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	6
EL-23-LN	Surgitie™ Ligating Loop with Delivery System with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	6

MATERIALS OF MANUFACTURE

Code	Description
ED-20-L	Nylon, clear polycarbonate, thermoplastic elastomer, stainless steel and PET suture
ED-21-L	Nylon, clear polycarbonate, thermoplastic elastomer, stainless steel and PET suture
EL-20-L	Nylon, clear polycarbonate, thermoplastic elastomer, stainless steel and polysorb suture
EL-21-L	Nylon, clear polycarbonate, thermoplastic elastomer, stainless steel and polysorb suture
EL-22-LN	Nylon, clear polycarbonate, thermoplastic elastomer, stainless steel and polysorb suture
EL-23-LN	Nylon, clear polycarbonate, thermoplastic elastomer, stainless steel and polysorb suture

SHIPPING INFORMATION

Code	Inner Box Size (mm)
ED-20-L	L 197 × W 152 × H 470
ED-21-L	L 197 × W 152 × H 470
EL-20-L	L 156 × W 152 × H 470
EL-21-L	L 156 × W 152 × H 470
EL-22-LN	L 156 × W 152 × H 470
EL-23-LN	L 156 × W 152 × H 470

STERILIZATION METHOD

Code	
ED-20-L	Radiation
ED-21-L	Radiation
EL-20-L	Ethylene Oxide
EL-21-L	Ethylene Oxide
EL-22-LN	Ethylene Oxide
EL-23-LN	Ethylene Oxide

ADDITIONAL INFORMATION

Latex Free	✓
DEHP	✓
Shelf Life	5 years
CE Marked	Class III Directive 93/42/EEC on Medical Devices (MDD), Annex II

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 17-weu-surgitie-tech-sheet-1515966

Medtronic
Further, Together