



Thoracoport™

Single Use Trocars with Non conductive Sleeves

PRODUCT INFORMATION

Code	Description	Qty Per Case
179301	Thoracoport™ 10.5 mm for instrument up to 11 mm, Single Use Trocar, non Conductive Sleeve	12 per box
179303	Thoracoport™ 11.5 mm for instrument up to 12 mm, Single Use Trocar, non Conductive Sleeve	12 per box
179305	Thoracoport™ 5.5 mm for instrument up to 6 mm, Single Use Trocar, non Conductive Sleeve	12 per box
179307	Thoracoport™ 15 mm for instrument up to 16,5 mm, Single Use Trocar, non Conductive Sleeve	12 per box

MATERIALS OF MANUFACTURE

Code	Description
179301	ABS, GE White, 302, 304, and 430 Stainless Steel, Polycarbonate, GE Lexan, PVC Alpha, Dow Polycarbonate Megarad, Polypropylene white, Colorite Unichem PVC compound, Elastomer Santoprene and Fraton
179303	ABS GE, Cicolac, Acetal %15 glass filled PTFE, C-Flex 30 Durometer
179305	ABS GE, Cicolac, Acetal %15 glass filled PTFE, C-Flex 30 Durometer
179307	ABS GE, Cicolac, Acetal %15 glass filled PTFE, C-Flex 30 Durometer

SHIPPING INFORMATION

	Shipped Case Dimensions (mm)
Inner box	L 193 × W 152 × H 195

ADDITIONAL INFORMATION

Latex Free	✓
DEHP free	✓
Sterilization method	Ethylene Oxide
Shelf Life	5 years
CE Marked	Class IIa According to Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 17-weu-thoracoport-trocars-tech-sheet-en-1515172

Medtronic
Further, Together