



Endo Babcock™

laparoscopic hand instrument 10 mm

PRODUCT INFORMATION

Code	Description	Qty Per Case
174001	Endo Babcock™ Laparoscopic Hand Instruments 10 mm, Single Use Clamp – 31 cm Length	1

MATERIALS OF MANUFACTURE

302, 303, 304, 316L AND 410 Stainless Steel, Polycarbonate 20% Glass, Polycarbonate 20% Glass Filled, Polycarbonate 30% Carbon, Polycarbonate GE, Polycarbonate 10% Glass Filled, PVC Grade 6k10, Style Vw-1, Acetal Dupont Delrin, LNP, Rubber/Silicon Durometer, Aluminum 6061-T6

SHIPPING INFORMATION

	Inner Box Size (mm)	Shipper Box Size (mm)
174001	L 162 × W 27 × H 560	L 181 × W 168 × H 572

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Ethylene Oxide
Shelf Life	5 years
CE Marked	Class IIa According to Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 18-weu-endo-babcock-tech-sheet-2341954

Medtronic
Further, Together