



# Endo Grasp™

## 5 mm Single Use Instrument

### PRODUCT INFORMATION

Code	Description	Qty Per Case
173030	Endo Grasp™ 5 mm, Single Use Instrument with On/Off Ratchet Switch – 31 cm Length	6

### MATERIALS OF MANUFACTURE

301, 302, 303, Stainless Steel, 17-4 Stainless Steel, 304 Sintered Stainless Steel, Polycarbonate 20% Glass, 10% Glass Filled Polycarbonate Lnp, Polycarbonate 20% Glass Filled, Pvc Grade 6k10, Style Vw-1, Aluminum, 2011-T3 Or 6061-T6, Rubber/Silicon Durometer, Polyurethane

### SHIPPING INFORMATION

	Inner Box Size (mm)
173030	L 610 × W 191 × H 70

### ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Radiation
Shelf Life	5 years
CE Marked	Class IIa Accordign to Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

**IMPORTANT:** Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 18-weu-endo-grasp-tech-sheet-2324090

**Medtronic**  
Further, Together