



# Endo Peanut™

## single use device

### 5 mm

#### PRODUCT INFORMATION

Code	Description	Qty Per Case
173019	Endo Peanut™ 5 mm, Single Use Device – 45 cm Length	3 units per pouch – 4 pouches per box

#### MATERIALS OF MANUFACTURE

Polycarbonate 20% Glass Fiber, Polycarbonate 20% Glass Filled, Polycarbonate, Ge, Aluminum, 6061-T6, Pvc Grade 4k10 Std Fw:21, Elastomer Santoprene, Lnp Df

#### SHIPPING INFORMATION

	Inner Box Size (mm)
173019	L 83 × W 19 × H 489

#### ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Radiation
Shelf Life	5 years
CE Marked	Class IIa According to Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

**IMPORTANT:** Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.18-weu-endo-peanut-tech-sheet-2342877

**Medtronic**  
Further, Together