



Endo Retract™

10 mm laparoscopic retractor

PRODUCT INFORMATION

Code	Description	Qty Per Case
176613	Endo Retract™ 10 mm Single Use Instrument – 31 cm Length, 3 Finger	1
176647	Endo Retract™ II 10 mm Single Use Instrument – 32 cm Length, 5 Finger articulating	1

MATERIALS OF MANUFACTURE

176613	Polycarbonate, PVC, Acetal, Santoprene, Stainless Steel, Nylon
176647	301, 302, 303, 304, 316L and 17-7 P.H Stainless Steel, Polycarbonate 30% carbon, PVC Grade 6K10, LNP, Acetal 15% glass, Acetal 10% PTFE, Santoprene

SHIPPING INFORMATION

	Inner Box Size (mm)	Shipper Box Size (mm)
176613	L 109 × W 28 × H 559	L 183 × W 114 × H 564
176647	L 109 × W 41 × H 559	L 240 × W 133 × H 572

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Ethylene Oxide
Shelf Life	5 years
CE Marked	Class IIb Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 18-weu-endo-retract-10-mm-tech-sheet-2342616

Medtronic
Further, Together