



# EndoCatch™ Gold specimen retrieval pouch

## PRODUCT INFORMATION

Code	Description	Qty Per Case
173050G	Endo Catch™ Gold 10 mm 2.5" × 6" Specimen Pouch, Volume capacity 190 cc, 29,5 cm Shaft Length	6

## MATERIALS OF MANUFACTURE

303, 304 Stainless Steel, Polycarbonate 20% Glass, Ptfе, Polycarbonate Ge Lexan, Rubber/Silicon Durometer, Polypropylene Red, Tinel Bb, Rnf 150, Polyurethan

## SHIPPING INFORMATION

	Inner Box Size (mm)
173050G	L 152 × W 111 × H 562

## ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Ethylene Oxide
Shelf Life	3 years
CE Marked	Class IIa According to Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

**IMPORTANT:** Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 18-weu-endocatch-gold-tech-sheet-2652114