



OPERATING ROOM **ESSENTIALS**

THE COMPLETE
SOLUTION
YOU CAN TRUST



Surgical Instruments
Catalogue

September 2020

Medtronic
Further, Together

WHY MEDTRONIC?

Medtronic has an extensive heritage in partnerships with Public Hospitals and Private Groups. Delivering value-based healthcare, clinical expertise and innovation through products and services that are focused on patient care and clinical outcomes.

OR ESSENTIALS


We pride ourselves on being a brand you can trust. With our proven track record of **quality products and services**, you can feel reassured that Medtronic's OR Essentials will help you to continue providing **the very best patient care** from the first mark to the last stitch. Let's embrace the future together.



Laparoscopic Access



Surgical Instruments



Wound Closure

INTRODUCTION

A comprehensive, versatile line of hand instruments and ligation devices built on more than 40 years of investment in the science and engineering of minimally invasive surgical products, hand instruments and ligation devices from Medtronic deliver the choice and flexibility to meet your needs.

Surgical Clips

Endoscopic clip appliers include ligation devices as Endo Clip™ and Lapro-Clip™ instruments.

The Lapro-Clip™ ligation system offers a reusable single clip applier that delivers absorbable clips in laparoscopic procedures.

Open ligation devices include a wide range of Premium Surgiclip™ clip appliers.

Laparoscopic hand instruments

Wide range of instruments to cope with a large variety of requirements.

Hand Instruments include: Endo Dissect™, Endo Shears™, Endo Mini Shears™, Endo Grasp™, Endo Clinch™, Endo Babcock™.

Surgical retractors

Laparoscopic retractors include, Endo Paddle Retract™, Endo Mini Retract™ and several Endo Retract™ instruments.

Suction and irrigation devices instrumentation includes Surgiwand™ II devices for endoscopic surgery – and Vital Vue™ instruments for open procedures.

Specimen retrieval products include EndoCatch™ and ReliaCatch™ specimen systems and Endo Bag™ specimen retrieval system that isolate and contain specimens during laparoscopic surgery.

Automatic suturing devices

Endoscopic suturing devices include Endo Stitch™ and Endo Close™ instruments.

Our endoscopic suturing accessories include Endo Stitch™ and V-Loc™ reloads, as well as Surgitie™ and Surgiwip™ suturing products.

CONTENTS

CLIP APPLIERS	 Endoscopic Clip Appliers.....	06
	 Open Clip Appliers.....	08
LAPAROSCOPIC HAND INSTRUMENTS	 Dissectors, Graspers, Shears.....	10
	 Retractors.....	12
	 Suction and Irrigation.....	13
	 Specimen Retrieval Products.....	15
AUTOMATIC SUTURING DEVICES	 Endo Stitch™ Suturing Devices.....	17
	 Suturing accesories.....	20
	 Order Information.....	22

CLIP APPLIERS

Endo Clip™
clip applicator



Surgiclip™
clip applicator

ENDOSCOPIC CLIP APPLIERS

Clip Logic™ Technology

Endo Clip™ III 5 mm

Single Use Clip Applier with Clip Logic™ Technology



ORDER CODE	DESCRIPTION	BOX QTY
176630	Endo Clip™ III 5 mm 16 M/L Titanium Clips, 33 cm Shaft Length	6
176630B	Endo Clip™ III 5 mm Without clip display (no battery inside)	6

Endo Clip™ II 10 mm Pistol Grip

Single Use Clip Applier with Super Interlock™ Titanium Clips



ORDER CODE	DESCRIPTION	BOX QTY
176657	Endo Clip™ II 10 mm Pistol Grip 20 M/L Super Interlock™ Titanium Clips, 29 cm Shaft Length	6

Endo Clip™ 10 mm Pistol Grip

Single Use Clip Applier with Super Interlock™ Titanium Clips



ORDER CODE	DESCRIPTION	BOX QTY
176625	Endo Clip™ L 10 mm Pistol Grip 15 Large Super Interlock™ Titanium Clips, 33 cm Shaft Length	6

Lapro-Clip™ Clip Applier

Reusable stainless steel laparoscopic clip applier with absorbable Lapro-Clip™ clips



ORDER CODE	DESCRIPTION	BOX QTY
8886848700	Lapro-Clip™ 26,5 cm Shaft Length	1
8886848701	Lapro-Clip™ Long 33 cm Shaft Length	1

Lapro-Clip™ Absorbable Loading Units

Disposable loading units for use with the Lapro-Clip™ clip applicator code 8886848700. The 2 piece absorbable clip is made of an external polyglycolic acid body that runs over a flexible internal polyglactone part.



ORDER CODE	DESCRIPTION	BOX QTY
8886848808	Lapro-Clip™ Clip Disposable loading unit with 6 absorbable Medium/Large (8 mm) clips per sterile box.	6
8886848508	Lapro-Clip™ Multi Loading Unit Automatic disposable loading unit with 8 absorbable clips per sterile box.	5
8886848812	Lapro-Clip™ Clip Disposable loading unit with 6 absorbable Large (12 mm) clips per sterile box.	6
8886848882	Lapro-Clip™ Clip Disposable loading unit with 2 absorbable Medium-Large (8 mm) clips per sterile box.	10
8886848813	Lapro-Clip™ Clip Disposable loading unit with 2 absorbable Large (12 mm) clips per sterile box.	10

OPEN CLIP APPLIERS

Premium Surgiclip™ III

Single Use Automatic Clip Appliers with Super Interlock™ Titanium Clips



ORDER CODE	DESCRIPTION	BOX QTY
133650	Premium Surgiclip™ III S-9.0 20 Small Titanium Clips	6

Premium Surgiclip™ II

Single Use Automatic Clip Appliers with Super Interlock™ Titanium Clips



ORDER CODE	DESCRIPTION	BOX QTY
134051	Premium Surgiclip™ II M-9.75 20 Medium Titanium Clips	6
134053	Premium Surgiclip™ II M-11.5 30 Medium Titanium Clips	6

Premium Surgiclip™

Single Use Automatic Clip Appliers with Super Interlock™ Titanium Clips



ORDER CODE	DESCRIPTION	BOX QTY
134048	Premium Surgiclip™ L-13.0 15 Large Titanium Clips	6

LAPAROSCOPIC HAND INSTRUMENTS






DISSECTORS, GRASPERS, SHEARS

Roticator™ 5 mm Laparoscopic Hand Instruments

Single Use – 31 cm Length, 80° articulation








ORDER CODE	COLOR CODE	DESCRIPTION	BOX QTY
174213		Roticator™ Endo Dissect™ 5 mm Single Use Instrument with Monopolar Cautery	6
174309		Roticator™ Endo Mini-Shears™ 5 mm Single Use Instrument with Monopolar Cautery	6
174233		Roticator™ Endo Grasp™ 5 mm Single Use Instrument with Spin Lock™	6

Laparoscopic Hand Instruments 5 mm

Single Use – 31 cm Length



ORDER CODE	COLOR CODE	DESCRIPTION	BOX QTY
176643		Endo Shears™ 5 mm Single Use Instrument with Monopolar Cautery	6
174301		Endo Mini-Shears™ 5 mm Single Use Instrument with Monopolar Cautery	6
176645		Endo Dissect™ 5 mm Single Use Instrument with Monopolar Cautery	6
173030		Endo Grasp™ 5 mm Single Use Instrument with On/Off Ratchet Switch – 31 cm Length	6
174317		Endo Clinch™ II 5 mm Single Use Instrument with On/Off Ratchet Switch – 31 cm Length	6

Laparoscopic Hand Instrument Long 5 mm

Single Use – 45 cm Length



ORDER CODE	COLOR CODE	DESCRIPTION	BOX QTY
174601	■	Endo Shears™ Long 5 mm Single Use Instrument with Monopolar Cautery	6










Laparoscopic Hand Instruments 10 mm



ORDER CODE	COLOR CODE	DESCRIPTION	BOX QTY
174001	■	Endo Babcock™ 10 mm Single Use Clamp – 31 cm Length	6

RETRACTORS

Laparoscopic Retractors

	ORDER CODE	COLOR CODE	DESCRIPTION	BOX QTY
	176613		Endo Retract™ 10 mm Single Use Instrument – 31 cm Length, 3 Finger	6
	176647		Endo Retract™ II 10 mm Single Use Instrument – 32 cm Length, 5 Finger articulating	6
	173046		Endo Paddle Retract™ 12 mm Single Use Instrument – 46 cm Length	3
	174209		Endo Mini-Retract™ 5 mm Single Use Instrument – 31 cm Length	6
	173019		Endo Peanut™ 5 mm Single Use Device – 45 cm Length (3 units per pouch – 4 pouches per box)	3

SUCTION AND IRRIGATION

Surgiwand™ II 5 mm

Single Use Suction/Irrigation Devices without Cautery,
34 cm Shaft Length



ORDER CODE	DESCRIPTION	BOX QTY
178081	Surgiwand™ II 5 mm Suction /Irrigation Device Without Tubing For Use with the Nitrogen Pump System	3



178083	Surgiwand™ II 5 mm Suction/Irrigation Device With Tubing For Gravity Feed System	3
--------	--	---

Surgiwand™ II 5 mm with Cautery

Single Use Suction/Irrigation Devices, 34 cm Shaft Length



ORDER CODE	DESCRIPTION	BOX QTY
178091	Surgiwand™ II 5 mm Suction/Irrigation Device With Cautery and L-Hook Tip without Tubing For use with the Nitrogen Pump System	3



178099	Surgiwand™ II 5 mm Suction/Irrigation Device With Cautery and Spatula Tip without Tubing For use with the Nitrogen Pump System	3
--------	---	---



178093	Surgiwand™ II 5 mm Suction/Irrigation Device With Cautery, L-Hook Tip and Tubing For Gravity Feed System	3
--------	---	---



178094	Surgiwand™ II 5 mm Suction/Irrigation Device With Cautery, Spatula Tip and Tubing For Gravity Feed System	3
--------	--	---

Vital Vue™ Surgical Instrument



ORDER CODE	DESCRIPTION	BOX QTY
8886827906	Bulb tip Used in many gynaecological and plastic surgery cases as well as other shallow procedures where a bulb tip is appropriate	6
8886828006	Extended bulb tip For lower abdominal and pelvic cases where extra reach is required and a bulb tip is appropriate	6
8886828106	Extended slender tip The tip of choice for surgeons who prefer precise suction with a slender, 1,35 cm profile tip	6
8886828206	Extended metal tip The 12 Fr., 5 cm extended metal tip is ideal where reach and precise aspiration is required, such as oral surgery and some vascular cases	6
8886828306	Extra long extended metal tip Designed specifically for orthopaedic cases, this 16 Fr., 8 cm metal tip is ideal where extra extension is necessary particularly in hip arthroplasty	6

Vital Vue™ Light Power Source



CODE	DESCRIPTION	REGION	BOX QTY
8886828421	Vital Vue™ light power source	UK	1
8886828422	Vital Vue™ light power source	Spain	1
8886828423	Vital Vue™ light power source	Germany	1
8886828424	Vital Vue™ light power source	Italy	1
8886828425	Vital Vue™ light power source	France	1

SPECIMEN RETRIEVAL PRODUCTS

Single Use Specimen Retrieval Products (Latex Free)



ORDER CODE	DESCRIPTION	BOX QTY
CATCH10	Reliacatch™ Specimen Retrieval Bag 10 mm 6.2 × 16.6 cm Volume capacity 275 mL	5
173050G	Endo Catch™ Gold 10 mm 5.8 × 15 cm Specimen Pouch, Volume capacity 190 cc,	6
173049	Endo Catch™ II 15 mm 9.7 × 23.6 cm Specimen Pouch, Volume capacity 1500 cc,	3

Single Use Specimen Retrieval Pouches



ORDER CODE	DESCRIPTION	BOX QTY
250-30	Endobag™ 7.6 × 15.2 cm Specimen Retrieval System – 20.9 cm Shaft Length	20
250-40	Endobag™ 12.7 × 20.3 cm Specimen Retrieval System – 20.9 cm Shaft Length	20

AUTOMATIC SUTURING DEVICES

Suturing Made Easy
At Every Angle



Endo Stitch™
suturing device

Endo Stitch™
triple-stitch reload

ENDO STITCH™ SUTURING DEVICES



Endo Stitch™ 10 mm Suturing Device

ORDER CODE	DESCRIPTION	BOX QTY
173016	Endo Stitch™ 10 mm Suturing Device, 34 cm Shaft Length	3

Endo Stitch™ Triple-Stitch Loading Unit

Single Use Loading Units for use with SILS™ Stitch and Endo Stitch™ 10 mm Suturing Devices. Each Loading Unit contains three strands of suture, each attached to a 9 mm needle.



ORDER CODE	SUTURE MATERIAL	USP SIZE	METRIC SIZE	LENGTH	BOX QTY
170043	Surgidac™ Ctd, Brd Polyester non absorbable	0	3.5	120 cm / 48"	3
170040	Surgidac™ Ctd, Brd Polyester non absorbable	0	3.5	18 cm / 7"	3
170041	Surgidac™ Ctd, Brd Polyester non absorbable	2-0	3	18 cm / 7"	3
170044	Surgidac™ Ctd, Brd Polyester non absorbable	2-0	3	120 cm / 48"	3

Endo Stitch™ Single-Stitch Loading Unit

Single Use Loading Units for use with SILS™ Stitch and Endo Stitch™ 10 mm Suturing Devices.



ORDER CODE	SUTURE MATERIAL	USP SIZE	METRIC SIZE	LENGTH	BOX QTY
170001-	Sofsilks™ Ctd, Brd Silk	0	3.5	18 cm / 7"	12
170003	Sofsilks™ Ctd, Brd Silk	0	3.5	120 cm / 48"	12
170002	Sofsilks™ Ctd, Brd Silk	2-0	3	18 cm / 7"	12
170004	Sofsilks™ Ctd, Brd Silk	2-0	3	120 cm / 48"	12
170010	Bralon™ Ctd, Brd Nylon	0	3.5	18 cm / 7"	12
170012	Bralon™ Ctd, Brd Nylon	0	3.5	120 cm / 48"	12
170011	Bralon™ Ctd, Brd Nylon	2-0	3	18 cm / 7"	12
170013	Bralon™ Ctd, Brd Nylon	2-0	3	120 cm / 48"	12
173026	Surgidac™ Ctd, Brd Polyester	0	3.5	18 cm / 7"	12
173024	Surgidac™ Ctd, Brd Polyester	0	3.5	120 cm / 48"	12
173021	Surgidac™ Ctd, Brd Polyester	2-0	3	18 cm / 7"	12
173023	Surgidac™ Ctd, Brd Polyester	2-0	3	120 cm / 48"	12
170050	Polysorb™ Ctd, Brd Syn. Abs.	0	3.5	18 cm / 7"	12
170052	Polysorb™ Ctd, Brd Syn. Abs.	0	3.5	120 cm / 48"	12
170054	Polysorb™ Ctd, Brd Syn. Abs.	0	3.5	18 cm / 7"	12
170051	Polysorb™ Ctd, Brd Syn. Abs.	2-0	3	18 cm / 7"	12
170053	Polysorb™ Ctd, Brd Syn. Abs.	2-0	3	120 cm / 48"	12
170055	Polysorb™ Ctd, Brd Syn. Abs.	2-0	3	18 cm / 7"	12
170070	Polysorb™ Ctd, Brd Syn. Abs.	3-0	2	18 cm / 7"	12
170071	Polysorb™ Ctd, Brd Syn. Abs.	3-0	2	120 cm / 48"	12
170072	Polysorb™ Ctd, Brd Syn. Abs.	3-0	2	18 cm / 7"	12
170073	Polysorb™ Ctd, Brd Syn. Abs.	3-0	2	120 cm / 48"	12
170092	Polysorb™ Ctd, Brd Syn. Abs.	4-0	1.5	18 cm / 7"	12
170090	Polysorb™ Ctd, Brd Syn. Abs.	4-0	1.5	120 cm / 48"	12
170094	Polysorb™ Ctd, Brd Syn. Abs.	4-0	1.5	18 cm / 7"	12

V-Loc™ Wound Closure Reload

Use with Endo Stitch™ and SILS™ Stitch Suturing Devices.



ORDER CODE	DESCRIPTION	USIZE	LENGTH	BOX QTY
VLOCA004	V-Loc™ 180 absorbable reload	0	10 cm / 4"	6
VLOCA006L	V-Loc™ 180 absorbable reload	0	15 cm / 6"	6
VLOCA008L	V-Loc™ 180 absorbable reload	0	20 cm / 8"	6
VLOCA204L	V-Loc™ 180 absorbable reload	2-0	10 cm / 4"	6
VLOCA206L	V-Loc™ 180 absorbable reload	2-0	15 cm / 6"	6
VLOCA208L	V-Loc™ 180 absorbable reload	2-0	20 cm / 8"	6
VLOCA304L	V-Loc™ 180 absorbable reload	3-0	10 cm / 4"	6
VLOCA306L	V-Loc™ 180 absorbable reload	3-0	15 cm / 6"	6
VLOCA308L	V-Loc™ 180 absorbable reload	3-0	20 cm / 8"	6
VLOCN004L	V-Loc™ PBT non-absorbable reload	0	10 cm / 4"	6
VLOCN006L	V-Loc™ PBT non-absorbable reload	0	15 cm / 6"	6
VLOCN008L	V-Loc™ PBT non-absorbable reload	0	20 cm / 8"	6
VLOCN204L	V-Loc™ PBT non-absorbable reload	2-0	10 cm / 4"	6
VLOCN206L	V-Loc™ PBT non-absorbable reload	2-0	15 cm / 6"	6
VLOCN208L	V-Loc™ PBT non-absorbable reload	2-0	20 cm / 8"	6
VLOCN304L	V-Loc™ PBT non-absorbable reload	3-0	10 cm / 4"	6
VLOCN306L	V-Loc™ PBT non-absorbable reload	3-0	15 cm / 6"	6
VLOCN308L	V-Loc™ PBT non-absorbable reload	3-0	20 cm / 8"	6

SUTURING ACCESSORIES

Endo Close™ Trocar Site Closure Device



ORDER CODE	DESCRIPTION	BOX QTY
173022	Endo Close™ Trocar Site Closure Device	12

Surgitie™ Single Use Ligating Loops with Delivery System



ORDER CODE	COLOR CODE	SUTURE MATERIAL	USP SIZE	METRIC SIZE	LENGTH UNITS	IN.	CM	UNITS/BOX
ED-20-L	Green	Surgitie™ Ligating Loop with Surgidac™ Coated, Braided, Polyester Suture	2-0	4	1 ×	21	53	12
ED-21-L	Green	Surgitie™ Ligating Loop with Surgidac™ Coated, Braided, Polyester Suture	0	3.5	1 ×	21	53	12
EL-20-L	Violet	Surgitie™ Ligating Loop with Delivery System with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	2-0	3	1 ×	21	53	6
EL-21-L	Violet	Surgitie™ Ligating Loop with Delivery System with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	0	3.5	1 ×	21	53	6
EL-22-LN	Violet	Surgitie™ Ligating Loop with Delivery System with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	2-0	3	1 ×	21	53	6
EL-23-LN	Violet	Surgitie™ Ligating Loop with Delivery System with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	0	3.5	1 ×	21	53	6





Surgiwip™


Single Use Suture Ligatures with Delivery System

ORDER CODE	COLOR CODE	SUTURE MATERIAL	USP SIZE	METRIC SIZE	LENGTH UNITS	IN.	CM	UNITS/ BOX
------------	------------	-----------------	----------	-------------	--------------	-----	----	------------

ELW-30	Violet	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	0	3.5	1 ×	48	120	6
---------------	---------------	--	----------	------------	------------	-----------	------------	----------

ELW-31	Violet	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	2-0	3	1 ×	48	120	6
---------------	---------------	--	------------	----------	------------	-----------	------------	----------

EST STRAIGHT TAPER



19 mm

EV-20 TAPER



1/2 26 mm

ELW-42	Violet	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	0	3.5	1 ×	48	120	6
---------------	---------------	--	----------	------------	------------	-----------	------------	----------

ELW-43	Violet	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	2-0	3	1 ×	48	120	6
---------------	---------------	--	------------	----------	------------	-----------	------------	----------

ESK SKI TAPER



23 mm

ELW-52	Violet	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	0	3.5	1 ×	48	120	6
---------------	---------------	--	----------	------------	------------	-----------	------------	----------

ELW-53	Violet	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	2-0	3	1 ×	48	120	6
---------------	---------------	--	------------	----------	------------	-----------	------------	----------

EV-20 TAPER




1/2 26 mm

EDW-50	Green	Surgiwip™ Suture Ligature with Surgidac™ Coated, Braided, Polyester Suture	0	3.5	1 ×	48	120	12
---------------	--------------	---	----------	------------	------------	-----------	------------	-----------

EDW-51	Green	Surgiwip™ Suture Ligature with Surgidac™ Coated, Braided, Polyester Suture	2-0	3	1 ×	48	120	12
---------------	--------------	---	------------	----------	------------	-----------	------------	-----------

ESK SKI TAPER



23 mm

EDW-52	Green	Surgiwip™ Suture Ligature with Surgidac™ Coated, Braided, Polyester Suture	0	3.5	1 ×	48	120	12
---------------	--------------	---	----------	------------	------------	-----------	------------	-----------

ORDER INFORMATION

Clip Appliers – page 5

Endoscopic Clip Appliers – page 6

Endo Clip™ III 5 mm Single Use Clip Applier with Clip Logic™ Technology

ORDER CODE	DESCRIPTION	BOX QTY
176630	Endo Clip™ III 5 mm 16 M/L Titanium Clips, 33 cm Shaft Length	6
176630B	Endo Clip™ III 5 mm without clip display (no battery inside)	6

Endo Clip™ II 10 mm Pistol Grip Single Use Clip Applier with Super Interlock™ Titanium Clips

ORDER CODE	DESCRIPTION	BOX QTY
176657	Endo Clip™ II 10 mm Pistol Grip 20 M/L Super Interlock™ Titanium Clips, 29 cm Shaft Length	6

Endo Clip™ 10 mm Pistol Grip Single Use Clip Applier with Super Interlock™ Titanium Clips

ORDER CODE	DESCRIPTION	BOX QTY
176625	Endo Clip™ L 10 mm Pistol Grip 15 Large Super Interlock™ Titanium Clips, 33 cm Shaft Length	6

Lapro-Clip™ Clip Applier Reusable stainless steel laparoscopic clip applier with absorbable Lapro-Clip™ clips

ORDER CODE	DESCRIPTION	BOX QTY
8886848700	Lapro-Clip™ 26,5 cm Shaft Length	1
8886848701	Lapro-Clip™ Long 33 cm Shaft Length	1

Lapro-Clip™ Absorbable Loading Units

Disposable loading units for use with the Lapro-Clip™ clip applier code 8886848700. The 2 piece absorbable clip is made of an external polyglycolic acid body that runs over a flexible internal polyglyconate part.

ORDER CODE	DESCRIPTION	BOX QTY
8886848808	Lapro-Clip™ Clip Disposable loading unit with 6 absorbable Medium/Large (8 mm) clips per sterile box.	6
8886848508	Lapro-Clip™ Multi Loading Unit Automatic disposable loading unit with 8 absorbable clips per sterile box.	5
8886848812	Lapro-Clip™ Clip Disposable loading unit with 6 absorbable Large (12 mm) clips per sterile box.	6
8886848882	Lapro-Clip™ Clip Disposable loading unit with 2 absorbable Medium-Large (8 mm) clips per sterile box.	10
8886848813	Lapro-Clip™ Clip Disposable loading unit with 2 absorbable Large (12 mm) clips per sterile box.	10

Open Clip Appliers – page 8

Premium Surgiclip™ III Single Use Automatic Clip Appliers with Super Interlock™ Titanium Clips

ORDER CODE	DESCRIPTION	BOX QTY
133650	Premium Surgiclip™ III S-9.0 20 Small Titanium Clips	6

Premium Surgiclip™ II Single Use Automatic Clip Appliers with Super Interlock™ Titanium Clips

ORDER CODE	DESCRIPTION	BOX QTY
134051	Premium Surgiclip™ II M-9.75 20 Medium Titanium Clips	6
134053	Premium Surgiclip™ II M-11.5 30 Medium Titanium Clips	6

Premium Surgiclip™ Single Use Automatic Clip Appliers with Super Interlock™ Titanium Clips

ORDER CODE	DESCRIPTION	BOX QTY
134048	Premium Surgiclip™ L-13.0 15 Large Titanium Clips	6






Laparoscopic Hand Instruments – page 9

Dissectors, Graspers, Shears – page 10


Rotulator™ 5 mm Laparoscopic Hand Instruments Single Use – 31 cm Length, 80° articulation

ORDER CODE	COLOR CODE	DESCRIPTION	BOX QTY
174213		Rotulator™ Endo Dissect™ 5 mm Single Use Instrument with Monopolar Cautery	6
174309		Rotulator™ Endo Mini-Shears™ 5 mm Single Use Instrument with Monopolar Cautery	6
174233		Rotulator™ Endo Grasp™ 5 mm Single Use Instrument with Spin Lock™	6


Laparoscopic Hand Instruments 5 mm Single Use – 31 cm Length

ORDER CODE	COLOR CODE	DESCRIPTION	BOX QTY
176643		Endo Shears™ 5 mm Single Use Instrument with Monopolar Cautery	6
174301		Endo Mini-Shears™ 5 mm Single Use Instrument with Monopolar Cautery	6
176645		Endo Dissect™ 5 mm Single Use Instrument with Monopolar Cautery	6
173030		Endo Grasp™ 5 mm Single Use Instrument with On/Off Ratchet Switch – 31 cm Length	6
174317		Endo Clinch™ II 5 mm Single Use Instrument with On/Off Ratchet Switch – 31 cm Length	6

Laparoscopic Hand Instrument Long 5 mm Single Use – 45 cm Length






ORDER CODE	COLOR CODE	DESCRIPTION	BOX QTY
174601		Endo Shears™ Long 5 mm Single Use Instrument with Monopolar Cautery	6

Laparoscopic Hand Instruments 10 mm

ORDER CODE	COLOR CODE	DESCRIPTION	BOX QTY
174001		Endo Babcock™ 10 mm Single Use Clamp – 31 cm Length	6

Retractors – page 12

Laparoscopic Retractors

ORDER CODE	COLOR CODE	DESCRIPTION	BOX QTY
176613		Endo Retract™ 10 mm Single Use Instrument – 31 cm Length, 3 Finger	6
176647		Endo Retract™ II 10 mm Single Use Instrument – 32 cm Length, 5 Finger articulating	6
174311		Endo Retract™ Maxi 10 mm Single Use Instrument – 32 cm Length	6
173046		Endo Paddle Retract™ 12 mm Single Use Instrument – 46 cm Length	3
174209		Endo Mini-Retract™ 5 mm Single Use Instrument – 31 cm Length	6
173019		Endo Peanut™ 5 mm Single Use Device – 45 cm Length (3 units per pouch – 4 pouches per box)	3

Suction and Irrigation – page 13

Surgiwand™ II 5 mm Single Use Suction/Irrigation Devices without Cautery, 34 cm Shaft Length

ORDER CODE	DESCRIPTION	BOX QTY
178081	Surgiwand™ II 5 mm Suction/Irrigation Device without Tubing, For use with the Nitrogen Pump System	3
178083	Surgiwand™ II 5 mm Suction/Irrigation Device without Tubing, For Gravity Feed System	3

Surgiwand™ II 5 mm with Cautery Single Use Suction/Irrigation Devices, 34 cm Shaft Length

ORDER CODE	DESCRIPTION	BOX QTY
178091	Surgiwand™ II 5 mm Suction/Irrigation Device With Cautery and L-Hook Tip without Tubing For use with the Nitrogen Pump System	3
178099	Surgiwand™ II 5 mm Suction/Irrigation Device With Cautery and Spatula Tip without Tubing For use with the Nitrogen Pump System	3
178093	Surgiwand™ II 5 mm Suction/Irrigation Device With Cautery, L-Hook Tip and Tubing For Gravity Feed System	3
178094	Surgiwand™ II 5 mm Suction/Irrigation Device With Cautery, Spatula Tip and Tubing For Gravity Feed System	3

Vital Vue™ Surgical Instrument

ORDER CODE	DESCRIPTION	BOX QTY
8886827906	Bulb tip – Used in many gynaecological and plastic surgery cases as well as other shallow procedures where a bulb tip is appropriate	6
8886828006	Extended bulb tip – For lower abdominal and pelvic cases where extra reach is required and a bulb tip is appropriate	6
8886828106	Extended slender tip – The tip of choice for surgeons who prefer precise suction with a slender, 1,35 cm profile tip	6
8886828206	Extended metal tip – The 12 Fr., 5 cm extended metal tip is ideal where reach and precise aspiration is required, such as oral surgery and some vascular cases	6
8886828306	Extra long extended metal tip – Designed specifically for orthopaedic cases, this 16 Fr., 8 cm metal tip is ideal where extra extension is necessary particularly in hip arthroplasty	6

Vital Vue™ Light Power Source

ORDER CODE	DESCRIPTION	REGION	BOX QTY
8886828421	Vital Vue™ light power source	UK	1
8886828422	Vital Vue™ light power source	Spain	1
8886828423	Vital Vue™ light power source	Germany	1
8886828424	Vital Vue™ light power source	Italy	1
8886828425	Vital Vue™ light power source	France	1

Specimen Retrieval Products – page 15

Single Use Specimen Retrieval Products (Latex Free)

ORDER CODE	DESCRIPTION	BOX QTY
CATCH10	Reliacatch™ Specimen Retrieval Bag 10 mm 6.2 × 16.6 cm, Volume capacity 275 mL	5
173050G	Endo Catch™ Gold 10 mm 5.8 × 15 cm Specimen Pouch, Volume capacity 190 cc, 29,5 cm Shaft Length	6
173049	Endo Catch™ II 15 mm 9.7 × 23.6 cm Specimen Pouch, Volume capacity 1500 cc, 29,5 cm Shaft Length	3

Single Use Specimen Retrieval Pouches

ORDER CODE	DESCRIPTION	BOX QTY
250-30	Endobag™ 7.6 × 15.2 cm Specimen Retrieval System – 20.9 cm Shaft Length	20
250-40	Endobag™ 12.7 × 20.3 cm Specimen Retrieval System – 20.9 cm Shaft Length	20

Automatic Suturing Devices – page 16

Endo Stitch™ Suturing Devices – page 17

Endo Stitch™ 10 mm Suturing Device

ORDER CODE	DESCRIPTION	BOX QTY
173016	Endo Stitch™ 10 mm Suturing Device, 34 cm Shaft Length	3

Endo Stitch™ Triple-Stitch Loading Unit Single Use Loading Units for use with SILS™ Stitch and Endo Stitch™ 10 mm Suturing Devices. Each Loading Unit contains three strands of suture, each attached to a 9 mm needle

ORDER CODE	SUTURE MATERIAL	USP SIZE	METRIC SIZE	LENGTH	BOX QTY
170043	Surgidac™ Ctd, Brd Polyester non absorbable	0	3.5	120 cm / 48"	3
170040	Surgidac™ Ctd, Brd Polyester non absorbable	0	3.5	18 cm / 7"	3
170041	Surgidac™ Ctd, Brd Polyester non absorbable	2-0	3	18 cm / 7"	3
170044	Surgidac™ Ctd, Brd Polyester non absorbable	2-0	3	120 cm / 48"	3

Endo Stitch™ Single-Stitch Loading Unit

Single Use Loading Units for use with SILS™ Stitch and Endo Stitch™ 10 mm Suturing Devices.

ORDER CODE	SUTURE MATERIAL	USP SIZE	METRIC SIZE	LENGTH	BOX QTY
170001-	Sofsilk™ Ctd, Brd Silk	0	3.5	18 cm / 7"	12
170003	Sofsilk™ Ctd, Brd Silk	0	3.5	120 cm / 48"	12
170002	Sofsilk™ Ctd, Brd Silk	2-0	3	18 cm / 7"	12
170004	Sofsilk™ Ctd, Brd Silk	2-0	3	120 cm / 48"	12
170010	Bralon™ Ctd, Brd Nylon	0	3.5	18 cm / 7"	12
170012	Bralon™ Ctd, Brd Nylon	0	3.5	120 cm / 48"	12
170011	Bralon™ Ctd, Brd Nylon	2-0	3	18 cm / 7"	12
170013	Bralon™ Ctd, Brd Nylon	2-0	3	120 cm / 48"	12
173026	Surgidac™ Ctd, Brd Polyester	0	3.5	18 cm / 7"	12
173024	Surgidac™ Ctd, Brd Polyester	0	3.5	120 cm / 48"	12
173021	Surgidac™ Ctd, Brd Polyester	2-0	3	18 cm / 7"	12
173023	Surgidac™ Ctd, Brd Polyester	2-0	3	120 cm / 48"	12
170050	Polysorb™ Ctd, Brd Syn. Abs.	0	3.5	18 cm / 7"	12
170052	Polysorb™ Ctd, Brd Syn. Abs.	0	3.5	120 cm / 48"	12
170054	Polysorb™ Ctd, Brd Syn. Abs.	0	3.5	18 cm / 7"	12
170051	Polysorb™ Ctd, Brd Syn. Abs.	2-0	3	18 cm / 7"	12
170053	Polysorb™ Ctd, Brd Syn. Abs.	2-0	3	120 cm / 48"	12
170055	Polysorb™ Ctd, Brd Syn. Abs.	2-0	3	18 cm / 7"	12
170070	Polysorb™ Ctd, Brd Syn. Abs.	3-0	2	18 cm / 7"	12

Endo Stitch™ Single-Stitch Loading Unit

Single Use Loading Units for use with SILS™ Stitch and Endo Stitch™ 10 mm Suturing Devices.

ORDER CODE	SUTURE MATERIAL	USP SIZE	METRIC SIZE	LENGTH	BOX QTY
170071	Polysorb™ Ctd, Brd Syn. Abs.	3-0	2	120 cm / 48"	12
170072	Polysorb™ Ctd, Brd Syn. Abs.	3-0	2	18 cm / 7"	12
170073	Polysorb™ Ctd, Brd Syn. Abs.	3-0	2	120 cm / 48"	12
170092	Polysorb™ Ctd, Brd Syn. Abs.	4-0	1.5	18 cm / 7"	12
170090	Polysorb™ Ctd, Brd Syn. Abs.	4-0	1.5	120 cm / 48"	12
170094	Polysorb™ Ctd, Brd Syn. Abs.	4-0	1.5	18 cm / 7"	12

V-Loc™ Wound Closure Reload

Use with Endo Stitch™ and SILS™ Stitch Suturing Devices..

ORDER CODE	SUTURE MATERIAL	USP SIZE	METRIC SIZE	LENGTH	BOX QTY
VLOCA004	V-Loc™ 180 absorbable reload	0	3.5	10 cm / 4"	6
VLOCA006L	V-Loc™ 180 absorbable reload	0	3.5	15 cm / 6"	6
VLOCA008L	V-Loc™ 180 absorbable reload	0	3	20 cm / 8"	6
VLOCA204L	V-Loc™ 180 absorbable reload	2-0	3	10 cm / 4"	6
VLOCA206L	V-Loc™ 180 absorbable reload	2-0	3.5	15 cm / 6"	6
VLOCA208L	V-Loc™ 180 absorbable reload	2-0	3.5	20 cm / 8"	6
VLOCA304L	V-Loc™ 180 absorbable reload	3-0	3	10 cm / 4"	6
VLOCA306L	V-Loc™ 180 absorbable reload	3-0	3.5	15 cm / 6"	6
VLOCA308L	V-Loc™ 180 absorbable reload	3-0	3.5	20 cm / 8"	6
VLOCN004L	V-Loc™ PBT non-absorbable reload	0	3	10 cm / 4"	6
VLOCN006L	V-Loc™ PBT non-absorbable reload	0	3	15 cm / 6"	6
VLOCN008L	V-Loc™ PBT non-absorbable reload	0	3.5	20 cm / 8"	6
VLOCN204L	V-Loc™ PBT non-absorbable reload	2-0	3.5	10 cm / 4"	6
VLOCN206L	V-Loc™ PBT non-absorbable reload	2-0	3.5	15 cm / 6"	6
VLOCN208L	V-Loc™ PBT non-absorbable reload	2-0	3	20 cm / 8"	6
VLOCN304L	V-Loc™ PBT non-absorbable reload	3-0	3	10 cm / 4"	6
VLOCN306L	V-Loc™ PBT non-absorbable reload	3-0	3	15 cm / 6"	6
VLOCN308L	V-Loc™ PBT non-absorbable reload	3-0	3	20 cm / 8"	6

Suturing Accessories – page 20

Endo Close™ Trocar Site Closure Device

ORDER CODE	DESCRIPTION	BOX QTY
173022	Endo Close™ Trocar Site Closure Device	12

Surgitie™ Single Use Ligating Loops with Delivery System

ORDER CODE	COLOR CODE	SUTURE MATERIAL	USP SIZE	METRIC SIZE	LENGTH UNITS	IN.	CM	UNITS/BOX
ED-20-L	Green	Surgitie™ Ligating Loop with Surgidac™ Coated, Braided, Polyester Suture	2-0	4	1 ×	21	53	12
ED-21-L	Green	Surgitie™ Ligating Loop with Surgidac™ Coated, Braided, Polyester Suture	0	3.5	1 ×	21	53	12
EL-20-L	Violet	Surgitie™ Ligating Loop with Delivery System with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	2-0	3	1 ×	21	53	6
EL-21-L	Violet	Surgitie™ Ligating Loop with Delivery System with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	0	3.5	1 ×	21	53	6
EL-22-LN	Violet	Surgitie™ Ligating Loop with Delivery System with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	2-0	3	1 ×	21	53	6
EL-23-LN	Violet	Surgitie™ Ligating Loop with Delivery System with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	0	3.5	1 ×	21	53	6

Surgiwip™ Single Use Suture Ligatures with Delivery System

ORDER CODE	COLOR CODE	SUTURE MATERIAL	USP SIZE	METRIC SIZE	LENGTH UNITS	IN.	CM	UNITS/BOX
ELW-30	Violet	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	0	3.5	1 ×	48	120	6
ELW-31	Violet	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	2-0	3	1 ×	48	120	6
ELW-33	Violet	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture without Delivery System	2-0	3	1 ×	48	120	6
ELW-42	Violet	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	0	3.5	1 ×	48	120	6
ELW-43	Violet	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	2-0	3	1 ×	48	120	6
ELW-52	Violet	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	0	3.5	1 ×	48	120	6
ELW-53	Violet	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	2-0	3	1 ×	48	120	6
EDW-50	Violet	Surgiwip™ Suture Ligature with Surgidac™ Coated, Braided, Polyester Suture	0	3.5	1 ×	48	120	12
EDW-51	Violet	Surgiwip™ Suture Ligature with Surgidac™ Coated, Braided, Polyester Suture	2-0	3	1 ×	48	120	12
EDW-52	Violet	Surgiwip™ Suture Ligature with Surgidac™ Coated, Braided, Polyester Suture	0	3.5	1 ×	48	120	12
EDW-53	Violet	Surgiwip™ Suture Ligature with Surgidac™ Coated, Braided, Polyester Suture	2-0	3	1 ×	48	120	12

Medtronic

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 20-weu-instruments-catalogue-may20-4046158-v2

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)