

# VERSAONE™ UNIVERSAL CANNULA



## PRODUCT INFORMATION

Code	Description	Ship Case Qty
UNVCA5SHF	VersaOne™ universal cannula 5 mm Short (70 mm) Fixation	6 per box
UNVCA5STF	VersaOne™ universal cannula 5 mm Standard (100 mm) Fixation	6 per box
UNVCA11STF	VersaOne™ universal cannula 11 mm Standard (100 mm) Fixation	6 per box
UNVCA12STF	VersaOne™ universal cannula 12 mm Standard (100 mm) Fixation	6 per box
UNVCA12STS	VersaOne™ universal cannula 12 mm Standard (100 mm) Smooth	6 per box
UNVCA5LGF	VersaOne™ universal cannula 5 mm Long (150 mm) Fixation	6 per box
UNVCA12LGF	VersaOne™ universal cannula 12 mm Long (150 mm) Fixation	6 per box

## MATERIALS OF MANUFACTURE

- Cannula
- Transparent polyamide
  - Polysoprene (synthetic rubber)
  - C-Flex Thermoplastic Elastomer
  - Marlex 9006 HDPE
  - Silicone oil

## DIMENSIONS

Trocar length	Short	70 mm
	Standard	100 mm
	Long	150 mm

## PACKAGING (APPLIES TO ALL CODES)

		Material	Qty	Size (cm)
Out packaging	5 mm	Carton box	1	L 23,0 x W 12,0 x D 26,5
Out packaging	11 and 12 mm	Carton box	1	L 25,5 x W 10,0 x D 7,5
Inner packaging	5, 11 and 12 mm	Soft breather Tyvek® pouch	6	

## SHIPPING INFORMATION

	Shipped case dimensions (cm)
UNVCA5SHF, UNVCA5STF	L 30,5 x W 12,7 x D 12,7
UNVCA5LGF	L 36,8 x W 12,1 x D 14,0
UNVCA11STF, UNVCA12STS, UNVCA12STF	L 22,9 x W 12,3 x D 26,7
UNVCA12LGF	L 28,2 x W 12,3 x D 26,7

## ADDITIONAL INFORMATION

Latex Free	✓
DEHP free	✓
Sterilization method	Ethylene Oxide
Shelf Life	5 years
CE Marked	Class IIa According to EU directive 93/42/EEC, Annex II.3

**IMPORTANT:** Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2017 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 16-eu-universal-cannula-tech-sheet-en-1477541

**Medtronic**  
Further, Together