

Medtronic

VERSAONE™ FAMILY TROCARS:

The access without the excess

VersaOne™ Bladeless Trocar
VersaOne™ Access System



Operating room
essentials: the complete
solution you can trust.

Efficacy & simplicity.
Every time.

The VersaOne™ Bladeless Trocar is available in 5 mm, 11 mm, 12 mm and 15 mm sizes

Smoother insertion

Unique “dolphin nose” tip facilitates easier trocar placement²

Port fixation and security

Advanced fixation ribs keep cannula from moving during surgical procedure

Outstanding clarity

Universal, clear cannula provides visibility to all tissue layers




Ergonomic design

provides control and comfort during laparoscopic procedures




Enhanced seal system

allows for easy instrument exchange¹



Unique dolphin nose tip

facilitates smooth insertion for easy trocar placement¹



Proprietary ribbed cannula

provides better port fixation and enhanced security compared to the Applied Kii™ 12 mm trocar²

VersaOne™ bladeless trocar

Sold six units per box

	5 mm Fixation	11 mm Fixation	12 mm Smooth	12 mm Fixation	15 mm Fixation
Short Length Cannula (70 mm)	NONB5SHF	-	-	NONB12SHF	-
Standard Length Cannula (100 mm)	NONB5STF	NONB11STF	NONB12STS	NONB12STF	NONB15STF
Long Length Cannula (150 mm)	NONB5LGF	-	-	NONB12LGF	NONB15LFG

VersaOne™ universal cannula

Sold six units per box

	5 mm Fixation	11 mm Fixation	12 mm Smooth	12 mm Fixation
Short Length Cannula (70 mm)	UNVCA5SHF	-	-	-
Standard Length Cannula (100 mm)	UNVCA5STF	UNVCA11STF	UNVCA12STS	UNVCA12STF
Long Length Cannula (150 mm)	UNVCA5LGF	-	-	UNVCA12LGF

Contact your Medtronic representative today for more information about the VersaOne™ Access System or visit [medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)

1. Based on internal test report#2143-123. Applies to Medtronic 12 mm trocars, when compared to the Versaport™ Bladeless 12 mm trocar. March 2013.

2. Based on internal test report #2143-114. Applies to 12 mm trocars, when compared to Kii™* (12 mm, Z-thread cannula). March 2013.

Medtronic

Important: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2022 Medtronic. Medtronic, Medtronic logo and Engineering the extraordinary are trademarks of Medtronic.

™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. EMEA-AI-2200005 - 21-weu-versaone-bladeless-info-sheet-5501643

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)