



# VERSAONE™ BLADELESS TROCAR

## PRODUCT INFORMATION

Code	Description	Ship Case Qty
NONB5SHF	VersaOne™ bladeless trocar 5 mm Short (70 mm) Fixation Cannula	6 per box
NONB12SHF	VersaOne™ bladeless trocar 12 mm Short (70 mm) Fixation Cannula	"
NONB5STF	VersaOne™ bladeless trocar 5 mm Standard (100 mm) Fixation Cannula	"
NONB11STF	VersaOne™ bladeless trocar 11 mm Standard (100 mm) Fixation Cannula	"
NONB12STS	VersaOne™ bladeless trocar 12 mm Standard (100 mm) Smooth Cannula	"
NONB12STF	VersaOne™ bladeless trocar 12 mm Standard (100 mm) Fixation Cannula	"
NONB15STF	VersaOne™ bladeless trocar 15 mm Standard (100 mm) Fixation Cannula	"
NONB5LGF	VersaOne™ bladeless trocar 5 mm Long (150 mm) Fixation Cannula	"
NONB12LGF	VersaOne™ bladeless trocar 12 mm Long (150 mm) Fixation Cannula	"
NONB15LGF	VersaOne™ bladeless trocar 15 mm Long (150 mm) Fixation Cannula	"

## MATERIALS OF MANUFACTURE

Cannula	<ul style="list-style-type: none"> <li>▪ Transparent polyamide</li> <li>▪ Polysoprene (synthetic rubber)</li> <li>▪ C-Flex Thermoplastic Elastomer</li> <li>▪ Marlex 9006 HDPE</li> <li>▪ Silicone oil</li> </ul>
Obturator	<ul style="list-style-type: none"> <li>▪ Transparent polyamide</li> <li>▪ Acrylonitrile-butadiene-styrene</li> <li>▪ Versaflex Thermoplastic Elastomer</li> </ul>

## DIMENSIONS

Trocar length	Short	70 mm
	Standard	100 mm
	Long	150 mm

## PACKAGING (APPLIES TO ALL CODES)

	Material	Qty
Out packaging	Carton box	1
Inner packaging	Soft breather Tyvek® pouch	6

## SHIPPING INFORMATION

	Shipped case dimensions (cm)
NONB5STF, NONB5SHF	L 30,5 x W 12,7 x D 12,7
NONB5LGF	L 36,8 x W 12,1 x D 14,0
NONB11STF, NONB12STF, NONB12STS, NONB12SHF, NONB15STF	L 22,9 x W 12,3 x D 26,7
NONB12LGF, NONB15LGF	L 28,2 x W 12,3 x D 26,7

**ADDITIONAL INFORMATION**

Latex Free	✓
DEHP free	✓
Sterilization method	Ethylene Oxide
Shelf Life	5 years
CE Marked	Class IIa According to EU directive 93/42/EEC, Annex II.3

**Medtronic**



Use scan app to read

For more information, please visit  
[medtronic.eu/product-catalog](http://medtronic.eu/product-catalog)

**IMPORTANT:** Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2017 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.  
All other brands are trademarks of a Medtronic company. 17-eu-versaone-bladeless-tech-sheet-en-1514856