



VERSAONE™ BLUNT TROCAR

PRODUCT INFORMATION

Code	Description	Ship Case Qty
BPT12STS	VersaOne™ blunt trocar with threaded anchor standard length smooth cannula (100 mm)	6 per box

MATERIALS OF MANUFACTURE

Cannula	<ul style="list-style-type: none"> Transparent polyamide Polysoprene (synthetic rubber) C-Flex Thermoplastic Elastomer Marlex 9006 HDPE Silicone oil
Obturator	<ul style="list-style-type: none"> Transparent polyamide Acrylonitrile-butadiene-styrene Versaflex Thermoplastic Elastomer
Threaded anchor screw	<ul style="list-style-type: none"> Nylon Acrylonitrile-butadiene-styrene Stainless steel

DIMENSIONS

Trocar length	Standard	100 mm
---------------	----------	--------

PACKAGING

	Material	Qty
Out packaging	Carton box	1
Inner packaging	Soft breather Tyvek® pouch	6

SHIPPING INFORMATION

	Shipped case dimensions (cm)
Case	L 24.8 x W 10.6 x D 16.2

ADDITIONAL INFORMATION

Latex Free	✓
DEHP free	✓
Sterilization method	Ethylene Oxide
Shelf Life	5 years
CE Marked	Class IIa According to EU directive 93/42/EEC, Annex II.3

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2017 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 17-eu-versaone-blunt-tech-sheet-en-1514777

Medtronic
Further, Together