

**Medtronic**

VERSAONE™ FAMILY TROCARS:

# The access without the excess

**VersaOne™ Blunt Trocar**  
VersaOne™ Access System



Operating room  
essentials: the complete  
solution you can trust.

Efficacy & simplicity.  
Every time.

The VersaOne™ Blunt Trocar with threaded anchor is designed for surgeons that use a Hasson or open technique. Available in 12 mm size.

#### Safe insertion\*

Blunt tip aids in reducing potential injury to internal structures upon insertion

#### Outstanding visualization

Universal, clear cannula provides visibility to all tissue layers

#### Secure placement

Threaded anchor ensures trocar placement and maintains pneumoperitoneum



#### Blunt obturator

Rounded distal tip designed to allow for a less traumatic entry

Low profile hub accommodates a variety of patients



#### Universal cannula

Increased inner diameter accommodates standard 12 mm laparoscopic instruments

Funneled entry, enhanced seal allows smooth instrument exchange<sup>1</sup>



#### Threaded anchor

Secure suture tie down with two suturing wings

Adjust anchor to desired location on cannula using a spring lever

### VersaOne™ blunt trocar with threaded anchor

Sold three units per box

Standard Length Cannula (100 mm)	12 mm Smooth
	BPT12STS

Contact your Medtronic representative today for more information about the VersaOne™ Access System or visit [medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)

\* Compared to bladed trocars

1. Based on internal test report#2143-123. Applies to Medtronic 12 mm trocars, when compared to the Versaport™ Bladeless 12 mm trocar. March 2013.

## Medtronic

Important: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2022 Medtronic. Medtronic, Medtronic logo and Engineering the extraordinary are trademarks of Medtronic.

™\* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. EMEA-AI-2200004 - 21-weu-versaone-blunt-info-sheet-5501745

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)