

Medtronic

VERSAONE™ FAMILY TROCARS:

The access without the excess

VersaOne™ Optical Trocar
VersaOne™ Access System



**Operating room
essentials: the complete
solution you can trust.**

Efficacy & simplicity.
Every time.

The VersaOne™ Optical Trocar is available in 5 mm, 11 mm, 12 mm, and 15 mm sizes

Smooth insertion

Unique 'dolphin nose' tip facilitates easier trocar placement¹

Port fixation and security

Advanced fixation ribs keep cannula from moving during surgical procedure¹

Outstanding visualization

Universal, clear cannula provides visibility to all tissue layers[†]



Low-profile housing

provides greater flexibility when placing multiple ports




Seal system

allows for easy instrument exchange²



Unique dolphin nose tip

facilitates smooth insertion for easy trocar placement¹



Proprietary ribbed cannula

provides better port fixation and enhanced security compared to the Applied Kii™ 12 mm trocar¹

VersaOne™ optical trocar

Sold six units per box

	5 mm Fixation	11 mm Fixation	12 mm Smooth	12 mm Fixation	15 mm Fixation
Short Length Cannula (70 mm)	ONB5SHF	-	-	ONB12SHF	-
Standard Length Cannula (100 mm)	ONB5STF	ONB11STF	ONB12STS	ONB12STF	ONB15STF
Standard Length Dual Pack Cannula (100 mm)	ONB5STF2C	-	-	-	-
Long Length Cannula (150 mm)	ONB5LGF	ONB11LGF	-	ONB12LGF	-

VersaOne™ universal cannula

Sold six units per box

	5 mm Fixation	11 mm Fixation	12 mm Smooth	12 mm Fixation
Short Length Cannula (70 mm)	UNVCA5SHF	-	-	-
Standard Length Cannula (100 mm)	UNVCA5STF	UNVCA11STF	UNVCA12STS	UNVCA12STF
Long Length Cannula (150 mm)	UNVCA5LGF	-	-	UNVCA12LGF

Contact your Medtronic representative today for more information about the VersaOne™ Access System or visit [medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)

[†] Images recorded on March 14, 2013, as part of a trocar clarity evaluation in an inanimate abdominal model; images representative of evaluation results

1. Based on internal test report #2143-114. Applies to 12 mm trocars, when compared to Applied Kii™* (12 mm, Z-thread cannula). March 2013.
2. Based on internal test report#2143-123. Applies to Medtronic 12 mm trocars, when compared to the Versaport™ Bladeless 12 mm trocar. March 2013.

Medtronic

Important: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2022 Medtronic. Medtronic, Medtronic logo and Engineering the extraordinary are trademarks of Medtronic.

™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. EMEA-AI-2200003 - 21-weu-versaone-optical-info-sheet-5501694

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)