



Endo Clip™ II

10 mm clip applier

PRODUCT INFORMATION

Code	Description	Qty Per Case
176657	Endo Clip™ II 10 mm Pistol Grip 20 M/L Super Interlock™ Titanium Clips, 29 cm Shaft Length	1

MATERIALS OF MANUFACTURE

301, 302, 303, 304, 305 and 17-4 P.H Stainless Steel, Polycarbonate 20% and 30% glass filled, Polycarbonate GE, Lexan, Polypropylene/Talc Filled, Polypropylene 20% Talc Rexene, Nylon 20% Glass, Titanium, Buna-N Rubber, Nylon, Polyvinylchloride film

SHIPPING INFORMATION

	Inner Box Size (mm)	Shipper Box Size (mm)
176657	L 184 × W 42 × H 542	L 551 × W 273 × H 192

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Ethylene Oxide
Shelf Life	5 years
CE Marked	Class IIb According to Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 17-weu-endo-clip-ii-10-mm-tech-sheet-1516759