



Endo-Close™

trocar site closure device

PRODUCT INFORMATION

Code	Description	Qty per case
173022	Endo Close™ Trocar Site Closure Device	12

MATERIALS OF MANUFACTURE

Polycarbonate GE LEXAN, 302, 304 Stainless Steel, 17-7 PH Stainless Steel, Nylon Natural color

SHIPPING INFORMATION

	Shipped case dimensions (mm)
Inner box	L 257 × W 118 × H 149

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Ethylene Oxide
Shelf Life	5 years
CE Marked	Class IIa According to Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 17-weu-endo-close-trocar-tech-sheet-1515887