



Endo Mini-Retract™

5 mm laparoscopic retractor

PRODUCT INFORMATION

Code	Description	Qty Per Case
174209	Endo Mini-Retract™ 5 mm, Single Use Instrument – 31 cm Length	6

MATERIALS OF MANUFACTURE

303, 304 Stainless Steel, Polycarbonate 20% Glass, Acetal A 5% Glass Filled With Lubricant, Acetal 10% PTFE, PVC Grade 4K10 STD FW:21, Nylon 6/6 Black, Tinel BB

SHIPPING INFORMATION

	Inner Box Size (mm)
174209	L 152 × W 111 × H 562

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Ethylene Oxide
Shelf Life	5 years
CE Marked	Class IIa Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 18-weu-endo-mini-retract-tech-sheet-2342479

Medtronic
Further, Together