

Endo Mini-Shears™ 5 mm Single Use Instrument with Monopolar Cautery



PRODUCT INFORMATION

Code	Description	Qty Per Case
174301	Endo Mini-Shears™ 5 mm, Single Use Instrument with Monopolar Cautery	6

MATERIALS OF MANUFACTURE

301, 302, 303, 304, And 416 Stainless Steel, Polycarbonate 10% Glass, Polycarbonate 20% Glass, Polycarbonate 20% Glass Fiber, Pvc Grade 6k10, Style Vw-1, 2%, Union Carbide, Aluminum 6061-T6

SHIPPING INFORMATION

	Inner Box Size (mm)
174301	L 610 × W 191 × H 70

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Radiation
Shelf Life	5 years
CE Marked	Class IIa According to Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 18-weu-endo-mini-shears-5-mm-tech-sheet-2341693

Medtronic
Further, Together