



Endo Paddle Retract™

single use instrument

12 mm

PRODUCT INFORMATION

Code	Description	Qty Per Case
173046	Endo Paddle Retract™ 12 mm, Single Use Instrument – 46 cm Length	3

MATERIALS OF MANUFACTURE

ABS GE, ABS, 5% Glass, 410 Stainless Steel, 301, 302, 303, 304 Stainless Steel, Polycarbonate 20% Glass Filled, Acetal, 15% Glass, 10% PTFE, Aluminum, 2011-T3 OR 6061-T6, Rubber (Synthetic) Natsyn 2200 Durometer, Nylon 20% Glass, Santoprene, RNF 150, Made From 95005-03, Polyester Knit Adele Knit Style #D4931

SHIPPING INFORMATION

	Inner Box Size (mm)	Shipper Box Size (mm)
173046	L 125 × W 39 × H 630	L 132 × W 127 × H 638

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Ethylene Oxide
Shelf Life	5 years
CE Marked	Class IIa Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 20-weu-endo-paddle-retract-tech-sheet-4749627