



Endo Stitch™ suturing device

PRODUCT INFORMATION

Code	Description	Qty Per Case
173016	Endo Stitch™ 10 mm Suturing Device, 34 cm Shaft Length	3

MATERIALS OF MANUFACTURE

173016	302, 303, 304, 17-4 P.H. Stainless Steel, 17-4 Sintered Stainless Steel
--------	---

SHIPPING INFORMATION

	Shipper Box Size (mm)	Inner Box Size (mm)
173016	L 629 × W 142 × H 92	L 130 × W 27 × H 613

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Ethylene Oxide
Shelf Life	5 years
CE Marked	Class IIa According to Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 20-weu-endo-stitch-10-mm-tech-sheet-4751206