



Fred™ Anti-Fog Solution

PRODUCT INFORMATION

Code	Description	Box Qty
220-50	FRED™ 2 Anti-Fog Kit Bulk Pack	20 bottles

MATERIALS OF MANUFACTURE

LDPE, D-SURF, R.O Water, Foam

COMPOSITION OF THE SOLUTION

Isopropyl Alcohol	<15%
Surfactant	<2%
Water	85%

SHIPPING INFORMATION

Code	Inner box size (mm)
220-50	L 108 × W 165 × H 229

ADDITIONAL INFORMATION

Latex/DEHP Free	✓
Sterilization method	Radiation
Shelf Life	5 years
CE Marked	Class I According to EU directive Directive 93/42/EEC on Medical Devices (MMD), Annex V
Flammability / Inflammability	Fred Anti-Fog Solution is classified as a flammable liquid category number 3 and non-inflammable.

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2019 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 20-weu-fred-tech-sheet-3890027

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)

Medtronic
Further, Together