



Lapro-Clip™ loading units

PRODUCT INFORMATION

Code	Description	Qty Per Case
8886848508	Lapro-Clip™ Multi Loading Unit, Automatic disposable loading unit with 8 absorbable clips per sterile box	1
8886848808	Lapro-Clip™ Clip, Disposable loading unit with 6 absorbable Medium/Large (8 mm) clips per sterile box	6
8886848812	Lapro-Clip™ Clip, Disposable loading unit with 6 absorbable Large (12 mm) clips per sterile box	6
8886848813	Lapro-Clip™ Clip, Disposable loading unit with 2 absorbable Large (12 mm) clips per sterile box	10
8886848882	Lapro-Clip™ Clip, Disposable loading unit with 2 absorbable Medium-Large (8 mm) clips per sterile box	10

MATERIALS OF MANUFACTURE

Polyglycolic Acid Polymer, Stainless Steel, Polycarbonate, LCP

SHIPPING INFORMATION

	Inner Box Size (mm)	Shipper Box Size (mm)
8886848508	L 135 × W 19 × H 565	L 139 × W 102 × H 607
8886848808	L 226 × W 113 × H 113	N/A
8886848812	L 226 × W 113 × H 113	N/A
8886848813	L 226 × W 113 × H 113	N/A
8886848882	L 226 × W 113 × H 113	N/A

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Ethylene Oxide
Shelf Life	5 Years
CE Marked	Class III According to Directive 93/42/EEC on Medical Devices (MDD), Annex II

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2021 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. EMEA-AI-2000026 - 20-weu-lapro-clip-loading-units-tech-sheet-4757441

Medtronic
Further, Together