

# Mini Step™ Bladeless Trocars



## PRODUCT INFORMATION

Code	Description	Unit of Measure	Qty
MS101003	Mini Step™ 2/3 mm cannula and dilator with radially expandable sleeve. For use with Step™ and Mini Step™ instruments only.	Box	10
MS101008	Mini Step™ 7 mm / 8 mm cannula and dilator with radially expandable sleeve and 7/8 - 5 mm reducer	Box	10
MS101005	Mini Step™ 5 mm cannula and dilator with radially expandable sleeve and 5 - 2/3 mm reducer	Box	10
MS100703	Mini Step™ Short 2 mm/3 mm cannula and dilator with radially expandable sleeve	Box	10
MS100708	Mini Step™ Short 7/8 mm cannula and dilator with radially expandable sleeve and 7/8 - 5 mm reducer	Box	10
MS100705	Mini Step™ Short 5 mm cannula and dilator with radially expandable sleeve and 5 - 2/3 mm reducer	Box	10

## MATERIALS OF MANUFACTURE

Sleeve	ABS, Nylon
Stylet	Polyurethane
Cannula	Polyurethane

## SHIPPING INFORMATION

Code	Inner Box Size (mm)
MS101003	L 256 × W 146 × H 399
MS101008	L 256 × W 146 × H 399
MS101005	L 256 × W 146 × H 399
MS100703	L 256 × W 146 × H 399
MS100708	L 257 × W 117 × H 152
MS100705	L 257 × W 117 × H 152

## ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterile	✓
Sterilization Method	Ethylene Oxide, Normal conservation and storage conditions
Shelf Life	5 years

**IMPORTANT:** Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2017 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 18-weu-mini-step-bladeless-trocars-tech-sheet-2518825

**Medtronic**  
Further, Together