



Mini Step™

Single Use Radially Expandable Access System

PRODUCT INFORMATION

Code	Description	Unit of Measure	Qty
S110000	Step™ Short Insufflation/Access Needle, 14 Gauge Compatible with Step™ and VersaStep™ Access System	Box	10
S100000	Step™ Insufflation/Access Needle, 14 Gauge Compatible with Step™ and VersaStep™ Access System	Box	10
MS100703	Mini Step™ Short 2 mm/3 mm Cannula and Dilator with Radially Expandable Sleeve	Box	10
MS100708	Mini Step™ Short 7 mm/8 mm Cannula and Dilator with Radially Expandable Sleeve and 7 mm/8 mm-5 mm Reducer	Box	10
S101000	Step™ Radially Expandable Sleeve	Box	10

MATERIALS OF MANUFACTURE

S110000	ABS, Stainless Steel, Polycarbonate
S100000	ABS, Stainless Steel, Polycarbonate
MS100703	ABS, Silicone, Polyethylene, Urethane, Polyisoprene, Polycarbonate
MS100708	ABS, Silicone, Polyethylene, Urethane, Polyisoprene, Polycarbonate

SHIPPING INFORMATION

Code	Inner Case Size (mm)
S110000	L 257 × W 117 × H 152
S100000	L 257 × W 117 × H 152
MS100703	L 256 × W 146 × H 399
MS100708	L 256 × W 146 × H 399

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterile	✓
Sterilization Method	Ethylene Oxide / Radiation
Shelf Life	5 years: S110000 and S100000 3 years: MS100703 and MS100708
CE Marked	Class IIa According to Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2017 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 18-weu-mini-step-single-use-tech-sheet-2518423