



ReliaCatch™ 10 mm specimen retrieval bag

PRODUCT INFORMATION

Code	Description	Ship Case Qty
CATCH10	Specimen Retrieval Bag 10 mm	5

MATERIALS OF MANUFACTURE

Bag	Ripstop nylon
Device tube	Polycarbonate
Device handle and body	ABS (Acrylonitrile butadiene styrene)
Device forks	301SS
Suture	Nylon
Suture loop tubing	LDPE (low density polyethylene)

DIMENSIONS

	Length	Bag Size	Outer Diameter
CATCH10	35.56 cm (working length below housing to end of fork)	275 mL	max 1 cm

PACKAGING (5 UNITS PER BOX)

Individual Unit Packaging Type	Sealed Nylon Tyvek breather soft pouch
Packaging Recycle Content (%):	Soft pouch is 100% recyclable Box is made from 100% recycled fibers

SHIPPING INFORMATION

Code	Dimensions (cm)
CATCH10	L 55.88 x W 15.24 x D 9.53

ADDITIONAL INFORMATION

Polyvinyl Chloride (PVC) Free Packaging and Device (If not does it have patient contact)	✓
Latex Free Packaging and Device (If not, does it have patient contact)	✓
DEHP Free	✓
Silicone Free	✓
Sterile – single use	✓
Is it Xray Dectactable/Radiopacity	No
Shelf Life (in months)	48 months
Sterilization Method	Ethylene oxide

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2017 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 18-weu-catch10-tech-sheet-2265313

Medtronic
Further, Together