



Surgineedle™

single use pneumoperitoneum needles

PRODUCT INFORMATION

Code	Description	Qty Per Case
172015	Surgineedle™ single use pneumoperitoneum needle 120 mm	12
172016	Surgineedle™ single use pneumoperitoneum needle 150 mm Long	12

MATERIALS OF MANUFACTURE

ABS GE, Polycarbonate GE LEXAN, 302, 304 stainless steel, 17-7 PH stainless steel, Clear ABS

SHIPPING INFORMATION

	Shipped Case Dimensions (mm)
Inner box	L 257 × W 117 × H 132

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Ethylene Oxide
Shelf Life	5 years
CE Marked	Class IIa According to Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 17-weu-surgineedle-tech-sheet-en-1515330

Medtronic
Further, Together