



# SURGISLEEVE™ WOUND PROTECTOR

## PRODUCT INFORMATION

Code	Description	Ship Case Qty
<a href="#">WPXSM24</a>	SurgiSleeve™ wound protector extra small	5 per box
<a href="#">WPSM256</a>	SurgiSleeve™ wound protector small	"
<a href="#">WPMD509</a>	SurgiSleeve™ wound protector medium	"
<a href="#">WPLGR914</a>	SurgiSleeve™ wound protector large with retraction ring	"
<a href="#">WPXLGR1117</a>	SurgiSleeve™ wound protector extra large with retraction ring	"

## MATERIALS OF MANUFACTURE

Part	Size	Material
Sleeve Film	Small, Medium and Large	PS-8011 Polyurethane
	Small	Pellethane 2363-80A
Proximal Ring	Medium	Pellethane 2363-80A
	Large	Pellethane 2363-90A
Distal Ring	Small	Pellethane 2363-90A
	Medium	Pellethane 2363-55D
	Large	Pellethane 2363-55D

## DIMENSIONS (INCISION SIZE)

Code	Incision Size (cm)
<a href="#">WPXSM24</a>	2 – 4
<a href="#">WPSM256</a>	2.5 – 6
<a href="#">WPMD509</a>	5 – 9
<a href="#">WPLGR914</a>	9 – 14
<a href="#">WPXLGR1117</a>	11 – 17

## PACKAGING (5 UNITS PER BOX)

Code	Inner Box (cm)	Outer Box (cm)
<a href="#">WPXSM24</a>	6,34 L x 5,25 W x 8,38 D	16,10 L x 63,50 W x 21,28 D
<a href="#">WPSM256</a>	6,34 L x 5,25 W x 8,38 D	16,10 L x 63,50 W x 21,28 D
<a href="#">WPMD509</a>	7,34 L x 5,25 W x 9,88 D	18,64 L x 63,50 W x 25,09 D
<a href="#">WPLGR914</a>	10,75 L x 5,25 W x 10,88 D	27,30 L x 63,50 W x 27,63 D
<a href="#">WPXLGR1117</a>	12,50 L x 6 W x 12 D	31,75 L x 15,24 W x 30,48 D

## ADDITIONAL INFORMATION

Latex Free	✓
Phthalate (DHP) free	✓
Sterilization method	Ethylene Oxide
Shelf Life	3 years
CE Marked	Class IIa According to EU directive 93/42/EEC, Annex II.3

**IMPORTANT:** Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2017 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 16-weu-surgisleeve-tech-sheet-en-1477383

**Medtronic**  
Further, Together