



# Surgiwip™

## single use suture ligatures with delivery system

### PRODUCT INFORMATION

Code	Description	Qty Per Case
<b>EDW-50</b>	Surgiwip™ Suture Ligature with Surgidac™ Coated, Braided, Polyester Suture	12
<b>EDW-51</b>	Surgiwip™ Suture Ligature with Surgidac™ Coated, Braided, Polyester Suture	12
<b>ELW-30</b>	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	6
<b>ELW-31</b>	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	6
<b>ELW-33</b>	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture without Delivery System	6
<b>ELW-42</b>	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	6
<b>ELW-43</b>	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	6
<b>ELW-53</b>	Surgiwip™ Suture Ligature with Polysorb™ Coated, Braided, Synthetic Absorbable Suture	6

### MATERIALS OF MANUFACTURE

Code	Description
<b>EDW-50</b>	Nylon, clear polycarbonate, thermoplastic elastomer, stainless steel and PET suture
<b>EDW-51</b>	Nylon, clear polycarbonate, thermoplastic elastomer, stainless steel and PET suture
<b>ELW-30</b>	Nylon, clear polycarbonate, thermoplastic elastomer, stainless steel and polysorb suture
<b>ELW-31</b>	Nylon, clear polycarbonate, thermoplastic elastomer, stainless steel and polysorb suture
<b>ELW-33</b>	Nylon, clear polycarbonate, thermoplastic elastomer, stainless steel and polysorb suture
<b>ELW-42</b>	Nylon, clear polycarbonate, thermoplastic elastomer, stainless steel and polysorb suture
<b>ELW-43</b>	Nylon, clear polycarbonate, thermoplastic elastomer, stainless steel and polysorb suture
<b>ELW-52</b>	Nylon, clear polycarbonate, thermoplastic elastomer, stainless steel and polysorb suture
<b>ELW-53</b>	Nylon, clear polycarbonate, thermoplastic elastomer, stainless steel and polysorb suture

**SHIPPING INFORMATION**

Code	Inner Box Size (mm)
EDW-50	L 197 × W 152 × H 470
EDW-51	L 197 × W 152 × H 470
ELW-30	L 156 × W 152 × H 470
ELW-31	L 156 × W 152 × H 470
ELW-33	L 156 × W 152 × H 470
ELW-42	L 156 × W 152 × H 470
ELW-43	L 156 × W 152 × H 470
ELW-52	L 156 × W 152 × H 470
ELW-53	L 156 × W 152 × H 470

**STERILIZATION METHOD**

Code	
EDW-50	Ethylene Oxide/Radiation
EDW-51	Ethylene Oxide/Radiation
ELW-30	Ethylene Oxide
ELW-31	Ethylene Oxide
ELW-33	Ethylene Oxide
ELW-42	Ethylene Oxide
ELW-43	Ethylene Oxide
ELW-52	Ethylene Oxide
ELW-53	Ethylene Oxide

**ADDITIONAL INFORMATION**

Latex Free	✓
DEHP	✓
Shelf Life	5 years
CE Marked	Class III According to Directive 93/42/EEC on Medical Devices (MDD), Annex II

**IMPORTANT:** Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2021 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. EMEA-AI-2100028 - 20-weu-surgiwip-tech-sheet-4750055

**Medtronic**