



VersaStep™

Insufflation / Access Needle

PRODUCT INFORMATION

Code	Description	Unit of Measure	Qty
VS070000	VersaStep™ Short Insufflation/Access Needle, 14 Gauge	Box	12
VS150000	VersaStep™ Long Insufflation/Access Needle 14 Gauge	Box	12

MATERIALS OF MANUFACTURE

Needle	Stainless Steel
--------	-----------------

SHIPPING INFORMATION

	Dimensions (mm)	Weight (g)
Outer Box	L 262 × B 120 × H 154	400

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterile	✓
Sterilization Method	Ethylene Oxide, Normal conservation and storage conditions
Shelf Life	5 years

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2017 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 18-weu-versastep-access-needle-tech-sheet-2519026

Medtronic
Further, Together