



VersaStep™ Plus

Single Use Radially Expandable Access System

PRODUCT INFORMATION

Code	Description	Unit of Measure	Qty
VS101000	VersaStep™ Radially Expandable Sleeve	Box	6
VSR100005	VersaStep™ Reposable Radially Expandable Sleeve with 5 mm Versaseal™	Box	6
VSR100812P	VersaStep™ Reposable Radially Expandable Sleeve 5 - 12 mm with 5-12 mm Versaseal Plus™	Box	6
VSR100800	VersaStep™ Short Reposable Radially Expandable Sleeve	Box	6
VSR107005	VersaStep™ Short Reposable 5 mm Versaseal™ with Radially Expandable Sleeve	Box	6
VSR107812P	VersaStep™ Reposable Short Radially Expandable Sleeve with 5-12 mm Versaseal Plus™	Box	6

MATERIALS OF MANUFACTURE

Sleeve	ABS, Nylon
Stylet	Polyurethane
Cannula	Polyurethane

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterile	✓
Sterilization Method	Ethylene Oxide, Normal conservation and storage conditions
Shelf Life	5 years

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2017 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 18-weu-versastep-plus-single-use-tech-sheet-2518624